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Date 1st Edition:
 7th 2009
 Date 5th Revision:
 7th 2024

DESCRIPTION

An inhibitory blood based medium for the isolation of gram (-) anaerobic bacteria.

PRINCIPLE OF THE METHOD

The combination of Vancomycin and Kanamycin enables the medium to selectively support the growth of gram-negative bacteria. The selected peptones provide optimal nutrients for the growth of anaerobic bacteria. Starch and sodium bicarbonate neutralize toxins, while haemin serves as both a key nutritional supplement and an enhancer for pigment production in *Porphyromonas melaninogenicus*. Cysteine promotes the growth of *Fusobacterium necrophorum*, *Propionibacterium acnes*, and *Bacteroides fragilis*. Arginine acts as an energy source for *Eubacterium spp.*, while soluble pyrophosphate supports the growth of *Porphyromonas gingivalis* and *Porphyromonas asaccharolyticus*. Pyruvate neutralizes hydrogen peroxide and serves as an energy source for *Veillonella spp.* Vitamin K and sodium succinate are essential factors for the development of several anaerobic microbes. The small amount of glucose (0.1%) acts as an energy source and, due to minimal acid production during its metabolism, does not significantly alter the medium's pH (low acid production).

FORMULA	g/litre
Peptone mix	23.0
Sodium chloride	5.0
Soluble starch	1.0
Sodium bicarbonate	0.4
Glucose	1.0
Sodium pyruvate	1.0
Cysteine HCl monohydrate	0.5
Haemin	0.01
Vitamin K	0.001
L-Arginine	1.0
Soluble pyrophosphate	0.25
Sodium succinate	0.5
Agar No. 2	12.0
Kanamycin	7,5mg
Vancomycin	10mg
Horse Blood	60ml

Appearance: Red - crimson, not clear due to the addition of blood

Final pH 7.2 ± 0.2 .

PRECAUTIONS

ANAEROBE CDC – ANAEROBE BLOOD AGAR W/ VANCOMYCIN & NEOMYCIN is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

Petri dishes should be stored at **2–12°C** in their packaging until ready for use. Prolonged storage at temperatures below **2°C** can cause excess moisture to accumulate in the medium, increasing the risk of contamination. Freezing, even briefly, destroys the medium. Avoid excessive heating as well. Petri dishes can be used up to the expiration date printed on the label. When you open the airtight package, if you have any leftover plates, store them in the bag until the expiration date. For transport, our stability studies indicate that the dishes can remain at **18–25°C** for **3 days** or at **25–40°C** for **24 hours** without affecting product performance.

USAGE

Unsealing in a clean environment (Laminar flow), wearing gloves. Drying the petri dish in the incubator (37°C) for 30–45 minutes. Inoculation of the sample as soon as possible after collection, with successive dilutions for isolated colonies. Incubation under anaerobic conditions at 35–37°C for 24 hours up to five days. Anaerobic conditions can be monitored with a specific anaerobiosis indicator, such as the GasPak anaerobic condition indicator.

INTERPRETATION OF RESULTS

Bacteroides fragilis grows relatively easily (24-48 hours) and forms white gray, occasionally mucous, circular, smooth colonies (1-3 mm) without a hemolysis zone.


Clostridium perfringens grows easily (24 hours) and forms medium-sized gray, mucous colonies (2-5 mm) with a β -hemolysis zone.

LIMITATIONS OF THE METHOD

Final identification must be performed using biochemical and serological tests (e.g., Latex Agglutination Test), which can be carried out directly from the suspected colonies. Anaerobic bacteria are highly sensitive to oxygen exposure. It is essential to always use an anaerobic indicator to ensure anaerobic conditions are maintained. Additionally, the sample must be placed under anaerobic conditions promptly.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Germ	ATCC	Growth	Colonies
<i>Bacteroides fragilis</i>	25285	Good	No haemolysis
<i>Clostridium perfringens</i>	13124	Good	β haemolysis
<i>Peptostreptococcus spp</i>	27337	It doesn't grow	

Clostridium perfringens

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

ANAEROBE BLOOD AGAR W/ KANAMYCIN & VANCOMYCIN - CND: W0104010499 - 

ITEM	CODE	PACKAGE	STORAGE	SELF LIFE
Petri dish 9cm	010220	10 pieces	2 – 12 °C	2 months

Produced in Greece by the company Bioprep in accordance with the requirements of the European Directive 2017/746.

BASIC UDI-DI: 5212037714010401WF. EDMA (14 01 04 01) non-chromogenic media (Plates).

The Bioprep company has been certified according to the standards: EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004

BIBLIOGRAPHY

Brazier, J.S. (1986). Yellow fluorescence of Fusobacteria Letters in Applied Microbiol. 2: 124-126.

Brazier, J.S. (1986). A note on ultraviolet red fluorescence of anaerobic bacteria in vitro. J. Appl. Bact. 60: 121-126.

Eley, A., Clarry, T., Bennett, K.W. (1989). Selective and differential medium for isolation of *Bacteroides ureolyticus* from clinical specimens. European Journal of Clinical Microbiology, Infectious Diseases. 8: 83-85.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr