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Esteve Z. 1984, Lab Med., 15:258

BACTEROIDES BILE ESCULINE AGAR (BBE)

Anderson, N.L., et al. Cumitech 3B; Quality Systems in the Clinical Microbiology Laboratory, Coordinating ed., A.S. Weissfeld. American Society for Microbiology, Washington, D.C.

Tille, P., et al. Bailey and Scott's Diagnostic Microbiology, C.V. Mosby Company, St. Louis, MO.

Facklam, R.R. and M.D. Moody. 1970. Appl. Microbiol.; 20:245

Isenberg, H.D. Clinical Microbiology Procedures Handbook, Vol. I & II. American Society for Microbiology, Washington, D.C.

MacFaddin, J.F. 1985. Media for Isolation, Cultivation, Identification, Maintenance of Bacteria, Vol. I. Williams & Wilkins, Baltimore, MD.

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6th 2024

DESCRIPTION

Divided Petri dish ANAEROBE C.N.A. – BACTEROIDES BILE ESCULINE AGAR (BBE) for the isolation and identification of anaerobic gram (+) cocci (C.N.A.) isolation and identification of *Bacteroides fragilis* (BBE).

PRINCIPLE OF THE METHOD

In ANAEROBE C.N.A., the combination of peptone mixes and beef extract provides a broad range of essential nutrients. Cornstarch also serves as an energy source enriched with vitamins. Sodium chloride maintains osmotic balance. Horse blood is an important growth factor and is used to assess the haemolytic properties of bacteria. Haemin and vitamin K1 supply additional growth factors and are essential components for anaerobic bacteria. L-cysteine and dithiothreitol act as reducing agents. The addition of colistin and nalidixic acid facilitates the isolation of gram-positive anaerobic cocci from samples containing mixed flora, specifically inhibiting the growth of gram-negative anaerobic bacteria.

FORMULA	g/litre
Peptone Mixture	20.00
Sodium Chloride	5.00
Beef Extract	3.00
Corn Starch	1.00
Nalidixic Acid	0.015
Colistin Sulfate	0.01
D-dithiothreitol	0,1
L-Cystein HCL	0,5
Vitamin K1	0.01
Hemin	0.01
Bacteriological Agar	15.0
Horse Blood Defibrinated	50ml

Appearance: Red - crimson, not clear due to the addition of blood.

Final pH 7.3 ± 0.2 in 25 °C.

B.B.E. combines simultaneous inhibition of facultative anaerobes as well as most gram-negative anaerobic microorganisms. This differentiation is achieved due to the presence of gentamicin and bovine bile included in the medium's components. *Bacteroides fragilis* hydrolyzes esculin, producing esculetin and dextrose. Esculetin reacts with the iron salts present in the medium, forming a dark brown to black complex that surrounds the colonies of *B. fragilis* on the medium.

FORMULA	g/litre
Peptone	8.0
Bile salts	20.0
Ferric citrate	0.5
Aesculin	1.0
Agar	15.0
Sodium Chloride	5.0
Beef Extract	5.0
Bacteriological Peptone	7.0
Lysed Horse Blood	10ml
Vitamin K1	5mg
Hemin	50mg
Gentamicin	70mg

Appearance: Agar light brown, red clear.

Final pH 7.1 ± 0.2 in 25°C.

PRECAUTIONS

AEROMONAS BLOOD AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of

deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

Petri dishes should be stored at **2 – 12 °C** in their packaging until ready for use. Prolonged storage at temperatures below **2°C** can cause excessive moisture in the medium, increasing the risk of contamination. Freezing, even briefly, destroys the medium. Avoid excessive heating as well. Petri dishes can be used until the expiration date indicated on the label. When you open the airtight package, if you have any leftover plates, store them in the bag until the expiration date. For transport, our stability studies have shown that the dishes can remain at **6 – 25°C** for **4 days** or at **25 – 40°C** for **48 hours** without affecting the product's performance.

USAGE

Unsealing should be done in a clean environment (Laminar flow) while wearing gloves. Place the petri dishes in the incubator at 37°C for 30–45 minutes to ensure no moisture remains on the surface of the medium. Inoculate the sample as soon as possible after collection, using successive dilutions to isolate single colonies. Incubate under anaerobic conditions at 35–37°C for 24 hours to five days. Anaerobic conditions can be monitored using a specific anaerobiosis indicator, such as the GasPak anaerobic condition indicator.

INTERPRETATION OF RESULTS

ANAEROBE CNA:

Peptostreptococcus anaerobius grows within 24–48 hours, forming small white-gray colonies (1–2 mm) with a β-hemolysis zone.


BACTEROIDES BILE ESCULINE AGAR:

After 48 hours of incubation, most anaerobes, except *Bacteroides fragilis*, are inhibited. Colonies of *B. fragilis* have a diameter greater than 1 mm, appear gray, circular, and raised, with blackening of the medium around the colonies due to esculin hydrolysis.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

ANAEROBE C.N.A.


Germ	ATCC	Colonies
<i>Peptostreptococcus anaerobius</i>	27337	B hemolysis
<i>Escherichia coli</i>	25922	Suspension (partial or total)

Peptostreptococcus anaerobius

BACTEROIDES BILE ESCULINE AGAR (BBE)

Germ	ATCC	Growth	Colonies
<i>Bacteroides fragilis</i>	25285	Good	Black color on the material.
<i>Clostridium perfringens</i>	13124	No growth	
<i>Proteus mirabilis</i>	12453	No growth	

Bacteroides fragilis ATCC 25285

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

ANAEROBE C.N.A. – BACTEROIDES BILE ESCULINE AGAR (BBE) - 

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Petri dish halved 9cm	020224	10 items	2 – 12 °C	2 months

It is produced in Greece from the Bioprep company in accordance with the requirement of the European guidance 98/79/EK. Code by EDMA 14 01 04 90. The Bioprep company has been certified according to the standards EN ISO 9001:2008 / ΔY86/1348/2004.

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45:502-504. Ruoff, k., I. 1995 Streptococcus.

American Public Health Association (1950). Diagnostic Procedures and Reagents. 3rd edn. A.P.H.A., New York.

American Public Health Association (1946). Standard Methods for the examination of Water and Sewage. 9th edn. A.P.H.A., New York.

Heginbotham M., Fitzgerald T.C., and Wade W.G. (1990). Comparison of solid media for the culture of anaerobes. J. Clin. Path. 43: 253-256.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr