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Eley, A., Clarry, T., Bennett, K.W. (1989). Selective and differential medium for isolation of Bacteriodes ureolyticus from clinical specimens. European Journal of Clinical Microbiology, Infectious Diseases. 8: 83-85.

Wade W. Griffiths, M. (1987). Comparison of Media for cultivation of subgingival bacteria. J. Dent. Res. 66: no. 4 abstract 334.

Heginbotham M., Fitzgerald T.C., and Wade W.G. (1990). Comparison of solid media for the culture of anaerobes. J. Clin. Path. 43: 253-256.

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**TECHNICAL DATA SHEET**

PRODUCT: **ANAEROBE C.N.A. BLOOD AGAR**
 REFERENCE: **010003**

Date 1st Edition:
 7th 2009
 Date 4th Revision:
 6th 2024

DESCRIPTION

ANAEROBE C.N.A. BLOOD AGAR is used for the selective isolation of anaerobic streptococci.

PRINCIPLE OF THE METHOD

The combination of peptone mixes and beef extract provides a broad range of essential nutrients. Cornstarch also serves as an energy source enriched with vitamins. Sodium chloride maintains osmotic balance. Horse blood is a significant growth factor and is also used to assess the hemolytic properties of bacteria. Haemin and vitamin K1 supply growth factors and are essential components for anaerobic bacteria. L-cysteine and dithiothreitol act as reducing agents. The addition of colistin and nalidixic acid facilitates the isolation of gram-positive anaerobic cocci from samples containing mixed flora. Specifically, they inhibit the growth of gram-negative anaerobic bacteria. The addition of colistin and nalidixic acid facilitates the isolation of gram-positive anaerobic cocci from samples containing mixed flora. Specifically, they inhibit the growth of gram-negative anaerobic bacteria.

COMPOSITION	g/litre
Peptone Mixture	20.00
Sodium Chloride	5.00
Beef Extract	3.00
Corn Starch	1.00
Nalidixic Acid	0.015
Colistin Sulfate	0.01
D-dithiothreitol	0,1
L-Cystein HCL	0,5
Vitamin K1	0.01
Hemin	0.01
Bacteriological Agar	15.0
Horse Blood Defibrinated	50ml

Appearance: red – crimson not clear due to the addition of blood.

Final pH 7.3 ± 0.2 in 25 °C.

PRECAUTIONS

ANAEROBE C.N.A. BLOOD AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

Petri dishes should be stored at **2–12°C** in their packaging until ready for use. Prolonged storage at temperatures below **2°C** can cause excess moisture to form in the medium, increasing the risk of contamination. Freezing, even momentarily, destroys the medium. Avoid excessive heating as well. Petri dishes can be used until the expiration date indicated on the label. When you open the airtight package, if you have any leftover plates, store them in the bag until the expiration date. For transport, stability studies have shown that the dishes can remain at **18–25°C** for **3 days** or at **25–40°C** for **24 hours** without affecting the product's performance.

USAGE

Place the petri dishes in the incubator at 37°C for 30–45 minutes to ensure no moisture remains on the surface of the medium. Inoculate the sample as soon as possible after collection, using serial dilutions to isolate individual colonies. Incubate under anaerobic conditions at 35–37°C for 24 hours to five days. Anaerobic conditions can be monitored using a specific anaerobic indicator, such as the GasPak anaerobic condition indicator.

INTERPRETATION OF RESULTS


Peptostreptococcus anaerobius grows within 24–48 hours, forming small white-gray colonies (1–2 mm) with a β -hemolysis zone.

LIMITATIONS OF THE METHOD

Final identification should be performed using biochemical and serological tests (e.g., Latex Agglutination Test), which can be carried out directly from the suspected colonies. Anaerobic bacteria are highly sensitive to oxygen exposure. It is essential to always use an anaerobic indicator to ensure anaerobic conditions are maintained. Additionally, the sample must be placed under anaerobic conditions promptly.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganisms	ATCC	Colonies
<i>Peptostreptococcus anaerobius</i>	27337	B hemolysis
<i>Escherichia coli</i>	25922	Suspension (partial or total)

Peptostreptococcus anaerobius

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

ANAEROBE C.N.A. BLOOD AGAR – CND: W0104010499 - 

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Petri dish 90mm	010003	10 pieces	2 – 12 °C	2 months

It is produced in Greece from the Bioprep company in accordance with the requirement of the European guidance 98/79/EK. Code by EDMA 14 01 04 90. The Bioprep company has been certified according to the standards EN ISO 9001:2008 / ΔY86/1348/2004.

BIBLIOGRAPHY

- Brazier, J.S. (1986). Yellow fluorescence of Fusobacteria Letters in Applied Microbiol. 2: 124-126.
- Brazier, J.S. (1986). A note on ultraviolet red fluorescence of anaerobic bacteria in vitro. J. Appl. Bact. 60: 121-126.
- Eley, A., Clarry, T., Bennett, K.W. (1989). Selective and differential medium for isolation of Bacteriodes ureolyticus from clinical specimens. European Journal of Clinical Microbiology, Infectious Diseases. 8: 83-85.
- Wade W. Griffiths, M. (1987). Comparison of Media for cultivation of subgingival bacteria. J. Dent. Res. 66: no. 4 abstract 334.
- Heginbotham M., Fitzgerald T.C., and Wade W.G. (1990). Comparison of solid media for the culture of anaerobes. J. Clin. Path. 43: 253-256.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr