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E-mail: bioprep1@otenet.gr www.bioprepare.gr

TECHNICAL DATA SHEET

PRODUCT: **ANAEROBE CDC BLOOD AGAR (FASTIDIOUS)**

REFERENCE: **010002**

Date 1st Edition:

7th 2009

Date 5th Revision:

8th 2025

DESCRIPTION

A primary isolation medium suitable for most clinically significant anaerobic organisms.

PRINCIPLE OF THE METHOD

ANAEROBE CDC BLOOD AGAR (FASTIDIOUS) has been proven superior to other formulations as a primary isolation medium for fastidious anaerobic bacteria. The selected peptones provide optimal nutrients for anaerobic bacterial growth. Starch and sodium bicarbonate neutralize toxins, while haemin, in addition to being a key nutritional supplement, aids pigment production in *Porphyromonas melaninogenicus*. Cysteine promotes the growth of *Fusobacterium necrophorum*, *Propionibacterium acne*, and *Bacteroides fragilis*. Arginine serves as an energy source for *Eubacterium* spp. Soluble pyrophosphate supports *Porphyromonas gingivalis* and *Porphyromonas asaccharolyticus*. Pyruvate helps neutralize hydrogen peroxide and serves as an energy source for *Veillonella* spp. Vitamin K and sodium succinate are essential for the growth of various anaerobic bacteria. The low glucose concentration (0.1%) provides energy without significantly altering the medium's pH due to minimal acid production.

COMPOSITION	g/litre
Peptone mix	23.0
Sodium chloride	5.0
Soluble starch	1.0
Sodium bicarbonate	0.4
Glucose	1.0
Sodium pyruvate	1.0
Cysteine HCl monohydrate	0.5
Haemin	0.01
Vitamin K	0.001
L-Arginine	1.0
Soluble pyrophosphate	0.25
Sodium succinate	0.5
Agar No. 2	12.0
Horse Blood	60ml

Appearance: Red to burgundy, opaque due to the addition of blood.

Final pH 7.2 ± 0.2 at 25 °C.

PRECAUTIONS

ANAEROBE CDC BLOOD AGAR is an in vitro diagnostic laboratory material and should only be handled by trained laboratory personnel. This material contains peptones and extracts of animal origin. Certificates regarding the source and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious and handled using standard safety precautions (avoid ingestion and inhalation). Always handle plates wearing gloves and within a Class II laminar flow cabinet to prevent contamination, especially by saprophytic fungi. Do not use plates that show signs of microbial contamination. The material should not be used after the expiration date. In case of skin contact, wash immediately with plenty of water and soap. Positive samples must be disposed of according to the appropriate hygiene regulations for infectious materials.

STORAGE AND TRANSPORT CONDITIONS

The plates must be stored at **2–12 °C** in their packaging until the time of use. Prolonged storage below **2 °C** generates excess moisture in the medium, posing a risk of contamination. Even brief freezing irreversibly damages the material. Excessive heating should also be avoided. Plates may be used until the expiration date indicated on the label. If the airtight packaging is accidentally opened, the plate may be stored in the refrigerator for **5–7 days**, sealed with parafilm or placed in a bag. Stability studies have shown that the plates can remain at **18–25 °C** for **3 days** or at **25–40 °C** for **48 hours** without affecting product performance.

INSTRUCTIONS FOR USE

Unseal in a perfectly clean area (Laminar flow), using gloves. Dry the plate in an incubator at 37 °C for 30–45 minutes. Inoculate the sample as soon as possible after collection, using serial dilutions to isolate single colonies. Incubate under anaerobic conditions at 35–37 °C for 24 hours up to five days. Anaerobic conditions should be verified with a special indicator such as the GasPak anaerobic indicator.

RESULT INTERPRETATION

Bacteroides fragilis grows relatively easily (24–48 hours), forming white-grey, sometimes mucoid, round smooth colonies (1–3 mm) without a hemolysis zone.

Clostridium perfringens grows rapidly (24 hours), forming medium-sized grey, mucoid colonies (2–5 mm) with a zone of β -hemolysis.

METHOD LIMITATIONS


Final identification must be confirmed with biochemical and serological tests (e.g., Latex agglutination test) and can be performed directly from suspected colonies. Anaerobic bacteria are highly sensitive to oxygen exposure. Always ensure an anaerobic indicator is present to verify anaerobic conditions. Samples must be placed under anaerobic conditions promptly.

GENERAL QUALITY CONTROL CHARACTERISTICS


Microorganism	Growth	Colonies
<i>Bacteroides fragilis</i> ATCC 25285	Good	No hemolysis
<i>Clostridium perfringens</i> ATCC 13124	Good	β -hemolysis

Bacteroides fragilis ATCC 25285

Clostridium perfringens ATCC 13124

Clostridium difficile ATCC 9689

DISPOSAL OF MATERIALS IN WASTE

Materials that show no growth can be considered non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of in accordance with guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must manage and dispose of (or assign the management and disposal of) it according to the applicable regulations.

SPECIFICATIONS

ANAEROBE CDC BLOOD AGAR – CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Petri Dish 9cm	010002	10 pieces	2 – 12 °C	2 months
Split Petri Dish 9cm	020002	10 pieces	2 – 12 °C	2 months

Manufactured in Greece by Bioprep in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 5212037714010490X8 EDMA Classification: (14 01 04 90) Other Prepared Media in Plates. Bioprep is certified according to the standards: EN ISO 9001:2015 / EAOT EN ISO 13485:2016 / DY8d/1348/2004..

REFERENCES

- Brazier, J.S. (1986). Yellow fluorescence of *Fusobacteria* Letters in Applied Microbiol. 2: 124-126.
Brazier, J.S. (1986). A note on ultra violet red fluorescence of anaerobic bacteria in vitro. J. Appl. Bact. 60: 121-126.
Eley, A., Clarry, T., Bennett, K.W. (1989). Selective and differential medium for isolation of *Bacteroides ureolyticus* from clinical specimens. European Journal of Clinical Microbiology, Infectious Diseases. 8: 83-85.
Wade W. Griffiths, M. (1987). Comparison of Media for cultivation of subgingival bacteria. J. Dent. Res. 66: no. 4 abstract 334.
Heginbotham M., Fitzgerald T.C., and Wade W.G. (1990). Comparison of solid media for the culture of anaerobes. J. Clin. Path. 43: 253-256.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr