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Eley, A., Clarry, T., Bennett, K.W. (1989). Selective and differential medium for isolation of *Bacteroides ureolyticus* from clinical specimens. European Journal of Clinical Microbiology, Infectious Diseases. 8: 83-85.

Wade W. Griffiths, M. (1987). Comparison of Media for cultivation of subgingival bacteria. J. Dent. Res. 66: no. 4 abstract 334.

Heginbotham M., Fitzgerald T.C., and Wade W.G. (1990). Comparison of solid media for the culture of anaerobes. J. Clin. Path. 43: 253-256.

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TECHNICAL DATA SHEET

PRODUCT: **BACTEROIDES BILE ESCULINE AGAR**
REFERENCE: **010005**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

B.B.E. is used for the initial isolation and identification of *Bacteroides fragilis*.

PRINCIPLE OF THE METHOD

The B.B.E. combines the simultaneous inhibition of the growth of potentially anaerobic as well as most gram (-) anaerobic microorganisms. This distinction is achieved thanks to the presence of gentamicin and bovine bile included in the components of the material. *B. fragilis* hydrolyzes esculin to produce esculetine and dechirose. Esculetin reacts with the iron salts present inside, producing a dark brown as a black complex that surrounds the *B. fragilis* colonies in the middle.

COMPOSITION	g/litre
Peptone	8.0
Bile salts	20.0
Ferric citrate	0.5
Aesculin	1.0
Agar	15.0
Sodium Chloride	5.0
Beef Extract	5.0
Bacteriological Peptone	7.0
Lysed Horse Blood	10ml
Vitamin K1	5mg
Hemin	50mg
Gentamicin	70mg

Appearance: Agar brown - red bright.

Final pH 7.1 ± 0.2 in 25°C.

PRECAUTIONS

The BACTEROIDES BILE ESCULINE AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

Petri dishes must be stored in **6 – 12 °C** in their packaging until the time of their use. Prolonged storage at a temperature below **6 °C** creates enough moisture in the material at risk of contamination. Freezing even instantly destroys the material. Avoid overheating also. Petri dishes must be used until the expiration date which is shown on the label. When you open the airtight package, if you have any leftover plates, store them in the bag until the expiration date. For transport, our stability studies showed that the plates may remain at **6 - 25 °C** for **4 days** or at **25 - 40 °C** for **48 hours**, without affecting the attribution of the product.


USAGE

Laminar flow sealing with gloves. Drying the plate in the incubator (37 °C) for 30 - 45 minutes. Vaccinate the sample as soon as possible after taking it, with successive dilutions for individual colonies. Incubation under anaerobic conditions at 35-37 °C for 24 -48 hours. Anaerobic conditions can be controlled by a specific anaerobic index such as the Gas and anaerobic conditions indicator. After a 48 hour incubation, the *B. fragilis* colonies should have a diameter greater than 1mm and they have grey colour, thh colonies are also circular and swollen. In the final identification also helps the preparation gram staining of the colonies. Most anaerobes, other than *B. fragilis*, have been suspended. The hydrolysis of esculin is indicated by the tanning of the nutrient mean around the colonies.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Germ	ATCC	Growth	Colonies
<i>Bacteroides fragilis</i>	25285	Good	Black color on the material
<i>Clostridium perfringens</i>	13124	It does not grow	

<i>Proteus mirabilis</i>	12453	It does not grow	
--------------------------	-------	------------------	--

Bacteroides fragilis ATCC 25285

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

BACTEROIDES BILE ESCULINE AGAR - CND: W0104010499 - 

ITEM	CODE	PACKAGE	KEEPING	LIFE SPAN
Petri dish 90mm	010005	10 pieces	2 – 12 °C	2 months

It is produced in Greece from the Bioprep company in accordance with the requirement of the European guidance 98/79/EK. Code by EDMA 14 01 04 90. The Bioprep company has been certified according to the standards EN ISO 9001:2008 / ΔY85/1348/2004.

BIBLIOGRAPHY

- Brazier, J.S. (1986). Yellow fluorescence of Fusobacteria Letters in Applied Microbiol. 2: 124-126.
 Brazier, J.S. (1986). A note on ultra violet red fluorescence of anaerobic bacteria in vitro. J. Appl. Bact. 60: 121-126.
 Eley, A., Clarry, T., Bennett, K.W. (1989). Selective and differential medium for isolation of *Bacteriodes ureolyticus* from clinical specimens. European Journal of Clinical Microbiology, Infectious Diseases. 8: 83-85.
 Wade W. Griffiths, M. (1987). Comparison of Media for cultivation of subgingival bacteria. J. Dent. Res. 66: no. 4 abstract 334.
 Heginbotham M., Fitzgerald T.C., and Wade W.G. (1990). Comparison of solid media for the culture of anaerobes. J. Clin. Path. 43: 253-256.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr