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**TECHNICAL DATA SHEET**

PRODUCT: **BILE ESCULINE AZIDE WITH VANCOMYCIN**

REFERENCE: **010274**

Date 1st Edition:

7th 2009

Date 4th Revision:

6th 2024

**DESCRIPTION**

BILE ESCULINE AZIDE WITH VANCOMYCIN is used to isolate and identify the Vancomycin-resistant *Enterococcus*.

**PRINCIPLE OF THE METHOD**

Bile Esculine Azide Agar is a modification of Bile Esculine Agar. The modification consists of adding Sodium azide and reducing the salt concentration. The result is greater selectivity and better growth of *Enterococcus*. The ability of hydrolysing esculin in the presence of salts is characteristic of the group B enterococcus and streptococci. Sodium Azide inhibits the growth of gram-negative bacteria as well as gram-positive bacteria. Vancomycin acts against gram (+) strains other than resistant *Enterococcus*.

| FORMULA                | g/litre |
|------------------------|---------|
| Tryptone               | 17.0    |
| Peptone                | 3.0     |
| Yeast extract          | 5.0     |
| Bile salts             | 10.0    |
| Sodium chloride        | 5.0     |
| Esculin                | 1.0     |
| Amonium-Ferric citrate | 0.5     |
| Sodium azide           | 0.15    |
| Vancomycin             | 5mg     |
| Agar                   | 15.0    |

Appearance: Agar gray - greenish yellow.

Final pH 7.2 ± 0.2 in 25 °C.

**PRECAUTIONS**

BILE ESCULINE AZIDE WITH VANCOMYCIN is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

**STORAGE AND TRANSPORTS CONDITIONS**

Petri dishes must be stored in **2 - 12°C** in their packaging until the time of their use. Prolonged storage at a temperature below **2 °C** creates enough moisture in the material at risk of contamination. Freezing even instantly destroys the material. Avoid overheating also.

Petri dishes must be used until the expiration date which is shown on the label. When you open the airtight package, if you have any leftover plates, store them in the bag until the expiration date. For transport, our stability studies showed that the plates may remain at **14-25°C** for **3 days** or at **27-40 °C** for **24 hours**, without affecting the attribution of the product.

**USAGE**

Allow the plates to warm to room temperature. The surface of the agar must be dry prior to vaccination. Inoculate the specimen at one end of the plate and then make consecutive liner coatings on parallel lines to create individual colonies. Incubate the aerobic plates at 35 - 37 °C for 24 to 48 hours.

**INTERPRETATION OF RESULTS**

*VRE faecium* & *VRE faecalis* form colonies of small gray with tan material around the colonies.

**LIMITATIONS OF THE METHOD**

Final identification may require additional testing.

Some rare varieties of *Lactobacilli* & *Pediococcus* sometimes develop purple colonies. However, they can be differentiated by PYR control (1-pyrrolidonyl β-naftylamide):

PYR (+) -> *VRE Enterococci*.

PYR (-) -> *Lactobacillus* & *Pediococcus*.

After 24 hours of incubation, some rare varieties of *E. gallinarum* may sometimes develop purple colonies.

The R-type of vancomycin should be determined and confirmed by further tests.

The Oxidase Test Strips-Code MID-61G is recommended in pure culture colonies for full recognition. Final identification may require additional tests such as biochemical or immunoassay (latex agglutination).

#### GENERAL CHARACTERISTICS OF QUALITY CONTROL

| Microorganism                                   | Growth - Colony Color                        |
|-------------------------------------------------|----------------------------------------------|
| <i>VRE Enterococcus faecalis</i><br>ATCC® 51299 | Brown-black, small<br>Esculin Hydrolysis (+) |
| <i>Enterococcus faecalis</i><br>ATCC® 29212     | Inhibited                                    |
| <i>Escherichia coli</i><br>ATCC® 25922          | Inhibited                                    |
| <i>Candida albicans</i><br>ATCC® 10231          | Inhibited                                    |
| <i>Staphylococcus aureus</i><br>ATCC® 25923     | Inhibited                                    |

#### WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

#### SPECIFICATIONS

##### BILE ESCULINE AZIDE WITH VANCOMYCIN – CE

| PRODUCT   | CODE   | PACKING   | STORE     | SELF LIFE |
|-----------|--------|-----------|-----------|-----------|
| Plate 9cm | 010274 | 10 pieces | 2 – 12 °C | 2 months  |

Produced in Greece by Bioprep in accordance with the requirements of the European directive 98/79/EK. FEK B2198/2-10-2009. Code EDMA 14 01 02 90. Bioprep company has been certified according to standards EN ISO 9001:2015 / EAOT EN ISO 13485:2016.

#### BIBLIOGRAPHY

ISO Standard 7899-2 (2000) water quality detection and enumeration of intnal enterococci. Part 2: Membrane method.

ATLAS, R.M. & L.C. PARKS (1993) Handbook oo microbiological media. CRC Press Boca Raton. Fla.

ISENBERG, H.D., D. GOLDBREG & J. SAMPSON (1970) Laboratory studies with a selective entroccoccus medium. Appl. Microbiol. 20:433.

BRODSKY M.H. & D.A. SCHIEMANN (1976) Evaluation of Pfizer Selective Enterococcus and KF media for recovery of fecal streptococci from water by membrane filtration. Appl. Environ. Microbiol. 31:695-699.

#### IN VITRO MANUFACTURER'S DATA

#### G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: [bioprep1@otenet.gr](mailto:bioprep1@otenet.gr)

[www.bioprep.gr](http://www.bioprep.gr)