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**TECHNICAL DATA SHEET**

PRODUCT: **BLOOD AGAR – CHROMagar™ ORIENTATION**
 REFERENCE: **020211**

Date 1st Edition:
7th 2009
 Date 4th Revision:
6th 2024

DESCRIPTION

The petri dish halved BLOOD AGAR – CHROMagar™ ORIENTATION is used to cultivate all aerobic bacteria and identify hemolytic gram (+) grains (BLOOD AGAR). Also, for the isolation and identification of gram (-) bacteria, gram (+) grains (Enterococcus, E. coli, Klebsiella, Proteus) (CHROMagar™ ORIENTATION).

PRINCIPLE OF THE METHOD

The Blood Agar Base contains high nutritional ingredients that provide vitamins, hydrocarbons and other organic elements. The Sodium chloride provides the necessary minerals and maintains the osmotic and electrolyte balance. Adding 5% of horse blood the nutrition of the material becomes stronger. Also, haemolysis and the type of haemolysis of red horse help in identifying certain gram (+) grains.

COMPOSITION BLOOD AGAR	g/litre
Columbia Peptone Mixture	25.1
Corn Starch	1.0
Sodium chloride	5.0
Agar No. 2	12.0
Horse Blood	50ml

Appearance: Red - crimson, not-clear due to the addition of blood.

Final pH 7.3 ± 0.2 in 25 °C.

The CHROMagar ORIENTATION is used to color the pathogens of the urogenital tract with E. Coli to give red colonies, Klebsiella metallic blue, Proteus mirabilis clear with brown hawthorn and Enterococcus faecalis blue-greenish blue with green halo.

This makes easy visual detection of crops with mixed flora.

COMPOSITION	g/litre
CHROMagar™ ORIENTATION	
Chromogenic mix	1,0
Peptone, yeast and extract	17,0
Agar	15,0

Appearance: light beige clear

Final pH 7.0 ± 0.2 in 25 °C.

PRECAUTIONS

BLOOD AGAR – CHROMagar™ ORIENTATION is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

Petri dishes must be stored in **2 – 12 °C** in their packaging until the time of their use. Prolonged storage at a temperature below **2 °C** creates enough moisture in the material at risk of contamination. Freezing even instantly destroys the material. Avoid overheating also. Petri dishes must be used until the expiration date which is shown on the label. When you open the airtight package, if you have any leftover plates, store them in the bag until the expiration date. For transport, our stability studies showed that the plates may remain at **6 - 25 °C for 4 days** or at **25 - 40 °C for 48 hours**, without affecting the attribution of the product.

USAGE

Place the dishes in the incubator (35-37 °C) for 30-45'.

Inoculate the specimen as soon as possible after it has been taken and overlay with successive dilutions for individual colonies.

Incubate in aerobic conditions at 35 - 37 °C for 24 hours. Consult the colony color chart.


LIMITATIONS OF THE METHOD

In CHROMagar™ ORIENTATION Some E. coli strains may show a pale pink color to colorless in their colonies, so we need to confirm with biochemical tests.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

BLOOD AGAR


Germ	ATCC	Growth	Haemolysis
<i>Escherichia coli</i>	25922	Good	-
<i>Pseudomonas aeruginosa</i>	27853	Good	-
<i>Staphylococcus aureus</i>	25923	Good	Beta
<i>Streptococcus pneumoniae</i>	6303	Good	ΑΙΠ
<i>Streptococcus pyogenes</i>	19615	Excellent	Beta

Streptococcus pyogenes ATCC 19615

CHROMagar™ ORIENTATION

Microorganism	ATCC	Characteristics of colonies
<i>Escherichia coli</i>	25922	Κόκκινες – κεραμιδί
<i>Enterococcus faecalis</i>	29212	blue-greenish with light blue halo
<i>Klebsiella pneumoniae</i>	13883	Metalic blue
<i>Proteus mirabilis</i>	12473	Colorless with brown halo
<i>Pseudomonas aeruginosa</i>	10145	Yellow with green material
<i>Salmonella typhimurium</i>	14028	Colorless
<i>Staphylococcus aureus</i>	25923	Light yellow
<i>Staphylococcus saprophyticus</i>	15305	Light pink
<i>Candida albicans</i>	60193	Creamy white

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

BLOOD AGAR – CHROMagar™ ORIENTATION - CE

ITEM	CODE	PACKAGE	KEEPING	LIFE SPAN
Petri dish halved 9cm	020211	10 pieces	2 – 12 °C	2 months

It is produced in Greece from the Bioprep company in accordance with the requirement of the European guidance 98/79/EK. Code by EDMA 14 01 04 90. The Bioprep company has been certified according to the standards EN ISO 9001:2008 / ΔY85/1348/2004.

BIBLIOGRAPHY

BLOOD AGAR

Ellner, P.D., Stoessel, C.J., Drakeford, E. and Vasi, F. (1966). A new culture medium for medical bacteriology. Amer J. Clin Pathol. 45. 502-504.

CHROMagar™ ORIENTATION

(1) Merlino, J. et al. 1996. Evaluation of CHROMagar Orientation for Differentiation and Presumptive Identification of Gram-Negative Bacilli and Enterococcus Species, J.C.M. 34: 1788-1793. (2) Samra, Z. et al. 1998. Evaluation of Use of a New Chromogenic Agar in Detection of Urinary Tract Pathogens. J.C.M. 36: 990-994.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr