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Ellner, P.D., Stoessel, C.J., Drakeford, E. and Vasi, F. (1966). A new culture medium for medical bacteriology. Amer J. Clin Pathol. 45. 502-504.
American Public Health Association (1950). Diagnostic Procedures and Reagents. 3rd edn. A.P.H.A., New York.

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6th 2024

DESCRIPTION

A medium recommended for the cultivation of a wide variety of pathogenic microorganisms. This medium is prepared with the addition of 5% sheep blood.

PRINCIPLE OF THE METHOD

The BLOOD AGAR (SHEEP BLOOD) contains components of high nutritional value, providing vitamins, carbohydrates, and other organic elements. Sodium chloride supplies the necessary minerals and maintains osmotic and electrolyte balance. The addition of 5% sheep blood enhances the medium's nutritional properties. Furthermore, the hemolysis and the type of hemolysis of sheep red blood cells assist in the identification of certain gram (+) cocci.

FORMULA	g/litre
Columbia peptone mixture	25.1
Soluble starch	1.0
Sodium chloride	5.0
Agar	12.0
Sheep Blood Defibrinated	50ml

Appearance: Red - maroon not clear, due to the addition of blood.

Final pH 7.1 ± 0.2 at 25 °C.

PRECAUTIONS

BLOOD AGAR (SHEEP BLOOD) is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

Petri dishes should be stored at **2–12°C** in their packaging until ready for use. Prolonged storage at temperatures below **2°C** can cause moisture to accumulate in the medium, increasing the risk of contamination. Freezing, even momentarily, destroys the medium. Avoid excessive heating as well. Petri dishes can be used until the expiration date indicated on the label. When you open the airtight package, if you have any leftover plates, store them in the bag until the expiration date. For transport, stability studies have shown that the dishes can remain at **18–25°C** for **4 days** or at **25–40°C** for **48 hours** without affecting product performance.

USAGE

Spread the petri dishes using the streaking technique to isolate single colonies.
Incubate them at 35–37°C for 24 hours under aerobic conditions.

INTERPRETATION OF RESULTS

Streptococcus pyogenes form small, gray, white colonies with β -hemolysis.

Streptococcus pneumoniae forms small, flat colonies with α -hemolysis.

Staphylococcus aureus forms relatively large colonies (1–3 mm after 24 hours of incubation and 3–8 mm if incubation is extended to 5 days).

Their color ranges from cream-yellow to orange depending on the incubation time. Finally, it causes α or β hemolysis.

LIMITATIONS OF THE METHOD

Final identification must be performed using biochemical and serological tests (e.g., Latex Agglutination Test), which can be conducted directly from the suspected colonies.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	Growth	Haemolysis
<i>Escherichia coli</i> ATCC 25922	Good	-
<i>Pseudomonas aeruginosa</i> ATCC 27853	Good	-
<i>Staphylococcus aureus</i> ATCC 25923	Good	b
<i>Streptococcus pneumoniae</i> ATCC 6303	Good	a
<i>Streptococcus pyogenes</i> ATCC 19615	Good	b

Streptococcus pyogenes ATCC 19615

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

BLOOD AGAR (SHEEP BLOOD) – GR/CA01/GRM5/O/2 - **CE**

PRODUCT	CODE	PACKING	STORAGE	SELF LIFE
Plate 9cm	010011	10 pieces	2 – 12 °C	2 months
Plate 6cm	050011	10 pieces	2 – 12 °C	2 months

Produced in Greece by Bioprep in accordance with the requirements of the European Directive 98/79/EK. FEK B2198/2-10-2009. CODE EDMA 14 01 04 90. Bioprep company has been certified according to standards EN ISO 9001:2008 / ΔΥ86/1348/2004.

BIBLIOGRAPHY

Ellener, P.C., C.J. Stoessel, E. Drakeford, and F. Vassi 1966. A new culture medium for medical bacteriology. Am J. Clin Pathol. 45:502-504. Ruoff, k., I. 1995 Streptococcus.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr