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Eley, A., Clarry, T., Bennett, K.W. (1989). Selective and differential medium for isolation of *Bacteroides ureolyticus* from clinical specimens. European Journal of Clinical Microbiology, Infectious Diseases. 8: 83-85.

Wade W. Griffiths, M. (1987). Comparison of Media for cultivation of subgingival bacteria. J. Dent. Res. 66: no. 4 abstract 334.

Heginbotham M., Fitzgerald T.C., and Wade W.G. (1990). Comparison of solid media for the culture of anaerobes. J. Clin. Path. 43: 253-256.

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TECHNICAL DATA SHEET

PRODUCT: **BRUCELLA BLOOD AGAR W/ Vit. K1 + Hemin**
REFERENCE: **010395**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

AEROMONAS BLOOD AGAR with ampicillin is recommended for the cultivation and selective isolation of Aeromonas spp.

PRINCIPLE OF THE METHOD

Tryptone Soy Agar serves as the nutrient base. Ampicillin inhibits the growth of most enterobacteria from normal fecal flora. However, a drawback is that it also inhibits the growth of Aeromonas trota, a pathogenic species that is sensitive to ampicillin. For more selective pathogen isolation, CIN Agar (cefsulodin-irgasan-novobiocin) is recommended as an alternative.

FORMULA	g/litre
Meat Peptone	10.0
Casein Peptone	10.0
Sodium chloride	5.0
Yeast Extract	2.0
Dextrose	1.0
Sodium Bisulfite	0.1
Agar No. 2	12.0
Hemin	0.005
Vitamin K1	0.01
Sheep Blood	50ml

Appearance: Red to dark red (burgundy), opaque due to the addition of blood.

Final pH 7.2 ± 0.2 at 25°C.

PRECAUTIONS

Brucella Blood Agar with Vitamin K1 & Hemin is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The petri dishes should be stored at **2 – 12°C** in their packaging until use. Prolonged storage at temperatures below **2°C** creates excess moisture within the material, posing a risk of contamination. Freezing, even momentarily, destroys the material. Excessive heating should also be avoided. The plates can be used until the expiration date indicated on the label. For transport, stability studies have shown that the plates can remain at **6 - 25°C** for **4 days** or at **25 - 40°C** for **48 hours** without affecting product performance.

INSTRUCTIONS FOR USE

Place the plates in the incubator (35 - 37°C) for 30 – 45 minutes. Inoculate the sample as soon as possible after collection and streak with successive dilutions to obtain isolated colonies. Incubate under anaerobic conditions at 35 - 37°C for 24 hours to five days. Additionally, inoculate on an aerobic blood medium as there may also be optional anaerobic bacteria present. Incubate the plate aerobically with 5-10% CO₂. Anaerobic conditions can be monitored using a specific anaerobic indicator, such as the Gas Pak anaerobic conditions indicator.

QUALITY CONTROL

Microorganism	ATCC	Growth	Colonies
<i>Bacteroides fragilis</i>	25285	Very Good	Gray without hemolysis
<i>Clostridium perfringens</i>	13124	Very Good	Gray, white with β-hemolysis
<i>Fusobacterium nucleatum</i>	25586	Very Good	Gray, white with dark gray zones
<i>Peptostreptococcus anaerobius</i>	27337	Very Good	Off-white colonies

Bacteroides fragilis ATCC 25285

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

BRUCELLA BLOOD AGAR W/ Vit. K1 + Hemin - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Petri dish 90mm	010395	10 pieces	2 – 8 °C	2 months

Produced in Greece by the company Bioprepare in accordance with the requirements of the European Directive 2017/746.

BASIC UDI-DI: 5212037714010401WF. EDMA (14 01 04 01) non-chromogenic media (Plates).

The Bioprepare company has been certified according to the standards: EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004

LITERATURE REFERENCES

Brazier, J.S. (1986). Yellow fluorescence of Fusobacteria Letters in Applied Microbiol. 2: 124-126.

Brazier, J.S. (1986). A note on ultra violet red fluorescence of anaerobic bacteria in vitro. J. Appl. Bact. 60: 121-126.

Eley, A., Clarry, T., Bennett, K.W. (1989). Selective and differential medium for isolation of *Bacteriodes ureolyticus* from clinical specimens. European Journal of Clinical Microbiology, Infectious Diseases. 8: 83-85.

Wade W. Griffiths, M. (1987). Comparison of Media for cultivation of subgingival bacteria. J. Dent. Res. 66: no. 4 abstract 334.

Heginbotham M., Fitzgerald T.C., and Wade W.G. (1990). Comparison of solid media for the culture of anaerobes. J. Clin. Path. 43: 253-256.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepare.gr