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www.bioprepare.gr

**TECHNICAL DATA SHEET**

PRODUCT: **CAMPYLOBACTER SELECTIVE
AGAR (SKIRROW)**

REFERENCE: **010224**

Date 1st Edition:

7th 2009

Date 4th Revision:

6th 2024

DESCRIPTION

This material is used for the selective isolation of *Campylobacter* spp from clinical specimens, food and environmental samples.

PRINCIPLE OF THE METHOD

The nutritional base is Columbia agar which is enriched with lousy red horses. Sodium thiosulphate, pyruvic acid and iron sulphate enhance microaerobes the *Campylobacter* spp. Peptones provide nitrogenous growth factors, carbon, sulfur and trace elements. Yeast extract is a rich source of vitamin B. The use of 4 antimicrobial agents inhibits the growth of rich microbial faecal flora resulting in isolation of the *C. jejuni*.

COMPOSITION	g/litre
Pancreatic Digest of Casein	10,0
Yeast Extract	5,0
Meat Peptic Digest	5,0
Sodium chloride	5.0
Heart Pancreatic Digest	3,0
Maize Starch	1,0
Agar No. 2	12.0
Lysed Horse Blood	60ml
Sodium Pyruvate	250mg
Sodium Metabisulphite	250mg
Ferrous Sulphate	250mg
Vancomycin	10mg
Polymyxin B	2500 U 0,3mg
Trimethoprim	5mg

Appearance: Red - crimson clear.

Final pH 7.3 ± 0.2 in 25 °C.

PRECAUTIONS

CAMPYLOBACTER SELECTIVE AGAR (SKIRROW) is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The agar plates should be stored at 2 – 12°C in their packaging until use. Prolonged storage below 2°C may create excessive moisture within the medium, increasing the risk of contamination. Freezing, even momentarily, will damage the medium. Excessive heating should also be avoided. The agar plates can be used until the expiration date indicated on the label. Once the airtight packaging is opened, any remaining plates should be stored in the provided storage bag until their expiration date. For transportation, stability studies have shown that the agar plates can remain at 6 - 25°C for up to 3 days or at 25 - 40°C for up to 24 hours without compromising product performance.

USAGE

Inoculate the sample (1/10ww) in *Campylobacter* Enrichment Broth (Bolton) with code 070404 (10ml) and 150404 (225ml). Incubate at 37°C ± 1°C for 4–6 hours and then at 41.5°C ± 1°C under microaerophilic conditions for 44 ± 4 hours. Subsequently, using a loop, transfer 10µl of Bolton broth to the edge of a Skirrow agar plate and perform serial dilutions with a sterile loop. Incubate at 41.5°C ± 1°C for 44 ± 4 hours under microaerophilic conditions (5–6% oxygen, 10% CO₂, and 84–85% nitrogen).

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Germ	ATCC	Development
<i>Campylobacter jejuni</i>	33291	Good. Small colonies, grey / do not show haemolysis.
<i>Escherichia coli</i>	25922	Partial suspension.

Campylobacter jejuni

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CAMPYLOBACTER SELECTIVE AGAR (SKIRROW)

GR/CA01/GRM5/O/73 - **CE**

SAMPLE	CODE	PACKAGE	KEEPING	LIFE SPAN
Petri dish 90mm	010636	10 pieces	2 – 12 °C	2 months

It is produced in Greece from the Bioprep company in accordance with the requirements of the European guidance 98/79/EK. FEK B2198/2-10-2009. Code by EDMA 14 01 04 01. The Bioprep company has been certified according to the standards EN ISO 9001:2015 / EΛOT EN ISO 13485:2016.ΔΥ85/1348/2004.

BIBLIOGRAPHY

Moyer, N. P., and L. A. Holcomb. 1995. *Brucella*, p. 549-555. In P. R. Murray, E. J. Tenover, M. A. Tenover, F. C. Tenover, and R. H. Tenover (ed.). Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.

MadFaddin, J. D. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. 1, p. 110-114. Williams & Wilkins, Baltimore, MD.

Isenberg, H. D. (ed.). 1992. Clinical microbiology procedures handbook. American Society for Microbiology, Washington, D.C.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr