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Lewis, J. S., and Wiesner, P. J. (1980) Gonorrhea: Current laboratory methods. Lab. Management, Sept., p.33.

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E-mail: bioprep1@otenet.gr www.bioprepare.gr

TECHNICAL DATA SHEET

PRODUCT: **CHOCOLATE AGAR – THAYER MARTIN AGAR**
REFERENCE: **020141**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

Divided CHOCOLATE AGAR – THAYER MARTIN AGAR Petri dish for the cultivation of microorganisms with special nutritional requirements, such as *Haemophilus* species and *Neisseria* species. (CHOCOLATE AGAR) – and for the isolation of *Neisseria gonorrhoeae* (THAYER MARTIN).

METHOD PRINCIPLE

Chocolate Agar is an enriched medium for qualitative procedures, used for the isolation and cultivation of pathogenic microorganisms, particularly *Neisseria* species and *Haemophilus* species, from a variety of clinical samples.

FORMULA	g/litre
CHOCOLATE AGAR	
Columbia Peptone Mixture	23.0
Corn Starch	1.0
Sodium chloride	5.0
Agar No. 2	12.0
Horse Blood	60ml

Appearance: Brown – chocolate, non-transparent, due to the addition of blood.

Final pH 7.3 ± 0.2 at 25 °C

Thayer Martin Agar is used for the isolation of pathogenic *Neisseria* from mixed flora samples (including bacteria and fungi). Its base contains nitrogenous nutrients in the form of casein and meat peptones. Corn starch neutralizes the toxic fatty acids that may be present in the agar. Horse blood, through its hemolysis, serves as an important nutrient, providing the X and V factors. Thayer Martin also contains antimicrobial agents, such as vancomycin, colistin, nystatin, and trimethoprim (V-C-N-T Inhibitor), which suppress the growth of normal flora. Vancomycin acts against gram-positive cocci, colistin inhibits gram-negative bacteria, including *Pseudomonas* species, but not *Proteus* species, which are inhibited by trimethoprim. Nystatin inhibits fungi.

FORMULA	g/litre
THAYER MARTIN AGAR	
Columbia Peptone Mixture	23.0
Corn Starch	1.0
Sodium chloride	5.0
Agar No. 2	12.0
Horse Blood	60ml
Vancomycin	3mg
Colistin	7.5mg
Nystatin	12.5mg
Trimethoprim	5mg

Appearance: Brown – chocolate, non-transparent, due to the addition of blood.

Final pH 7.2 ± 0.2 at 25 °C.

PRECAUTIONS

CHOCOLATE AGAR – THAYER MARTIN AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE AND TRANSPORT CONDITIONS

The petri dishes must be stored at **2 – 12 °C** inside their packaging until use. Prolonged storage at temperatures below **2 °C** creates excess moisture inside the material, posing a risk of contamination. Freezing, even temporarily, destroys the material. Excessive heating should also be avoided. The petri dishes can be used until the expiration date indicated on the label. For transport, stability studies have shown that the petri dishes can remain at **6 – 25 °C** for **5 days** or at **25 – 40 °C** for **48 hours** without affecting the product's performance.

METHOD OF USE

Inoculate and spread the sample as soon as possible after receiving it in the laboratory. Alternatively, if the sample is to be cultured after swabbing, the following procedure can be followed:

1. Pass the swab over the agar, drawing a large "Z" on its surface. This allows the swab to be exposed to the nutrient medium for enough time, ensuring better transfer of microorganisms.
2. Spread the material using a sterile loop, crossing over the initial "Z" direction. This procedure is best done during sample collection. If it was not done at that time, it should be performed in the laboratory.
3. Place cultures in an aerobic environment enriched with 10% CO₂ as soon as possible.

- Incubate at 35-37 °C and examine the cultures after overnight incubation as the initial phase, and again after approximately 48 hours for the final phase.
- Re-culture for identification of *N. gonorrhoeae* within the next 18-24 hours.

INTERPRETATION OF RESULTS

CHOCOLATE AGAR:

Haemophilus influenzae forms small, moist colonies with a slight odor.

Streptococcus pneumoniae forms small, flat colonies with a green color.

Neisseria meningitidis forms large, grayish-blue colonies.

Neisseria gonorrhoeae forms small, gray, white, colorless, mucous colonies.

THAYER MARTIN AGAR:

Neisseria gonorrhoeae forms small, gray, white, colorless, mucous colonies.

Neisseria meningitidis forms medium to large-sized gray colonies.

METHOD LIMITATIONS

Final identification should be done with biochemical and serological tests (e.g., Latex Agglutination Test), and can be performed directly from the suspected colonies.

QUALITY CONTROL

CHOCOLATE AGAR


Microorganism	ATCC	Colonies
<i>Haemophilus influenzae</i>	10211	Small, moist, with a slight odor
<i>Streptococcus pneumoniae</i>	6305	Small, flat colonies with a green color
<i>Neisseria meningitidis</i>	13090	Large colonies, grayish blue in color
<i>Neisseria gonorrhoeae</i>	43069	Small, gray, white, colorless, mucous

Haemophilus influenzae ATCC 10211

THAYER MARTIN

Microorganism	ATCC	Colonies
<i>Neisseria gonorrhoeae</i>	43069	Small, white - grayish
<i>Neisseria meningitidis</i>	13090	Medium to large-sized gray colonies
<i>Escherichia coli</i>	25922	Does not grow

Neisseria gonorrhoeae ATCC 43069

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHOCOLATE AGAR – THAYER MARTIN AGAR – CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Divided Petri dish 9 cm	020141	10 pieces	2 – 12 °C	3 months

Produced in Greece by the company Bioprep in accordance with the requirements of the European Directive 2017/746.

BASIC UDI-DI: 5212037714010401WF. EDMA (14 01 04 01) non-chromogenic media (Plates).

The Bioprep company has been certified according to the standards: EN ISO 9001:2015 / ΕΛΟΤ EN ISO 13485:2016 DY8d/1348/2004

LITERATURE REFERENCES

CHOCOLATE AGAR

Sng, E.H., Rajan, V.S., and Lim, A.L. (1977) Simplified media for isolating *Neisseria gonorrhoeae* J. Clin. Microbiol., 5 (4), 387.

Sottnek, F. O., Biddle, J. W., Kraus, S. J., Weaver, R. E., And Stewart, J. A. (1980) Isolation and Identification of *Haemophilus ducreyi* in a clinical study. J. Clin. Microbiol., 12 (2), 170.

Lewis, J. S., and Wiesner, P.J. (1980) Gonorrhea: Current laboratory methods. Lab. Management, Sept., p.33.

THAYER MARTIN AGAR

Young, H. 1978. Cultural diagnoses of gonorrhoea with modified New York City (MNYC) medium. Brit. Journ. Ven. Dis. 54: 36-40:

Thayer, J. D. and Martin, J. E. 1966. Improved medium selective for the cultivation of *N. gonorrhoeae* and *N. Meningitidis*: Public Health Rep. 81: 559-562.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr www.bioprep.gr