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E-mail: bioprep1@otenet.gr www.bioprepare.gr

**TECHNICAL DATA SHEET**

PRODUCT: **CHOCOLATE AGAR – THAYER MARTIN – GARDNERELLA AGAR**

REFERENCE: **030156**

Date 1st Edition:

7th 2009

Date 4th Revision:

6th 2024

DESCRIPTION

A tri-compartment petri dish containing CHOCOLATE AGAR – THAYER MARTIN – GARDNERELLA AGAR for the cultivation of bacteria with special nutritional requirements, such as *Haemophilus* and *Neisseria* (CHOCOLATE AGAR) – for the isolation of *Neisseria gonorrhoeae* (THAYER MARTIN) – for the selective isolation of *Gardnerella vaginalis* (GARDNERELLA AGAR).

PRINCIPLE OF THE METHOD

Chocolate Agar is an enriched medium used in qualitative procedures for the isolation and cultivation of pathogenic microorganisms, particularly *Neisseria* and *Haemophilus* species, from various clinical specimens.

FORMULA	g/litre
CHOCOLATE AGAR	
Columbia Peptone Mixture	23.0
Corn Starch	1.0
Sodium chloride	5.0
Agar No. 2	12.0
Horse Blood	60ml

Appearance: Brown-chocolate, non-transparent due to the addition of blood.

Final pH 7.3 ± 0.2 at 25 °C

Thayer Martin Agar is used for the isolation of pathogenic *Neisseria* from samples containing mixed flora (bacteria and fungi). Its base contains nitrogenous nutrients in the form of casein and meat peptones. Cornstarch neutralizes toxic fatty acids that may be present in the agar. Horse blood, through hemolysis, serves as an important nutrient, providing factors X and V. Thayer Martin also contains antimicrobial agents—vancomycin, colistin, nystatin, and trimethoprim (V-C-N-T Inhibitor)—which suppress the growth of normal flora. Vancomycin acts against Gram-positive cocci, colistin inhibits Gram-negative bacteria (including *Pseudomonas* species but not *Proteus* species, which are inhibited by trimethoprim), and nystatin suppresses fungi.

FORMULA	g/litre
THAYER MARTIN	
Columbia Peptone Mixture	23.0
Corn Starch	1.0
Sodium chloride	5.0
Agar No. 2	12.0
Horse Blood	60ml
Vancomycin	3mg
Colistin	7.5mg
Nystatin	12.5mg
Trimethoprim	5mg

Appearance: Brown-chocolate, non-transparent due to the addition of blood.

Final pH 7.2 ± 0.2 at 25 °C.

GARDNERELLA AGAR is a selective medium used for the isolation of *Gardnerella vaginalis*. It is incubated in 10% CO₂. The Columbia base is a nutrient-rich formula that provides the necessary nutrients for the growth of fastidious bacteria. The enzymes from casein, meat peptones, and beef extract provide essential amino acids and other nitrogenous complexes. Yeast extract primarily provides vitamin B. Cornstarch helps neutralize fatty acids, which could be toxic to *G. vaginalis*. NaCl maintains osmotic balance. Colistin, nalidixic acid, and amphotericin B are added to facilitate the selective growth of *G. vaginalis* from clinical samples. Colistin and nalidixic acid inhibit most Gram-negative bacteria, while amphotericin B inhibits most fungi. The addition of horse blood enhances the nutritional quality of the medium and aids in the identification of *G. vaginalis* due to β -hemolysis on red blood cells.

FORMULA	g/litre
GARDNERELLA AGAR	
Peptone Mixture	20.0
Beef Extract	3.0
Corn Starch	1.0
Sodium chloride	5.0
Bacteriological Agar	15.0
Horse Blood	100ml
Colistin	10mg
Nalidixic acid	15mg
Amphotericin B	3mg

Appearance: Red – dark cherry, non-transparent, due to the addition of blood.

Final pH 7.3 ± 0.2 at 25 °C.

PRECAUTIONS

CHOCOLATE AGAR – THAYER MARTIN – GARDNERELLA AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE AND TRANSPORT CONDITIONS

The petri dishes must be stored at **2 – 12 °C** inside their packaging until use. Prolonged storage at temperatures below **2 °C** creates excess moisture inside the material, posing a risk of contamination. Freezing, even temporarily, destroys the material. Excessive heating should also be avoided. The petri dishes can be used until the expiration date indicated on the label. For transport, stability studies have shown that the petri dishes can remain at **6 – 25 °C** for **4 days** or at **25 – 40 °C** for **48 hours** without affecting the product's performance.

METHOD OF USE

Inoculate and spread the sample as soon as possible after receiving it in the laboratory. Alternatively, if the sample is to be cultured after swabbing, the following procedure can be followed:

1. Pass the swab over the agar, drawing a large "Z" on its surface. This allows the swab to be exposed to the nutrient medium for enough time, ensuring better transfer of microorganisms.
2. Spread the material using a sterile loop, crossing over the initial "Z". This procedure is best done when the sample is collected. If not, it should be performed in the laboratory.
3. Place cultures in an aerobic environment enriched with 10% CO₂ as soon as possible.
4. Incubate at 35-37 °C and examine the cultures after overnight incubation as an initial phase, and again after approximately 48 hours for the final phase.
5. Re-culture for identification of *N. gonorrhoeae* within the next 18-24 hours.

INTERPRETATION OF RESULTS

CHOCOLATE AGAR:

Haemophilus influenzae forms small, moist colonies with a slight odor.

Streptococcus pneumoniae forms small, flat colonies with a green color.

Neisseria meningitidis forms large, grayish-blue colonies.

Neisseria gonorrhoeae forms small, gray, white, colorless, mucous colonies.

THAYER MARTIN:

Neisseria gonorrhoeae forms small, gray, white, colorless, mucous colonies.

Neisseria meningitidis forms medium to large-sized gray colonies.

GARDNERELLA AGAR:

Gardnerella vaginalis forms small, gray, white colonies with beta-hemolysis.

METHOD LIMITATIONS

Final identification must be carried out with biochemical and serological tests (e.g., Latex Agglutination Test), and can be performed directly from the suspected colonies.

QUALITY CONTROL

CHOCOLATE AGAR


Microorganism	ATCC	Colonies
<i>Haemophilus influenzae</i>	10211	Small, moist, with a slight odor.
<i>Streptococcus pneumoniae</i>	6305	Small, flat, with a green color.
<i>Neisseria meningitidis</i>	13090	Large colonies, grayish blue in color.
<i>Neisseria gonorrhoeae</i>	43069	Small, gray, white, colorless, mucous.

Haemophilus influenzae ATCC 10211

THAYER MARTIN

Microorganism	ATCC	Colonies
<i>Neisseria gonorrhoeae</i>	43069	Small white - grayish.
<i>Neisseria meningitidis</i>	13090	Medium to large-sized gray colonies
<i>Escherichia coli</i>	25922	Does not grow.

Neisseria gonorrhoeae ATCC 43069

GARDNERELLA AGAR

Microorganism	ATCC	Growth	Hemolysis
<i>Gardnerella vaginalis</i>	14018	Excellent	B
<i>Candida albicans</i>	10231	Inhibited	-
<i>Proteus mirabilis</i>	12453	Inhibited	-

Gardnerella vaginalis ATCC 14018

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHOCOLATE AGAR – THAYER MARTIN – GARDNERELLA AGAR - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Petri dish, 9 cm divided	030156	10 pieces	2 – 12 °C	2 months

Produced in Greece by the company Bioprepere in accordance with the requirements of the European Directive 2017/746.

BASIC UDI-DI: 5212037714010401WF. EDMA (14 01 04 01) non-chromogenic media (Plates).

The Bioprepere company has been certified according to the standards: EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004

LITERATURE REFERENCES

CHOCOLATE AGAR

Sng, E.H., Rajan, V.S., and Lim, A.L. (1977) Simplified media for isolating *Neisseria gonorrhoeae* J. Clin. Microbiol., 5 (4), 387.

Sottnek, F. O., Biddle, J. W., Kraus, S. J., Weaver, R. E., And Stewart, J. A. (1980) Isolation and Identification of *Haemophilus ducreyi* in a clinical study. J. Clin. Microbiol., 12 (2), 170.

Lewis, J. S., and Wiesner, P.J. (1980) Gonorrhea: Current laboratory methods. Lab. Management, Sept., p.33.

THAYER MARTIN

Young, H. 1978. Cultural diagnoses of gonorrhoea with modified New York City (MNYC) medium. Brit. Journ. Ven. Dis. 54: 36-40:

Thayer, J. D. and Martin, J. E. 1966. Improved medium selective for the cultivation of *N. gonorrhoeae* and *N. Meningitidis*: Public Health rep. 81: 559-562.

GARDNERELLA AGAR

Ellner, P.D., Stoessel, C.J., Drakeford, E and Vasi, F. (1966). A new culture medium for medical bacteriology. Amer. J. Clin Pathol., 45:502-504.

Goldberg, R.L., and Washington, J.A., (1976). Comparison of isolation of *Haemophilus vaginalis* (*Corynebacterium vaginale*) from Peptone-Starch-Dextrose Agar and Columbia Colistin-Nalidixic Acid Agar. J. Clin. Microbiol., 4:245-247.

Thayer, D.D. and Martin, H. E. (1966). An improved medium for the cultivation of *N. gonorrhoeae* and *N. meningitidis*. Publ. Hlth. Report, 81:559-562.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepare.gr