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CHOCOLATE AGAR W/ POLIENRICHMENT SUPPLEMENT - CE

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E-mail: bioprep1@otenet.gr www.bioprep.gr

**TECHNICAL DATA SHEET**

PRODUCT: **CHOCOLATE AGAR**
W/ POLIENRICHMENT SUPPLEMENT
 REFERENCE: **010026 – 060026**

Date 1st Edition:
 7th 2009
 Date 4th Revision:
 6th 2024

DESCRIPTION

The CHOCOLATE AGAR W/ POLIENRICHMENT SUPPLEMENT is used for the cultivation of microaerophilic bacteria (*Haemophilus*, *Neisseria*).

PRINCIPLE OF THE METHOD

The Chocolate Agar is an enriched substrate for quality processes, for the isolation and cultivation of pathogenic microorganisms, particularly *Neisseria* and *Haemophilus* species, from a variety of clinical specimens. Adding GC Growth Supplement strengthens factors X and V.

CONFIGURATION	g/litre
Columbia Peptone Mixture	25.1
Soluble starch	1.0
Sodium chloride	5.0
Agar	12.0
Hemoglobin	10
GC Growth Supplement	5ml

Εμφάνιση: Brown – chocolate unclear.

Final pH 7.3 ± 0.2.

PRECAUTIONS

CHOCOLATE AGAR W / BVA is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE AND TRANSPORTS CONDITIONS

The agar plates should be stored at **2 – 12°C** in their packaging until use. Prolonged storage at temperatures below **2°C** creates excessive moisture inside the material, increasing the risk of contamination. Even momentary freezing damages the material. Excessive heating should also be avoided. The agar plates can be used until the expiration date stated on the label. If the airtight packaging is accidentally opened, the plate can be stored in the refrigerator for **5 – 7 days**, sealed with parafilm or placed in a storage bag. Stability studies have shown that for transport, the plates can remain at **6 – 25°C** for up to **4 days** or at **25 – 40°C** for **48 hours** without affecting product performance.

USAGE

Ready-to-use agar plate: Spread the sample on the edge of the medium and perform serial dilutions using a loop. The sample may include vaginal or cervical swabs, wound exudates, bronchial secretions, etc. Incubate the plates at 36 ± 1 °C for 24-48 hours, and if necessary, up to 72 hours under microaerophilic conditions with 8-10% CO₂ and 70% humidity.

100ml vial: Remove the vial from the refrigerator one hour before use to allow it to reach room temperature (20-25 °C). Heat the water bath until the water begins to boil gently. Place the vial in boiling water for 20-25 minutes, after unsealing and slightly loosening the cap. Remove the vial from the water bath and let it cool for 30-40 minutes at room temperature (20-25 °C), so the material reaches a temperature of 54-50 °C. Dispense the medium into agar plates: 90mm plates with 20ml and 60mm plates with 10ml. Spread the sample on the edge of the medium and perform serial dilutions using a loop. The sample may include vaginal or cervical swabs, wound exudates, bronchial secretions, etc. Incubate the plates at 36 ± 1 °C for 24-48 hours, and if necessary, up to 72 hours under microaerophilic conditions with 8-10% CO₂ and 70% humidity.

Note: Do not heat the liquefied vial for more than three hours at 45-50 °C. The vial can only be liquefied once.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	ATCC	Colonies
<i>Haemophilus influenzae</i>	10211	Small, wet, with a light smell.
<i>Streptococcus pneumoniae</i>	6305	Small, flat, appears green color.
<i>Neisseria meningitidis</i>	13090	Large colonies, gray, blue color.
<i>Neisseria gonorrhoeae</i>	43069	Small, gray, white, colorless, mucous.

Haemophilus influenzae

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHOCOLATE AGAR W/ POLIENRICHMENT SUPPLEMENT - **CE**

ITEM	CODE	PACKAGE	KEEPING	LIFE SPAN
Petri dish 90mm	010026	10 pieces	2 – 12 °C	3 months
Tube 100ml	060026	10 pieces	6 – 12 °C	8 months

Produced in Greece by the company Bioprepare in accordance with the requirements of the European Directive 2017/746.

BASIC UDI-DI: 5212037714010401WF. EDMA (14 01 04 01) non-chromogenic media (Plates).

The Bioprepare company has been certified according to the standards: EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004

BIBLIOGRAPHY

Sng, E. H., Rajan, V. S., and Lim, A. L. (1977) Simplified media for isolating *Neisseria gonorrhoeae*. J. Clin. Microbiol., 5 (4), 387.

Sottnek, F. O., Biddle, J. W., Kraus, S. J., Weaver, R. E., And Stewart, J. A. (1980) Isolation and Identification of *Haemophilus ducreyi* in a clinical study. J. Clin. Microbiol., 12 (2), 170.

Lewis, J. S., and Wiesner, P. J. (1980) Gonorrhea: Current laboratory methods. Lab. Management, Sept., p.33.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepare.gr