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Detection of *Acinetobacter baumannii* in Surveillance Cultures

2011. J. Kristie Johnson et al. University of Maryland Poster ASM 2010.

Modifications to CHROMagar *Acinetobacter* for improved selective growth of multi-drug-resistant *Acinetobacter baumannii*. 2011.

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TECHNICAL DATA SHEET

PRODUCT: **CHROMagar™ Acinetobacter**
REFERENCE: **010455 – 050455**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

CHROMagar™ Acinetobacter is used for the isolation and differentiation of Acinetobacter and other significant multidrug-resistant (MDR) bacteria.

INTRODUCTION

Acinetobacter is an organism with a high survival capacity on environmental surfaces. Its ability to acquire antimicrobial resistance is a major concern, particularly for hospital-acquired infections. In hospitals, Acinetobacter baumannii, for example, can enter the body through open wounds, catheters, and respiratory tubes. Any effective infection control policy should include stool testing as part of its protocol. CHROMagar™ Acinetobacter is specifically designed to facilitate this process, allowing the development of colonies with a distinctive bright red color.

PRINCIPLE OF THE METHOD

Meat peptones provide nitrogen, vitamins, minerals, and amino acids essential for growth. Yeast extract is the source of vitamins, particularly those of the B group. Sodium chloride provides essential electrolytes for transport and regulates the osmotic balance of the medium. The liquid supplement gives the medium selectivity and enhancement for growth. The chromogenic mix supplies chromogenic substances for amino acid binding. Bacteriological agar acts as the solidifying agent.

FORMULA	g/litre
Chromogenic mix	1,8
Peptone, yeast and extract	12,0
Agar	15,0
Salts	4,0
Liquid supplement	1ml

Appearance: Light pink clear

Final pH 7.0 ± 0.2 in 25 °C.

PRECAUTIONS

CHROMagar™ Acinetobacter is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE AND TRANSPORTS CONDITIONS

Petri dishes should be stored at **2–12 °C** in their packaging until use. Prolonged storage at temperatures below **2 °C** may cause significant moisture buildup in the medium, increasing the risk of contamination. Freezing, even briefly, will damage the medium. Avoid excessive heating as well. Petri dishes can be used until the expiration date indicated on the label. When you open the airtight package, if you have any leftover plates, store them in the bag until the expiration date. Stability studies have shown that for transportation, the dishes can remain at **6–25 °C** for up to **3 days** or at **25–40 °C** for up to **24 hours** without affecting product performance.

USAGE

Infectious material must be transported to the laboratory promptly without delay and should be protected from excessive heat and cold. If processing is delayed, the sample should be inoculated into an appropriate transport medium and maintained at 2–12 °C until further handling. Allow the Petri dishes to warm to room temperature before use. Ensure that the agar surface is dry before inoculation. Inoculate the sample at one edge of the dish and then perform successive streaks with a loop in parallel lines to create isolated colonies. Incubate the dishes aerobically at 35–37 °C for 20 to 28 hours.

INTERPRETATION OF RESULTS


Acinetobacter baumannii forms colonies with a bright red color, making it easy to distinguish from other bacteria. Other Gram-negative bacteria grow slowly and form blue colonies.

LIMITATIONS OF THE METHOD

Final identification may require additional tests, such as biochemical or immunological tests (e.g., latex agglutination).

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	Growth – Colony Colour
<i>Acinetobacter baumannii</i> ATCC 19906	Red
<i>Acinetobacter calcoaceticus</i> ATCC 19906	Red
Other Gram (-) Bacteria	Blue, violet colorless or inhibited
Gram (+) cocci and fungi	Mostly Inhibited

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

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Produced in Greece by the company Bioprep in accordance with the requirements of the European Directive 2017/746.

BASIC UDI-DI: 5212037714010401WF. EDMA (14 01 04 01) non-chromogenic media (Plates).

The Bioprep company has been certified according to the standards: EN ISO 9001:2015 / ΕΛΟΤ EN ISO 13485:2016 DY8d/1348/2004.

BIBLIOGRAPHY

Laboratory evaluation of different agar media for isolation of carbapenem-resistant *Acinetobacter* spp.. 2014. J.Moran-Gilad; A.Adler. D.Schwartz; S.Navon-Venezia; Y.Carmeli Eur J Clin Microbiol Infect Dis.

Detection of *Acinetobacter baumannii* in Surveillance Cultures

2011. J. Kristie Johnson et al. University of Maryland Poster ASM 2010.

Modifications to CHROMagar Acinetobacter for improved selective growth of multi-drug resistant *Acinetobacter baumannii*. 2011.

D W Wareham, N.C. Gordon JCP Online First 10.1136/jcp.2010.0834692010

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr