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E-mail: bioprep1@otenet.gr www.bioprepare.gr

DESCRIPTION

CHROMagar™ C. difficile is a chromogenic medium used for the color differentiation of *Clostridium difficile*.

INTRODUCTION

Clostridium difficile is the leading cause of hospital-associated infectious diarrhea in adults. These infections primarily occur in patients receiving medical care and antibiotic therapy. Symptoms of *C. difficile* infection include fever, abdominal cramps, and severe diarrhea, which can lead to death. In the United States, nearly 250,000 people are affected by *C. difficile* with at least 14,000 deaths (CDC estimate, 2013). Due to the emergence of the highly toxigenic strain of *C. difficile*, these infections have become more common and difficult to treat in recent years. While PCR has become the leading detection technique for *C. difficile*, culturing remains necessary for strain typing and antimicrobial resistance testing. CHROMagar™ C. difficile is a novel chromogenic medium, highly sensitive and selective, specifically designed to simplify and accelerate (24 hours) the cultivation of *C. difficile*.

METHOD PRINCIPLE

Peptones and yeast extract provide the medium with essential nutrients for bacterial growth. They are the primary source of organic nitrogen, particularly amino acids and polypeptides, and are also rich in vitamins, minerals, and trace elements. The chromogenic mix provides the medium with chromogenic substances for amino acid binding. Growth factors provide the necessary elements for *C. difficile* growth. The salts regulate the osmotic balance of the medium. The agar solidifies the medium.

FORMULA	g/litre
Peptones and yeast extract	25.0
Chromogenic mix	1.7
Salts	9.0
Growth factors	4.0
Agar	15.0

Appearance: Light beige, clear

Final pH 7.8 ± 0.2 at 25 °C.

PRECAUTIONS

CHROMagar™ C. difficile is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE AND TRANSPORT CONDITIONS


The Petri dishes must be stored at **2 – 12°C** in their packaging until use. Prolonged storage at temperatures below **2°C** may cause moisture buildup inside the material, posing a risk of contamination. The Petri dishes can be used until the expiration date indicated on the label. For transport, stability studies have shown that the tubes can remain at **6 – 25°C** for **4 days** or at **25 – 40°C** for **48 hours** without affecting product performance.

HOW TO USE

Inoculate the sample onto the surface of the material with serial dilutions. Incubate for 24 – 48 hours at 35 – 37°C under anaerobic conditions. Observe the colonies and consult the table.

READING AND INTERPRETATION

After 24 hours of incubation in an anaerobic atmosphere, *C. difficile* colonies appear large (approximately 2mm), colorless colonies that fluoresce under ultraviolet light at 365 nm. (Traditional media requires 48 hours). Other anaerobic bacteria that may grow display colorless colonies without fluorescence.

LIMITATIONS OF THE METHOD

Final identification should be performed with biochemical and serological tests (e.g., Microgen C. Difficile Latex Test, code: M41CE), which can be conducted directly from the suspected colorless *Clostridium* colonies that fluoresce under ultraviolet light. An oxidase test (Oxidase Test Strips-Code MID-61G) can also be performed directly on the suspected colony. Some Gram-negative bacteria resistant to the material's inhibitors may grow and produce false-positive colonies.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	Growth - Colony Color
<i>Clostridium difficile</i> ATCC 43255	Good - Colorless colonies that fluoresce under ultraviolet light
<i>Clostridium perfringens</i> ATCC® 13124	Inhibited
<i>Escherichia coli</i> ATCC 25922	Inhibited
<i>Candida albicans</i> ATCC® 10231	Inhibited

Other Anaerobes *Clostridium difficile*

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHROMagar™ C. difficile - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
9cm Plate, 20ml	010673	10 pieces	2 – 12 °C	2 months
6cm Plate, 9ml	050673	10 pieces	2 – 12 °C	2 months

Produced in Greece by Bioprep Company in accordance with Regulation (EU) 2017/746.

Basic UDI-DI: 5212037714010402WH. EDMA: (14 01 04 02) Chromogenic ID Media (Plates).

Bioprep Company is certified according to the following standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016, DY86/1348/2004.

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IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepare.gr