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TECHNICAL DATA SHEET

PRODUCT: **CHROMagar™ ENTEROCOCCUS**
REFERENCE: **010361 – 050361**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

CHROMagar™ ENTEROCOCCUS is used for the chromatic differentiation and isolation of *Enterococcus faecalis* and *Enterococcus faecium*.

PRINCIPLE OF THE METHOD

Meat peptones provide nitrogen, vitamins, minerals, and amino acids necessary for growth. Yeast extract is the source of vitamins, especially the B-group vitamins. Sodium chloride provides essential electrolytes for transport and regulates the osmotic balance of the medium. Bacteriological agar is the solidifying agent. In CHROMagar ENTEROCOCCUS, *Enterococcus faecalis* and *Enterococcus faecium* grow within 18-24 hours and are the only colonies that appear pink. The light-colored medium helps in easily identifying the pink colonies amidst the blue or colorless colonies of other Gram (+) cocci. Gram (-) bacteria do not grow.

FORMULA	g/litre
Chromogenic mix	27.3
Peptones and yeast extract	20.0
Sodium chloride	5,0
Agar	15,0

Appearance: White, non-transparent.

Final pH 6.8 ± 0.2 at 25 °C.

PRECAUTIONS

CHROMagar™ ENTEROCOCCUS is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

TORAGE AND TRANSPORT CONDITIONS

The plates should be stored at **2–12°C** in their packaging until use. Prolonged storage at temperatures below **2°C** causes moisture to accumulate within the material, which may lead to contamination. Freezing, even briefly, destroys the material. Excessive heating should also be avoided. The plates can be used up until the expiration date indicated on the label. For transportation, stability studies have shown that the plates can be stored at **17–25°C** for **4 days** or at **27–40°C** for **48 hours** without affecting the product's performance.

USAGE INSTRUCTIONS

The sample material should reach the laboratory quickly without delay and should be protected from excessive heat and cold. If there will be a delay in processing, the sample should be inoculated in an appropriate transport medium and kept at 2–12°C until inoculation. Allow the plates to warm to room temperature. The agar surface should be dry before inoculation. Inoculate the sample at one edge of the plate and then streak parallel lines with the loop to create isolated colonies. Incubate the plates aerobically at 35–37°C for 20 to 28 hours.

INTERPRETATION OF RESULTS

Enterococcus faecalis & *Enterococcus faecium* form colonies with a pink to purple color.

Streptococcus spp. form colonies with a blue color.


Gram-negative bacteria do not grow.

METHOD LIMITATIONS

Final identification should be carried out using biochemical and serological tests (e.g., Microgen Strep latex test, code: M47CE) and can be performed directly from the suspected *Enterococcus* colonies.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	Colony Development / Color
<i>Enterococcus faecalis</i> ATCC® 29212	Έντονο ροζ
<i>Streptococcus pyogenes</i> ATCC® 19615	Μπλε
<i>Escherichia coli</i> ATCC® 25922	Αναστέλλεται

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHROMagar™ ENTEROCOCCUS

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Petri Dish 9cm	010361	10 pieces	2 – 12 °C	3 months
Petri Dish 6cm	050361	10 pieces	2 – 12 °C	3 months

Produced in Greece by Bioprep Company in accordance with Regulation (EU) 2017/746.

Basic UDI-DI: 5212037714010402WH. EDMA: (14 01 04 02) Chromogenic ID Media (Plates).

Bioprep Company is certified according to the following standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016, DY86/1348/2004.

REFERENCES

Evaluation of CHROMagar compared with enterococcosel broth for the isolation of vancomycin resistant enterococci. 2012.

Hindley et al Department of microbiology and infectious diseases ASM 2012.

Evaluation of Three Chromogenic Media for Detection of Vancomycin-Resistant Enterococci in a tertiary-care Hospital M.L. Miller et al. 2011. Kingston General Hospital, ON, Canada. Poster P26, CACMID 2011.

Evaluation of Three Commercial Chromogenic Media and BEAA + van 6ug/mL for the Detection of Vancomycin-Resistant Enterococcus (VRE). 2011. Kornherr, Department of Microbiology Gamma Dynacare Medical Laboratories, Ottawa and Toronto, Ontario, Canada. ASM Meeting Poster 2010.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr