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Evaluation of three selective chromogenic media, CHROMagar ESB, CHROMagar CTX-M and CHROMagar KPC, for the detection of *Klebsiella pneumoniae* producing OXA-48 carbapenemase

2013. Michael Hornsey, Lynette Phee, Neil Woodford, Jane Turton, Daniele Meunier, Claire Thomas, David W Wareham.

Comparison of three selective media for the recovery of Extended Spectrum β -Lactamase (ESBL)-producing Enterobacteriaceae

2011. M. Jones, A. Sweeney, E. Stoeppler, M. Miller, and P. Gilligan Clinical Microbiology-Immunology Laboratories UNC Hospitals 2011

Evaluation of a chromogenic medium for extended-spectrum beta-lactamase (ESBL)-producing Enterobacteriaceae

2010. Philippe Lagacé-Wiens et al. University of Manitoba, Canada. ECCMID Poster 2010

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TECHNICAL DATA SHEET

PRODUCT: **CHROMagar™ ESBL**
REFERENCE: **010413– 050413**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

CHROMagar™ ESBL is used for the detection of Gram-negative bacteria with ESBL resistance (production of ESBL beta-lactamases).

INTRODUCTION

ESBLs are enzymes (beta-lactamases) produced by Gram-negative bacteria that inactivate beta-lactam antibiotics (e.g., ceftazidime, cefotaxime, and ceftriaxone). The presence of an organism that produces ESBL enzymes in a clinical infection can lead to treatment failure if one of these drugs is used. ESBL-producing bacteria are difficult to detect. CHROMagar™ ESBL enables the detection of bacteria that produce ESBL enzymes while inhibiting the growth of other bacteria, including those that produce AmpC beta-lactamases. This is an important feature because AmpC resistance has no clinical relevance but often leads to false-positive readings for ESBL bacteria in traditional diagnostic methods.

PRINCIPLE OF THE METHOD

The peptones in the medium provide the essential nutrients. Selective agents inhibit the growth of Gram-negative bacteria, yeasts, and some Gram-positive cocci. The chromogenic substrates are broken down by specific microbial enzymes, releasing insoluble colored compounds.

FORMULA	g/litre
Chromogenic mix	1,0
Peptones and yeast extract	17.0
Agar	15,0
ESBL supplement	570mg

Appearance: light beige clear

Final pH 7.0 ± 0.2 in 25 °C.

PRECAUTIONS

CHROMagar™ ESBL is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIOS

The plates should be stored at **2–12 °C** in their original packaging until ready for use. Prolonged storage at temperatures below **2 °C** can cause moisture buildup inside the material, risking contamination. Freezing, even briefly, destroys the material. Also, excessive heating should be avoided. The plates can be used until the expiration date indicated on the label. When you open the airtight package, if you have any leftover plates, store them in the bag until the expiration date. For transportation, stability studies have shown that the plates can remain at **6–25 °C** for **3 days** or at **25–40 °C** for **24 hours** without affecting the performance of the product.

USAGE

The infectious material must reach the laboratory quickly, without delay, and be protected from excessive heat and cold. If there is a delay in processing, the sample should be inoculated into an appropriate transport medium and stored at 2-12°C until inoculation. Allow the plates to warm to room temperature. The surface of the agar should be dry before inoculation. Inoculate the sample at one edge of the plate and then make successive streaks with the loop in parallel lines to isolate individual colonies. Incubate the plates aerobically at 35-37°C for 18 to 24 hours.

INTERPRETATION OF RESULTS

ESBL *E. coli* forms colonies with a dark pink to reddish color and is easily distinguishable from other bacteria.

ESBL KEC (*Klebsiella*, *Enterobacter*, *Citrobacter*) forms colonies with a metallic blue color.

Proteus BLSE forms colonies that are beige, brown with a light brown color in the medium.

ESBL *Acinetobacter* forms colonies with a cream color.





ESBL *Pseudomonas* forms colonies with a clear, natural cream color and a green hue in the medium.

LIMITATIONS OF THE METHOD

Some *Pseudomonas* spp and *Acinetobacter* spp, widely known as antibiotic-resistant bacteria, could grow with typical colonies as described. Final identification may require additional tests such as biochemical or immunological testing (latex agglutination). The confirmation test can be performed directly from the suspected colonies. Most bacteria producing AmpC do not grow, but some may grow poorly.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	Growth - Colony Color
<i>Pseudomonas aeruginosa</i> ATCC® 27853	Inhibit
<i>Escherichia coli</i> ATCC® 25922	Inhibit
<i>Candida albicans</i> ATCC® 10231	Inhibit
<i>Staphylococcus aureus</i> ATCC® 25923	Inhibit

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHROMagar™ ESBL - CE

ITEM	CODE	PACKAGE	KEEPING	LIFE SPAN
Τρυβλίο 9cm	010413	10 pieces	2 – 12 °C	2 months
Τρυβλίο 6cm	050413	10 pieces	2 – 12 °C	2 months

Produced in Greece by Bioprepare in accordance with the requirements of the European Directive 98/79/EK. FEK B2198/2-10-2009. CODE EDMA 14 01 04 02. Bioprepare company has been certified according to standards: EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004.

BIBLIOGRAPHY

Evaluation of three selective chromogenic media, CHROMagar ESBL, CHROMagar CTX-M and CHROMagar KPC, for the detection of *Klebsiella pneumoniae* producing OXA-48 carbapenemase

2013. Michael Hornsey, Lynette Phee, Neil Woodford, Jane Turton, Daniele Meunier, Claire Thomas, David W Wareham.

Comparison of three selective media for the recovery of Extended Spectrum β -Lactamase (ESBL)-producing Enterobacteriaceae

2011. M. Jones, A. Sweeney, E. Stoeppler, M. Miller, and P. Gilligan Clinical Microbiology-Immunology Laboratories UNC Hospitals 2011

Evaluation of a chromogenic medium for extended-spectrum beta-lactamase (ESBL)-producing Enterobacteriaceae

2010. Philippe Lagacé-Wiens et al. University of Manitoba, Canada. ECCMID Poster 201

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepare.gr