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TECHNICAL DATA SHEET

PRODUCT: **CHROMagar™ mSuperCARBA™**
 REFERENCE: **010710 - 050710**

Date 1st Edition:
 7th 2009
 Date 4th Revision:
 6th 2024

DESCRIPTION

CHROMagar™ mSuperCARBA™ is a selective and differential chromogenic culture medium, intended for use in the qualitative direct detection of gastrointestinal colonization with carbapenem-resistant Enterobacteria (CRE), including OXA-48 producers, to aid in the prevention and control of CRE in healthcare settings. The test is performed with rectal swab and stools from patients to screen for CRE colonization. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C.

INTRODUCTION

Carbapenem-resistant Enterobacteriaceae (CRE) are typically resistant to all β -lactams and most other classes of antimicrobial agents. Treatment options for patients infected with CRE are very limited. Healthcare-associated CRE outbreaks have been reported. Patients colonized with CRE are believed to serve as a source of transmission in healthcare settings. Identifying patients who are colonized with CRE and placing them in isolation with precautions can be a crucial step in preventing transmission. Alain Rambach and Patrice Nordmann have combined their efforts to develop a highly sensitive chromogenic medium, CHROMagar™ mSuperCARBA™, which belongs to a new generation of chromogenic materials achieving unprecedented performance: the detection of a wide range of carbapenems (KPC, NDM, VIM, IMP, OXA) with an impressive detection limit (10 CFU/ml), even for weakly expressed carbapenems like OXA-48, while maintaining a high level of selectivity.

PRINCIPLE OF THE METHOD

The peptones provide the necessary nutrients in the medium. A mixture of antibiotics allows for the selective growth of KPC, NDM, VIM, IMP, and OXA with a detection limit of 10 CFU/ml. Three chromogenic substrates are broken down by the microbial enzymes of the bacteria, releasing insoluble colored compounds that adhere to the microbial cell membrane during their multiplication, forming colonies with color. *Escherichia coli* produces the enzymes β -glucuronidase and β -galactosidase, and its colonies develop a pink to crimson color. *Klebsiella*, *Enterobacter*, *Serratia*, and *Citrobacter* produce β -glucuronidase, and their colonies turn blue green. If none of the substrates are used, there will be colorless or white colonies. The salts maintain the osmotic balance in the medium. The agar provides trace elements and electrolytes and solidifies the medium.

| FORMULA | g/litre |
|-------------------------------|---------|
| Chromogenic and selective mix | 0,8 |
| Peptones | 20,0 |
| Salts | 5,0 |
| Agar | 15,0 |
| Growth factors | 1,7 |

Appearance: Light beige, transparent.

Final pH 7.2 ± 0.2 at 25 °C.

PRECAUTIONS

CHROMagar™ mSuperCARBA™ is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The plates should be stored at **2 – 12 °C** inside their packaging until they are used. Prolonged storage at temperatures below **2 °C** can cause excessive moisture buildup in the material, leading to a risk of contamination. Freezing, even momentarily, destroys the material. Excessive heating should also be avoided. The plates can be used until the expiration date printed on the label. Once the airtight packaging is opened, if any plates remain, they should be stored in the provided bag until the expiration date. For transportation, stability studies have shown that the plates can remain at **6 – 25 °C** for **3 days** or at **25 – 40 °C** for **24 hours** without affecting the product's performance.

USAGE

Place the plates in the incubator (35 – 37 °C) for 30 – 45 minutes until they dry. Inoculate the plates with the sample to be tested as soon as possible after receiving it from the laboratory. Incubate the plates upside down at 35 – 37 °C in an aerobic atmosphere for 18 – 24 hours.

READING AND INTERPRETATION

After 18 - 24 hours of incubation, check the plates for possible growth of multi-drug-resistant bacteria.

- The pathogenic CPE *E. coli* shows pink to reddish colonies of 2 to 3 mm in size.
- *Klebsiella*, *Enterobacter*, *Serratia*, and *Citrobacter* show blue green to blue colonies.
- *Pseudomonas* shows green colonies.


- *Acinetobacter* shows beige to yellow colonies.

LIMITATIONS OF THE METHOD


Some *E. coli* strains may show pale pink to colorless colonies. In all cases, additional confirmation should be performed using biochemical tests or other methods.

GENERAL CHARACTERISTICS OF QUALITY CONTROL


| Microorganism | ATCC | Characteristics of colonies |
|------------------------------|----------------------|-----------------------------|
| <i>Escherichia coli</i> IMP | NCTC 13476 | dark pink to reddish |
| <i>Klebsiella pneumoniae</i> | ATCC® BAA 1705 | metallic blue |
| <i>Klebsiella pneumoniae</i> | KPC NCTC 13438 | Blue metallic |
| <i>Enterococcus faecalis</i> | 29212 | Inhibited |
| <i>K. pneumoniae</i> | ESBL ATCC® 700603 | mostly inhibited |


CPE *E. coli*

CPE Coliforms

Pseudomonas

Acinetobacter

WASTE DISPOSAL OF WASTE

Materials that do not show any growth can be considered non-hazardous waste and discarded accordingly.

Plates that develop colonies must be disposed of in accordance with the instructions for contaminants or potential contaminants.

The laboratory is responsible for the proper management of infectious waste in accordance with its nature and degree of risk and should handle and discard it (or assign its management and disposal) in accordance with the applicable regulations.

SPECIFICATIONS

CHROMagar™ mSuperCARBA™ - CE

| PRODUCT | CODE | PACKING | STORE | SELF LIFE |
|-----------|--------|-----------|-----------|-----------|
| Plate 9cm | 010710 | 10 pieces | 2 – 12 °C | 45 days |
| Plate 6cm | 050710 | 10 pieces | 2 – 12 °C | 45 days |

Produced in Greece by the company Bioprepere in accordance with the requirements of the European Directive 2017/746.

BASIC UDI-DI: 5212037714010401WF. EDMA (14 01 04 01) non-chromogenic media (Plates).

The Bioprepere company has been certified according to the standards: EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004

BIBLIOGRAPHY

García-Fernández S., Hernández-García M., Valverde A., Ruiz-Garbajosa P., Morosini MI., Cantón R.- Servicio de Microbiología, Hospital Universitario Ramón y Cajal and Instituto Ramón y Cajal de Investigación Sanitaria (IRYCIS), Madrid, Spain.: CHROMagar mSuperCARBA performance in carbapenem-resistant Enterobacteriaceae isolates characterized at molecular level and routine surveillance rectal swab specimens.

M. Gaskin, D.Yamamura, J. Korver Hamilton Regional Laboratory Medicine Program ECCMID Madrid April 2018.: Validation of Colorex™ ESBL/mSuperCARBA™ bi-plate on WASP®/WASPLab® to screen for ESBL and CPE.

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IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr www.bioprepare.gr