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The evaluation of chromogenic medium for quantitative analysis of B.cereus in Korean traditional sauce. 2012. Sangmo Kang, Bumsun Lee, Junwhaun Yang, Chyungte Kim of Nongshim Foods Safety Center Ltd. Poster, 2011.

Preliminary study on enumeration medium for presumptive Bacillus cereus

2013. ADRIA Normandie, 2011.

Bacillus cereus food poisoning and testing for Bacillus cereus 2013.

Professor Shigeko Ueda, PhD Kagawa Nutrition university laboratory of Hygiene The Chemical Times 2013 No.2

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TECHNICAL DATA SHEET

PRODUCT: **CHROMagar™ Mycoplasma**
REFERENCE: **050865**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

Chromogenic medium for the color differentiation of *Mycoplasma bovis*.

INTRODUCTION

Infections in cattle can lead to various clinical symptoms, such as pneumonia, arthritis, mastitis, and keratoconjunctivitis, none of which are specific to *M. bovis*. Therefore, laboratory diagnosis is crucial. The identification of *M. bovis* from nasal swabs or bronchoalveolar lavage samples can be done by directly plating the sample on CHROMagar™ Mycoplasma plates. The detection of *M. bovis* with a red coloration allows easy identification, resembling a fried egg appearance under a binocular microscope after 3-7 days of incubation at 37°C in a CO₂ atmosphere.

PRINCIPLE OF THE METHOD

Peptones and yeast extract provide the necessary nutrients to the medium. The chromogenic substrates are broken down by specific microbial enzymes that release insoluble colored compounds, which adhere to the microbial cell membrane. *Bacillus cereus* uses only one of the chromogenic substrates and forms colonies with a purple color. If none of the substrates are used, there will be natural or white colonies. The salts maintain the osmotic balance in the medium. Agar provides trace elements and electrolytes and solidifies the medium.

FORMULA	g/litre
Chromogenic and selective mix	0,3
Peptone and yeast extract	25,0
Mineral Salts	6,3
Growth factors	2,2
Agar	10,0

Appearance: Light beige, transparent

Final pH 7.6 ± 0.2 at 25 °C.

PRECAUTIONS

CHROMagar™ B. cereus is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE AND TRANSPORT CONDITIONS

The plates should be stored at **2 – 12°C** in their packaging until use. Prolonged storage at temperatures below **2°C** can cause excessive moisture buildup in the medium, posing a risk of contamination. Freezing, even momentarily, will destroy the medium. Excessive heating should also be avoided. The plates can be used until the expiration date indicated on the label. For transportation, stability studies have shown that the plates can remain at **6 – 25°C for 4 days** or at **25 – 40°C for 48 hours** without affecting the performance of the product.

INSTRUCTIONS FOR USE

Place the plates in the incubator (35 – 37°C) for 30 – 45 minutes until dry. Plate the samples to be examined as soon as possible after collection from the laboratory. Incubate at 28 – 30°C in an aerobic atmosphere. After 18 – 24 hours, check the plates for the following color changes.

READING AND INTERPRETATION


After 24 hours of incubation, check the plates for possible growth of *Yersinia*. *Bacillus cereus* forms blue colonies with a white halo in the medium, measuring 2 to 4 mm. Other bacilli are inhibited or form blue colonies without halos in the medium. Most Gram-negative bacteria are inhibited.

LIMITATIONS OF THE METHOD

It is recommended to conduct biochemical tests on colonies from pure cultures for full identification.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	Colony Growth / Color
<i>Bacillus cereus</i> ATCC® 9634	Blue colonies with a white halo, size 2 to 4 mm.
<i>Bacillus subtilis</i> ATCC® 6633	Inhibited or forms blue colonies without halos in the medium.
<i>Candida albicans</i> ATCC® 10231	Inhibited
<i>Escherichia coli</i> ATCC® 25922	Inhibited

Mycoplasma bovis

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHROMagar™ B. cereus - **CE**

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
9cm Plate	010517	10 pieces	2 – 12 °C	2 months
6cm Plate	050517	10 pieces	2 – 12 °C	2 months

Produced in Greece by the company Bioprep, in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 5212037714010402WH. EDMA Code: (14 01 04 02) Chromogenic ID Media (Plates). Bioprep is certified according to the standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 / DY8d/1348/2004.

REFERENCES

The evaluation of chromogenic medium for quantitative analysis of B.cereus in Korean traditional sauce. 2012. Sangmo Kang, Bumsun Lee, Junwhaun Yang, Chyungte Kim of Nongshim Foods Safety Center Ltd. Poster, 2011.

Preliminary study on enumeration medium for presumptive Bacillus cereus
2013. ADRIA Normandie, 2011.

Bacillus cereus food poisoning and testing for Bacillus cereus 2013.
Professor Shigeko Ueda, PhD Kagawa Nutrition university laboratory of Hygiene The Chemical Times 2013 No.2

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr