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6th 2024

DESCRIPTION

The tripartite plate CHROMagar ORIENTATION – MAC CONKEY AGAR No2 – CHROMagar CANDIDA is used for the isolation and identification of gram(-) bacteria, gram(+) cocci (*Enterococcus*, *E. coli*, *Klebsiella*, *Proteus*) (CHROMagar™ ORIENTATION). It is used for the culture, isolation, and identification of enterobacteria and enterococci (MAC CONKEY AGAR No2). Finally, it is used for the isolation and color-based identification of *Candida* strains (CHROMagar CANDIDA).

METHOD PRINCIPLE

CHROMagar ORIENTATION is used for the color differentiation of pathogens in the urogenital tract, with *E. Coli* producing red colonies, *Klebsiella* producing metallic blue, *Proteus mirabilis* showing transparent colonies with a brown halo, and *Enterococcus faecalis* producing blue-green colonies with a blue halo. This enables easy visual detection of cultures with mixed flora.

FORMULA	g/litre
CHROMagar™ ORIENTATION	
Chromogenic mix	1,0
Peptone, yeast and extract	17,0
Agar	15,0

Appearance: Light beige, clear

Final pH 7.0 ± 0.2 at 25 °C.

MAC CONKEY AGAR No. 2 is a nutrient medium in which, with the addition of bile salts No. 2, both enterobacteria and enterococci can grow. Gram-negative enterobacteria that ferment lactose produces red or pink colonies. The bile salt No. 2 and crystal violet inhibit the growth of gram-positive cocci, except for *Enterococcus*. The pH indicator, neutral red, changes color with the fermentation of lactose.

FORMULA	g/litre
MAC CONKEY AGAR No. 2	
Peptone	20.0
Lactose	10.0
Bile Salts No. 2	1.5
Sodium chloride	5.0
Neutral red	0.05
Crystal violet	0.001
Agar No. 2	15.0

Appearance: Pink-red, clear.

Final pH 7.0 ± 0.2 °C at 25 °C.

CHROMagar CANDIDA is a new differential chromogenic medium that allows for the selective isolation and color differentiation of *Candida fungi*. In addition to the reliable detection of *Candida albicans*, the detection of *Candida glabrata*, *Candida tropicalis*, and *Candida krusei*, which tend to become more resistant to antifungal agents, is important. This enhances antifungal treatment, especially in immunocompromised patients.

FORMULA	g/litre
CHROMagar CANDIDA	
Peptone	10.2
Chromogenic mix	22.0
Chloramphenicol	0.1
Agar	15,0

Appearance: Light beige, clear

Final pH 6.1 ± 0.2 at 25 °C.

PRECAUTIONS

CHROMagar ORIENTATION – MAC CONKEY AGAR No2 – CHROMagar CANDIDA is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The plates should be stored at **2–12°C** inside their packaging until they are ready to be used. Prolonged storage at temperatures below **2°C** can create excessive moisture within the medium, posing a risk of contamination. Freezing, even temporarily, destroys the medium. Excessive heating should also be avoided. The plates can be used up until the expiration date printed on the label. For transport, stability studies have shown that the plates can remain at **6–25°C** for **4 days** or at **25–40°C** for **48 hours** without affecting the product's performance.

HOW TO USE

Place the plates in an incubator (35–37°C) for 30–45 minutes.

Inoculate the sample as soon as possible after collection and streak with successive dilutions for isolated colonies.

Incubate under aerobic conditions at 35–37°C for 24 hours. Consult the table for colony color.


LIMITATIONS OF THE METHOD

In CHROMagar™ ORIENTATION, some strains of *E. coli* may produce colonies with a pale pink to colorless appearance. Therefore, confirmation with biochemical tests should be performed.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

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WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHROMagar ORIENTATION – MAC CONKEY AGAR No2 – CHROMagar CANDIDA - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Divided 9 cm plate	030121	10 pieces	2 – 8 °C	3 months

It is produced in Greece by the company Bioprepare in accordance with Regulation (EU) 2017/746.

Basic UDI-DI: 5212037714010490X8. EDMA: (14 01 04 90) Other Prepared Media in Plates. Bioprepare is certified according to the standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 / DY86/1348/2004.

REFERENCES

CHROMagar™ ORIENTATION

(1) Merlino, J. et al. 1996. Evaluation of CHROMagar Orientation for Differentiation and Presumptive Identification of Gram-Negative Bacilli and Enterococcus Species, J.C.M. 34: 1788-1793. (2) Samra, Z. et al. 1998. Evaluation of Use of a New Chromogenic Agar in Detection of Urinary Tract Pathogens. J.C.M. 36: 990-994.

MAC CONKEY AGAR No2

American Public Health Association (1950). Diagnostic Procedures and Reagents. 3rd edn. A.P.H.A., New York.

American Public Health Association (1946). Standard Methods for the examination of Water and Sewage. 9th edn. A.P.H.A., New York.

CHROMagar CANDIDA

CHROMagar Candida, a new differential isolation medium for the presumptive identification of clinically important Candida species.

Odds F. et al.

1994, Journal of Clinical Microbiology, 32 : 1923-1929. Use of CHROMagar Candida medium for isolation of yeasts from dental samples.

Beighton D. et al.

1995. Journal of Clinical Microbiology, 33 : 3025-3027.

Application of CHROMagar Candida for rapid screening of clinical specimens for Candida albicans, Candida tropicalis, Candida krusei, and Candida (Torulopsis) glabrata.

Pfaller M.A. et al.

1996. Journal of Clinical Microbiology, 34 : 58-61. Routine use of CHROMagar Candida medium for presumptive identification of Candida yeast species and detection of mixed fungal populations.

Bouchara J.P. et al.

1996. Technical Report, 202-208.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepare.gr