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CHROMagar™ ORIENTATION

(1) Merlino, J. et al. 1996. Evaluation of CHROMagar Orientation for Differentiation and Presumptive Identification of Gram-Negative Bacilli and Enterococcus Species, J.C.M. 34: 1788-1793. (2) Samra, Z. et al. 1998. Evaluation of Use of a New Chromogenic Agar in Detection of Urinary Tract Pathogens. J.C.M. 36: 990-994.

MAC CONKEY AGAR

American Public Health Association (1950). Diagnostic Procedures and Reagents. 3rd edn. A.P.H.A., New York.

American Public Health Association (1946). Standard Methods for the examination of Water and Sewage. 9th edn. A.P.H.A., New York.

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6th 2024

DESCRIPTION

The split CHROMagar™ ORIENTATION – MAC CONKEY AGAR plate is used for the chromatic differentiation of urogenital pathogens. *E. coli* forms red colonies, *Klebsiella* forms metallic blue colonies, *Proteus mirabilis* forms clear colonies with a brown halo, and *Enterococcus faecalis* forms blue-green colonies with blue halos (CHROMagar™ ORIENTATION). It is also used for the isolation and identification of Enterobacteriaceae (MAC CONKEY AGAR).

PRINCIPLE OF THE METHOD

Meat peptones provide nitrogen, vitamins, minerals, and amino acids necessary for growth. Yeast extract is the source of vitamins, especially B vitamins. The chromogenic substrates are cleaved by specific microbial enzymes that release insoluble colored compounds which attach to the microbial cell membrane. Sodium chloride provides essential electrolytes for transport and regulates the osmotic balance of the medium. Bacteriological agar is the solidifying agent.

FORMULA	g/litre
CHROMagar™ ORIENTATION	
Chromogenic mix	1,0
Peptone, yeast and extract	17,0
Agar	15,0

Appearance: Light beige, clear

Final pH 7.0 ± 0.2 at 25 °C.

MacConkey Agar is a nutrient medium used for the isolation of Enterobacteriaceae. Gram-negative enterobacteria that ferment lactose produces red or pink colonies. Bile salts No. 3 and crystal violet inhibit the growth of Gram-positive cocci. The neutral red indicator changes color with the breakdown of lactose.

FORMULA	g/litre
MAC CONKEY AGAR	
Peptone	20.0
Lactose	10.0
Bile Salts No. 3	1.5
Sodium chloride	5.0
Neutral red	0.03
Crystal violet	0.001
Agar No. 2	15.0

Appearance: Pink-red, clear.

Final pH 7.1 ± 0.2 at 25°C.

PRECAUTIONS

CHROMagar™ ORIENTATION – MAC CONKEY AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The plates must be stored at **2–12 °C** in their packaging until the time of use. Prolonged storage at temperatures below **2 °C** creates excessive moisture inside the medium, increasing the risk of contamination. Freezing, even momentarily, destroys the medium. Excessive heating should also be avoided. The plates can be used up to the expiration date indicated on the label. For transportation, our stability studies have shown that the plates can remain at **17–25 °C** for **4 days** or at **27–40 °C** for **48 hours** without affecting the performance of the product

INSTRUCTIONS FOR USESample collection:

Urine should be first-morning and midstream.

Careful cleaning of the external genital area is required.

Use a sterile urine collection container and collect 20–25 ml of urine.

Semi-quantitative culture:

Urine should be cultured immediately after collection.

Otherwise, it can be stored for 24 hours at 4–6 °C.

Place the plates in an incubator, lid side down, at 35–37 °C for 30–45 minutes until dry.

Using a sterile calibrated loop of known diameter, either 0.01 ml (10 µl) or 0.001 ml (1 µl), transfer a drop of urine onto the divided agar plate by dipping the loop vertically, ensuring it does not touch the container walls.
 Spread the drop on each side of the plate using the loop.
 Incubate under aerobic conditions at 35–37 °C for 18–24 hours.

INTERPRETATION OF RESULTS

The number of colonies per ml (CFU/ml) is calculated as follows:

For urine samples plated with 0.01 ml (10 µl): number of colonies × 100

For urine samples plated with 0.001 ml (1 µl): number of colonies × 1000

CHROMagar™ ORIENTATION:

Escherichia coli → red to brick-red colonies

Klebsiella, *Enterobacter*, *Citrobacter*, *Serratia* → blue colonies, sometimes with a pink halo

Proteus, *Morganella*, *Providencia* → colorless to beige colonies with a brown halo

Proteus vulgaris → blue colonies with a brown halo

Pseudomonas → yellow colonies with a green discoloration in the medium

Acinetobacter → cream-colored colonies

Stenotrophomonas → colorless colonies

Enterococcus → small light blue colonies with a blue halo

Staphylococcus aureus → golden, opaque small colonies

Staphylococcus epidermidis → distinct cream colonies

Staphylococcus saprophyticus → pink, opaque small colonies

Streptococcus Group B (Strep B) → blue colonies

MACCONKEY AGAR:

Gram-positive cocci do not grow.

Gram-negative bacteria that ferment lactose produce pink to red colonies due to acid production.

Escherichia coli forms pink to red colonies surrounded by a hazy zone. This occurs due to bile salt precipitation as a result of pH decrease.

Gram-negative bacteria that do not ferment lactose form colorless colonies.

Proteus does not swarm.

LIMITATIONS OF THE METHOD


Some strains of *E. coli* may appear as light pink to colorless colonies.

Final identification may require additional tests, such as biochemical assays or immunological tests (e.g., latex agglutination).

QUALITY CONTROL

CHROMagar™ ORIENTATION

Microorganism	ATCC	Colony Characteristics
<i>Escherichia coli</i>	25922	Red – brick-red
<i>Enterococcus faecalis</i>	29212	Bluish green with a blue halo
<i>Klebsiella pneumoniae</i>	13883	Metallic blue
<i>Proteus mirabilis</i>	12473	Colorless with a brown halo
<i>Pseudomonas aeruginosa</i>	10145	Yellow with green pigmentation in the medium
<i>Salmonella typhimurium</i>	14028	Colorless
<i>Staphylococcus aureus</i>	25923	Light yellow
<i>Staphylococcus saprophyticus</i>	15305	Light pink
<i>Candida albicans</i>	60193	White, creamy

MAC CONKEY AGAR

Microorganism	ATCC	Growth	Colonies
<i>S. typhimurium</i>	14028	Good	Colorless – semi-transparent
<i>P. mirabilis</i>	12453	Good	Colorless, does not swarm
<i>Escherichia coli</i>	25922	Good	Pink or red, surrounded by a hazy zone
<i>E. faecalis</i>	29212	Inhibited	

Escherichia coli

P. mirabilis

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

REFERENCES

CHROMagar™ ORIENTATION – MAC CONKEY AGAR - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Divided Agar Plate 9 cm	020142	10 pieces	2 – 12 °C	3 months

Manufactured in Greece by the company Bioprep, in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 521203771414010405SJ. EDMA: (14 01 04 05) Bi-Plates – including Non-Chromogenic & Chromogenic Media combination. Bioprep is certified according to the standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 / Ministerial Decision ΔΥ86/1348/2004

REFERENCES

CHROMagar™ ORIENTATION

(1) Merlino, J. et al. 1996. Evaluation of CHROMagar Orientation for Differentiation and Presumptive Identification of Gram-Negative Bacilli and Enterococcus Species, J.C.M. 34: 1788-1793. (2) Samra, Z. et al. 1998. Evaluation of Use of a New Chromogenic Agar in Detection of Urinary Tract Pathogens. J.C.M. 36: 990-994.

MAC CONKEY AGAR

American Public Health Association (1950). Diagnostic Procedures and Reagents. 3rd edn. A.P.H.A., New York.

American Public Health Association (1946). Standard Methods for the examination of Water and Sewage. 9th edn. A.P.H.A., New York.\

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr