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6th 2024

DESCRIPTION

The bipartite CHROMagar™ SALMONELLA PLUS – HEKTOEN ENTERIC AGAR is used for the selective isolation and chromogenic identification of *Salmonella* spp. (CHROMagar™ SALMONELLA PLUS). Additionally, it facilitates the isolation and identification of *Salmonella* and *Shigella* species (HEKTOEN ENTERIC AGAR). Recommended for stool culture.

PRINCIPLE OF THE METHOD

CHROMagar™ SALMONELLA PLUS is a new chromogenic medium for the selective isolation of *Salmonella* spp., including *S. typhi*, *S. paratyphi*, and lactose (+) strains, from other Gram (-) bacteria. This medium complies with the ISO 6579:2002 standard. The white color of the medium helps in the identification of pink colonies of *Salmonella*.

COMPOSITION CHROMagar SALMONELLA PLUS	g/litre
Peptone and yeast extract	8.0
Chromogenic mix	1.3
Supplement	6,0
Salts	8,5
Agar	15,0

Appearance: White opaque

Final pH 7.6 ± 0.2 in 25 °C.

HEKTOEN ENTERIC AGAR is a differential nutrient medium used for the isolation and identification of pathogenic Enterobacteriaceae from clinical samples. The bile salts in the medium inhibit the growth of Gram (+) bacteria. *Enterobacteriaceae* that ferment lactose, sucrose, and salicin, producing acid, form yellow orange to salmon-coloured colonies. On the other hand, *Salmonella* and *Shigella* species form green to blue-green colonies.

COMPOSITION HEKTOEN ENTERIC AGAR	g/litre
Meat Peptone	12.0
Yeast Extract	3.0
Lactose	12.0
Sucrose	12.0
Salicin	2.0
Bile Salts No. 3	9.0
Sodium chloride	5.0
Sodium thiosulphate	5.0
Ammonium ferric citrate	1.5
Acid fuchsin	0.10
Bromothymol blue	0.064
Bacteriological Agar	14.0

Appearance: Light green transparent

Final pH 7.5 ± 0.2 at 25 °C.

PRECAUTIONS

CHROMagar™ SALMONELLA PLUS - HEKTOEN ENTERIC AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The plates should be stored at **2–12 °C** within their packaging until the time of use. Prolonged storage at temperatures below **2 °C** can create excess moisture inside the material, which may lead to contamination. Freezing, even momentarily, destroys the material. Excessive heating should also be avoided. The plates can be used up until the expiration date indicated on the label. For transport, stability studies have shown that the plates can remain at **6–25 °C for 4 days** or at **25–40 °C for 48 hours** without affecting the product's performance.

USAGE

Place the plates in the incubator (35–37 °C) for 30–45 minutes until they dry. Inoculate the sample on one edge of the plate, then make successive streaks with the loop in parallel lines to create isolated colonies. Incubate at 35–37 °C in an aerobic atmosphere for 18–24 hours. After 18–24 hours, check the plates for the following color changes.

INTERPRETATION OF RESULTS

After 24 hours of incubation, check the plates for possible *Salmonella* growth.

CHROMagar™ SALMONELLA PLUS

- The pathogenic *Salmonella* produces purple colonies, ranging in size from 2 to 3mm.
- *Escherichia coli* produces colorless colonies.
- Other bacteria that develop or show limited growth produce blue metallic colonies (e.g., *Citrobacter*, *Enterobacter*, *Aeromonas*, etc.).

HEKTOEN ENTERIC AGAR

- *Salmonella spp.* produces blue-green colonies with a black center due to H₂S production.
- *Shigella flexneri* produces blue colonies.
- *E. coli* produces large yellow-orange colonies.

LIMITATIONS OF THE METHOD

It is recommended to perform biochemical and other confirmatory tests required for identification. *Proteus* is a bacterium that can produce black colonies on H.E. Agar, similar to *Salmonella*. This occurs because it can produce H₂S from the breakdown of ammonium ferric citrate.

GENERAL CHARACTERISTICS OF QUALITY CONTROL


CHROMagar™ SALMONELLA PLUS:

Microorganism	Growth	Colony Color
<i>Salmonella spp</i> including <i>S. typhi</i> <i>S. paratyphi</i>	Good	Pink
<i>Salmonella lactose</i> (+)	Good	Pink
<i>Escherichia coli</i>	Good	Colourless
<i>Klebsiella pneumoniae</i>	Good	Blue
<i>Shigella sonnei</i>	Good	Colourless
<i>Proteus mirabilis</i>	Inhibited	
Gram (+) κόκκοι	Inhibited	

GENERAL CHARACTERISTICS OF QUALITY CONTROL

HEKTOEN ENTERIC AGAR

Microorganism	ATCC	Growth	Colony shape and color
<i>Escherichia coli</i>	25922	Good	Large, color yellow - beige
<i>Enterococcus faecalis</i>	29212	Inhibited	
<i>Enterobacter cloacae</i>	23355	Good	Large, color yellow - beige
<i>Proteus vulgaris</i>	13315	Inhibited	Greenish
<i>Shigella flexneri</i>	12022	Good	Green
<i>Salmonella typhimurium</i>	14028	Good	Blue - green with a black center

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHROMagar™ SALMONELLA PLUS - HEKTOEN ENTERIC AGAR - CE

ITEM	CODE	PACKAGE	STORAGE	SELF LIFE
Bi-plate Petri dish 9cm	020523	10 pieces	2 – 12 °C	3 months

It is produced in Greece from the Bioprep company in accordance with the requirement of the European guidance 98/79/EK. Code by EDMA 14 01 04. Bioprep company has been certified according to standards EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004

BIBLIOGRAPHY

CHROMagar™ SALMONELLA PLUS

(1) de Beaumont C., Breuil J., Dedicova D. Tran Q. 2006. Poster presented during ECCMID meeting.

Evaluation of a new chromogenic medium CHROMagar™ Salmonella Plus for the detection of Salmonella spp including lactose positive Salmonella, S.Typhi and S.Paratyphi.

HEKTOEN ENTERIC AGAR

King, S. & Metzger Appl. Microbiol. 16:577. 1968. King, S. & Metzger Appl.

Microbiol. 16:579, 1968.

Isenberg, Kominos & Siegel. Appl. Microbiol. 18:656. 1969. Hoben, Aston & Peterson Appl. Microbiol. 26:126. 1973.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepare.gr