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CHROMagar™ SALMONELLA PLUS

(1) de Beaumont C., Breuil J., Dedicova D. Tran Q. 2006. Poster presented during ECCMID meeting.

Evaluation of a new chromogenic medium CHROMagar™ Salmonella Plus for the detection of *Salmonella* spp including lactose positive *Salmonella*, *S.Typhi* and *S.Paratyphi*.

S.S. AGAR

Isenberg, H.D., Kominos, S., and Siegel, M. (1969). Isolation of salmonellae and shigellae from an artificial mixture of fecal bacteria. Appl. Microbiol., 18: 4, 656-659.

Leifson, E. (1935). New culture media based on sodium desoxycholate for the isolation of intestinal pathogens and for the enumeration of colon bacilli in milk and water. J. Pathol. Bacteriol., 40: 581-589.

Taylor, W.I., Harris, B. (1965) Isolation of shigellae. II. Comparison of plating media and enrichment broths. Am. J. Clin. Pathol. 44: 4, 476-479.

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6th 2024

DESCRIPTION

The split plate CHROMagar™ SALMONELLA PLUS – S.S. AGAR is used for the selective isolation and chromogenic identification of *Salmonella* spp. species (CHROMagar™ SALMONELLA PLUS). It is also used for the isolation and identification of *Salmonella* and *Shigella* species (S.S. AGAR). RECOMMENDED FOR FECAL CULTURE.

METHOD PRINCIPLE

CHROMagar SALMONELLA PLUS is a new chromogenic medium for the selective separation of *Salmonella* spp. species, including *S. typhi*, *S. paratyphi*, and lactose (+) from other gram (-) bacteria. This medium meets the ISO 6579: 2002 standard. The white color of the medium helps in the recognition of the pink colonies of *Salmonella*.

FORMULA	g/litre
CHROMagar SALMONELLA PLUS	
Peptone and yeast extract	8.0
Chromogenic mix	1.3
Supplement	6,0
Salts	8,5
Agar	15,0

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In S.S. AGAR, Beef Extract and Balanced Peptone No. 1 provide nitrogen, vitamins, minerals, and amino acids essential for growth.

Lactose is the fermentable carbohydrate that provides carbon and energy. A mixture of Bile Salts No. 3, Sodium citrate, and Brilliant Green inhibits Gram-positive (+) bacteria, most coliforms, and some *Proteus* spp., while allowing the growth of *Salmonella* spp.

Neutral Red is the pH indicator. Sodium thiosulphate and Ferric citrate allow the detection of H₂S produced by bacteria such as *Proteus* and certain strains of *Salmonella*, as they form colonies with black centers and clear halos. Bacteria that do not ferment lactose (presumed pathogens) form clear, transparent, or colorless colonies, while coliforms are sufficiently inhibited and form small colonies ranging from pink to red in color. This formulation is highly selective and is not recommended for the primary isolation of *Shigella*. Some strains of *Shigella* spp. may not grow.

FORMULA	g/litre
S.S. AGAR	
Beef Extract	5.0
Balanced Peptone No. 1	5.0
Lactose	10.0
Bile Salts No. 3	8.5
Sodium citrate	8.5
Sodium thiosulphate	8.5
Ferric citrate	1.0
Brilliant Green	0.00033
Neutral Red	0.025
Agar No. 2	13.5

Appearance: Red - Pink, translucent

Final pH 7.4 ± 0.2. at 25 °C

PRECAUTIONS

CHROMagar™ SALMONELLA PLUS – S.S. AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The plates should be stored at **2 – 12°C** within their packaging until they are ready for use. Prolonged storage at temperatures below **2°C** may cause excessive moisture within the material, risking contamination. Freezing, even briefly, destroys the material. Also, excessive heating should be avoided. The plates can be used until the expiration date indicated on the label. For transportation, stability studies have shown that the plates can remain at **6 – 25°C** for **4 days** or at **25 – 40°C** for **48 hours** without affecting the product's performance.

METHOD OF USE

Place the plates in the incubator (35 – 37°C) for 30 – 45 minutes.

Inoculate the sample as soon as possible after collection and streak with successive dilutions to obtain isolated colonies.

Incubate under aerobic conditions at 35 – 37°C for 24 hours. Refer to the table for colony color.

READING AND INTERPRETATION

CHROMagar™ SALMONELLA PLUS:

After 24 hours of incubation, examine the plates for possible *Salmonella* growth. Pathogenic *Salmonella* appears as purple colonies, 2 to 3 mm in size. *Escherichia coli* appears as colorless colonies. Other bacteria that grow or exhibit limited growth present metallic blue colonies (*Citrobacter*, *Enterobacter*, *Aeromonas*, etc.).

S.S. AGAR:

Enterobacteriaceae are differentiated by their ability to ferment lactose. *Salmonella* and *Shigella* species that do not ferment lactose appear as colorless colonies. *Salmonella* species that produce H₂S show colonies with black centers. *E. coli* produces pink to red colonies. *Proteus* may grow on SS Agar, forming gray to black-centered colonies due to hydrogen sulfide production. *Enterococcus faecalis* is partially inhibited, forming colorless colonies.

LIMITATIONS OF THE METHOD

CHROMagar SALMONELLA PLUS:

Final identification should be performed using biochemical and serological tests (e.g., Microgen *Salmonella* Latex agglutination test, code: F42), which can be carried out directly from the suspicious purple *Salmonella* colonies. *Salmonella Dublin* may produce colorless colonies, but this *Salmonella* strain is very rare. Some *E. coli* strains may show slightly purple colonies. Certain *Pseudomonas* strains may also form purple colonies. This can be ruled out with an oxidase test (Oxidase Test Strips – Code: MID-61G), which can be performed directly on the suspect colony.


S.S. AGAR:

The incorporation of Brilliant Green makes this medium highly selective and has been shown to inhibit the growth of some *Shigella* species. Therefore, SS Agar is best used in combination with CHROMagar™ SALMONELLA PLUS (010315), HEKTOEN ENTERIC AGAR (010047), or MAC CONKEY AGAR (010065). Bile salts may crystallize over time and before the expiry date of the plate. This results in white spots on the agar surface resembling mold or cobwebs. These do not affect the performance of the medium and will reabsorb into the agar after incubation. Some *Shigella* strains, such as *S. sonnei* and *S. dysenteriae*, may ferment lactose slowly, and colonies may appear lactose-positive after incubation for two or more days. It is recommended that biochemical and other confirmatory tests be performed for proper identification.

QUALITY CONTROL

CHROMagar™ SALMONELLA PLUS

Microorganism	Colony Characteristics
<i>Salmonella typhimurium</i> ATCC® 14028	Purple
<i>Salmonella abaeutuba</i> ATCC® 35640	Purple
<i>Salmonella enteritidis</i> ATCC® 13076	Purple
<i>Salmonella arizonae</i> CIP 8230	Purple
<i>Escherichia coli</i> ATCC® 25922	Colorless
<i>Citrobacter freundii</i> ATCC® 8090	Metallic blue
<i>Staphylococcus aureus</i> ATCC® 25923	Inhibited
<i>Pseudomonas aeruginosa</i> ATCC® 9027	Inhibited

S.S. AGAR

Microorganism	ATCC	Colony Appearance and Color
<i>Escherichia coli</i>	25922	Partial inhibition, Red colonies
<i>Enterococcus faecalis</i>	29212	Partial inhibition
<i>Salmonella typhimurium</i>	14028	Colorless with black center
<i>Shigella flexneri</i>	12022	Colorless colonies

S.S. AGAR: *Salmonella typhimurium* ATCC 14028

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

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Manufactured in Greece by Bioprep in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 5212037714010490X8. EDMA: (14 01 04 90) Other Prepared Media in Plates. Bioprep is certified according to the standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 / DY86/1348/2004.

REFERENCES

CHROMagar™ SALMONELLA PLUS

(1) de Beaumont C., Breuil J., Dedicova D. Tran Q. 2006. Poster presented during ECCMID meeting.

Evaluation of a new chromogenic medium CHROMagar™ Salmonella Plus for the detection of Salmonella spp including lactose positive Salmonella, S.Typhi and S.Paratyphi.

S.S. AGAR

Isenberg, H.D., Kominos, S., and Siegel, M. (1969). Isolation of salmonellae and shigellae from an artificial mixture of fecal bacteria. Appl. Microbiol., 18: 4, 656-659.

Leifson, E. (1935). New culture media based on sodium desoxycholate for the isolation of intestinal pathogens and for the enumeration of colon bacilli in milk and water. J. Pathol. Bacteriol., 40: 581-589.

Taylor, W.I., Harris, B. (1965) Isolation of shigellae. II. Comparison of plating media and enrichment broths. Am. J. Clin. Pathol. 44: 4, 476-479.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr