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TECHNICAL DATA SHEET

PRODUCT: **CHROMagar™ SERRATIA**
REFERENCE: **010824 - 050824**

Date 1st Edition:
11oç 2022
Date 3rd Revision:
4oç 2024

DESCRIPTION

CHROMagar SERRATIA is a chromogenic medium used for the isolation and color-based detection of *Serratia marcescens*.

INTRODUCTION

Serratia species are involved in hospital-acquired infections. In several countries, *Serratia marcescens* is frequently associated with outbreaks in intensive care units, especially in neonatal and pediatric units. Surveillance of hospital-acquired infections requires effective recovery of clinical isolates from stool, wound exudates, and respiratory samples to prevent cross-contamination issues and potentially fatal infections. In this context, CHROMagar™ developed CHROMagar™ Serratia, a culture medium ideally suited for detecting *Serratia marcescens* in stool samples. Moreover, *S. marcescens* can survive from days to months on surfaces, in distilled water, and in hand soap, making it a significant pathogen in hospital infections and sporadic outbreaks.

PRINCIPLE OF THE METHOD

Meat peptones provide nitrogen, vitamins, minerals, and amino acids necessary for growth. Yeast extract is a source of vitamins, especially those of the B group. The chromogenic substrates are broken down by specific microbial enzymes, releasing insoluble colored compounds that bind to the cell membrane of the microbes. Inhibitory agents prevent the growth of Enterobacteria and Gram-positive cocci. Sodium chloride supplies essential electrolytes for transport and regulates the osmotic balance of the medium. Bacteriological agar is the solidifying agent.

FORMULA	g/litre
Peptones	20,0
Sodium Chloride	5,0
Growth factors	1,7
Chromogenic and selective mix	0,8
Agar	15,0

Appearance: Light beige, clear

Final pH: 7.1 ± 0.2 at 25 °C

PRECAUTIONS

CHROMagar™ SERRATIA is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE AND TRANSPORT CONDITIONS

The plates should be stored at **2 – 12 °C** inside their packaging until use. Prolonged storage at temperatures below **2 °C** can cause significant moisture buildup within the medium, increasing the risk of contamination. Freezing, even momentarily, will destroy the medium. Excessive heating should also be avoided. The plates can be used until the expiration date indicated on the label. For transport, our stability studies have shown that the plates can remain at **6 – 25 °C** for **5 days** or at **25 – 40 °C** for **48 hours** without affecting the product's performance.

INSTRUCTIONS FOR USE

CHROMagar™ SERRATIA can be used with the following samples: stool, rectal swabs, and surface swabs. The samples are inoculated directly with serial dilutions onto the plate. If the plate is frozen, allow it to warm to room temperature before inoculation. Incubate under aerobic conditions at 35-37 °C for 18-24 hours.

INTERPRETATION OF RESULTS

Serratia marcescens > Green blue to metallic blue.

Pseudomonas > Colorless to yellow colonies.

Yeasts, Gram (+) bacteria, and other Gram (-) bacteria > Do not grow.

LIMITATIONS AND ADDITIONAL TESTS

At incubation temperatures below 35 °C, rare strains of *Serratia marcescens* may stain due to their natural red pigmentation. Final identification may require additional tests such as biochemical or immunological tests (latex agglutination).

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	Growth / Colony Color
<i>Serratia marcescens</i> ATCC 8100	Green blue to metallic blue colony color
<i>Pseudomonas aeruginosa</i> ATCC 9027	Colorless to yellow colonies

<i>Staphylococcus aureus</i> ATCC® 25923	Inhibited
<i>Enterococcus faecalis</i> ATCC® 29212	Inhibited
<i>Escherichia coli</i> ATCC® 25922	Inhibited
<i>Candida albicans</i> ATCC® 10231	Inhibited

Serratia marcescens

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHROMagar™ SERRATIA - C€

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
9 cm plate	010824	10 pieces	2 – 12 °C	3 months
6 cm plate	050824	10 pieces	2 – 12 °C	3 months

Manufactured in Greece by the company Bioprepere in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 5212037714010402WH. EDMA: (14 01 04 02) Chromogenic ID Media (Plates). Bioprepere is certified according to the standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 / DY8d/1348/2004.

REFERENCES

(1) Merlino, J. et al. 1996. Evaluation of CHROMagar Orientation for Differentiation and Presumptive Identification of Gram-Negative Bacilli and Enterococcus Species, J.C.M. 34: 1788-1793. (2) Samra, Z. et al. 1998. Evaluation of Use of a New Chromogenic Agar in Detection of Urinary Tract Pathogens. J.C.M. 36: 990-994.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepere.gr