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CHROMagar™ VRE

Evaluation of CHROMagar compared with enterococcosel broth for the isolation of vancomycin resistant enterococci. 2012.

Hindley et al Department of microbiology and infectious diseases ASM 2012.

Evaluation of Three Chromogenic Media for Detection of Vancomycin-Resistant Enterococci in a tertiary-care Hospital M.L. Miller et al. 2011. Kingston General Hospital, ON, Canada. Poster P26, CACMID 2011.

Evaluation of Three Commercial Chromogenic Media and BEAA + van 6ug/mL for the Detection of Vancomycin-Resistant Enterococcus (VRE). 2011.

Kornherr, Department of Microbiology Gamma Dynacare Medical Laboratories, Ottawa and Toronto, Ontario, Canada. ASM Meeting Poster 2010.

CHROMagar™ mSuperCARBA™

García-Fernández S., Hernández-García M., Valverde A., Ruiz-Garbajosa P., Morosini M.L., Cantón R.- Servicio de Microbiología, Hospital Universitario Ramón y Cajal and Instituto Ramón y Cajal de Investigación Sanitaria (IRYCIS), Madrid, Spain.: CHROMagar mSuperCARBA performance in carbapenem-resistant Enterobacteriaceae isolates characterized at molecular level and routine surveillance rectal swab specimens.

M. Gaskin, D.Yamamura, J. Korver Hamilton Regional Laboratory Medicine Program ECCMID Madrid April 2018.: Validation of Colorex™ ESBL/mSuperCARBA™ bi-plate on WASP®/WASPLab® to screen for ESBL and CPE.

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TECHNICAL DATA SHEET

PRODUCT: **CHROMagar™ VRE –
CHROMagar™ mSuperCARBA™**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

REFERENCE: **020857**

DESCRIPTION

The bipartite CHROMagar™ VRE – CHROMagar™ mSuperCARBA™ plate is used for the detection of Vancomycin-Resistant Enterococci (CHROMagar™ VRE), as well as for the detection of Gram-negative bacteria with reduced susceptibility to most carbapenem agents (CHROMagar™ mSuperCARBA™).

PRINCIPLE OF THE METHOD

In CHROMagar™ VRE, meat peptones provide nitrogen, vitamins, minerals, and amino acids essential for microbial growth. Yeast extract serves as a source of vitamins, particularly from the B group. Sodium chloride supplies essential electrolytes for nutrient transport and helps regulate the osmotic balance of the medium. Bacteriological agar is the solidifying agent. Vancomycin inhibits the growth of vancomycin-sensitive enterococci as well as many other bacteria.

COMPOSITION	g/litre
CHROMagar™ VRE	
Chromogenic mix	27.3
Peptones and yeast extract	20.0
Sodium chloride	8.0
Agar	15.0
Vancomycin	60mg

Appearance: white non-transparent.

Final pH 6.8 ± 0.2 at 25 °C.

In CHROMagar™ mSuperCARBA™, peptones provide the necessary nutrients for the medium. A blend of antibiotics ensures the selective growth of carbapenemase-producing organisms such as KPC, NDM, VIM, IMP, and OXA, with a detection limit of 10 CFU/ml. Three chromogenic substrates are hydrolyzed by microbial enzymes, releasing insoluble colored compounds that bind to the bacterial cell membrane during proliferation, resulting in colonies with specific colors. *Escherichia coli* produces enzymes -glucuronidase and β -galactosidase, leading to pink to brick-red colonies. *Klebsiella*, *Enterobacter*, *Serratia*, and *Citrobacter* produce β -glucuronidase, resulting in bluish-green colonies. If none of the substrates are used by the organism, the colonies remain colorless or white. The salts maintain osmotic balance in the medium, while agar provides trace elements and electrolytes, as well as solidifying the medium.

COMPOSITION	g/litre
mSuperCARBA™	
Chromogenic and selective mix	0,8
Peptones	20,0
Salts	5,0
Agar	15,0
Growth factors	1,7

Appearance: light beige clear

Final pH 7.0 ± 0.2 at 25 °C.

PRECAUTIONS

To CHROMagar™ VRE – CHROMagar™ mSuperCARBA™ is an in vitro diagnostic laboratory material and must be handled only by qualified laboratory personnel. This product contains peptones and extracts of animal origin. Certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious, and standard safety precautions should be followed (they must not be ingested or inhaled). Plates should always be handled with gloves and within a Class II laminar flow hood to avoid contamination, especially from saprophytic fungi. If a plate is cracked or the packaging is damaged, do not use it. Do not use plates that show signs of microbial contamination. The agar layer should be 4–5 mm thick and the material should be free of cracks, dryness, or any other signs of deterioration. After the expiration date, the material is unsuitable for use. In case of contact with skin, wash immediately with plenty of water and soap. Positive samples must be disposed of according to hygiene regulations for the handling of infectious material.

STORAGE AND TRANSPORT CONDITIONS

The plates must be stored at **2 – 12 °C** in their original packaging until use. Prolonged storage below **2 °C** can create excess moisture in the medium, increasing the risk of contamination. Freezing, even briefly, will destroy the medium. Excessive heating should also be avoided. The plates can be used up to the expiration date indicated on the label. Stability studies have shown that for transportation, the plates can be kept at **6 – 25 °C** for up to **3 days** or at **25 – 40 °C** for **24 hours** without affecting the performance of the product.

INSTRUCTIONS FOR USE

Allow the plates to reach room temperature before use. The agar surface must be dry prior to inoculation. Inoculate the sample on one edge of the plate, then perform streaking with a loop in parallel lines to create isolated colonies. Incubate the plates aerobically at **35 – 37 °C** for **24 to 48 hours**.

RESULT INTERPRETATION

CHROMagar™ VRE

VRE faecium and VRE faecalis form colonies with a pink to mauve color.

VRE gallinarum and VRE casseliflavus show poor growth and form blue colonies.

CHROMagar™ mSuperCARBA™

Pathogenic CPE *E. coli* produces pink to brick-red colonies, 2 to 3 mm in size.

Klebsiella, *Enterobacter*, *Serratia*, and *Citrobacter* produce bluish green to blue colonies.

Pseudomonas produce green colonies.

Acinetobacter appears as beige to yellow colonies.

METHOD LIMITATIONS

CHROMagar™ VRE

Definitive identification may require additional testing. Some rare strains of *Lactobacilli* and *Pediococcus* may occasionally develop mauve colonies. However, they can be differentiated by a PYR test (l-pyrrolidonyl β-naphthylamide): PYR positive indicates VRE *Enterococci*, PYR negative indicates *Lactobacilli* and *Pediococcus*. After 24 hours of incubation, some rare strains of *E. gallinarum* may occasionally develop mauve colonies. The vancomycin resistance type (R type) should be determined and confirmed by further tests.

CHROMagar™ mSuperCARBA™



Some *E. coli* strains may display pale pink to colorless colonies.

In all cases, further confirmation is required using biochemical tests or other methods.

GENERAL QUALITY CONTROL CHARACTERISTICS


CHROMagar™ VRE

Μικρόβιο	Growth / Colony Color
VRE <i>Enterococcus faecalis</i> ATCC® 51299	Mauve, small colonies
<i>Enterococcus faecalis</i> ATCC® 29212	Inhibited
<i>Escherichia coli</i> ATCC 25922	Inhibited
<i>Staphylococcus aureus</i> ATCC 25923	Inhibited
<i>Candida albicans</i> ATCC® 10231	Inhibited


CHROMagar™ mSuperCARBA™


Microorganism	ATCC	Colony Characteristics
<i>Escherichia coli</i> IMP	NCTC 13476	Deep pink to brick-red
<i>Klebsiella pneumoniae</i>	ATCC® BAA 1705	Metallic blue
<i>Klebsiella pneumoniae</i>	KPC NCTC 13438	Metallic blue
<i>Enterococcus faecalis</i>	29212	Inhibited
<i>Candida albicans</i>	60193	Inhibited

CPE *E. coli*

CPE Coliforms

Pseudomonas

Acinetobacter

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

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Manufactured in Greece by Bioprepare in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 5212037714010404WM EDMA Code: (14 01 04 04) Chromogenic Resistance Marker Media (Plates) – MRSA, VRE, ESBL. Bioprepare is certified under the standards: EN ISO 9001:2015, ELOT EN ISO 13485:2016, DY8d/1348/2004.

REFERENCES

CHROMagar™ VRE

Evaluation of CHROMagar compared with enterococcosel broth for the isolation of vancomycin resistant enterococci. 2012.

Hindley et al Department of microbiology and infectious diseases ASM 2012.

Evaluation of Three Chromogenic Media for Detection of Vancomycin-Resistant Enterococci in a tertiary-care Hospital M.L. Miller et al. 2011. Kingston General Hospital, ON, Canada. Poster P26, CACMID 2011.

Evaluation of Three Commercial Chromogenic Media and BEAA + van 6µg/mL for the Detection of Vancomycin-Resistant Enterococcus (VRE). 2011.

Kornherr, Department of Microbiology Gamma Dynacare Medical Laboratories, Ottawa and Toronto, Ontario, Canada. ASM Meeting Poster 2010.

CHROMagar™ mSuperCARBA™

García-Fernández S., Hernández-García M., Valverde A., Ruiz-Garbajosa P., Morosini M., Cantón R.- Servicio de Microbiología, Hospital Universitario Ramón y Cajal and Instituto Ramón y Cajal de Investigación Sanitaria (IRYCIS), Madrid, Spain.: CHROMagar mSuperCARBA performance in carbapenem-resistant Enterobacteriaceae isolates characterized at molecular level and routine surveillance rectal swab specimens.

M. Gaskin, D.Yamamura, J. Korver Hamilton Regional Laboratory Medicine Program ECCMID Madrid April 2018.: Validation of Colorex™ ESBL/mSuperCARBA™ bi-plate on WASP®/WASPLab® to screen for ESBL and CPE.

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IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepare.gr