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Evaluation of CHROMagar compared with enterococcosel broth for the isolation of vancomycin resistant enterococci. 2012.

Hindley et al Department of microbiology and infectious diseases ASM 2012.

Evaluation of Three Chromogenic Media for Detection of Vancomycin-Resistant Enterococci in a tertiary-care Hospital M.L. Miller et al. 2011. Kingston General Hospital, ON, Canada. Poster P26, CACMID 2011.

Evaluation of Three Commercial Chromogenic Media and BEAA + van 6ug/mL for the Detection of Vancomycin-Resistant Enterococcus (VRE). 2011.

Kornherr, Department of Microbiology Gamma Dynacare Medical Laboratories, Ottawa and Toronto, Ontario, Canada. ASM Meeting Poster 2010.

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TECHNICAL DATA SHEET

PRODUCT: **CHROMagar™ VRE**
REFERENCE: **010360 – 050360**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

CHROMagar VRE allows easy and visible detection of vancomycin-resistant *Enterococcus faecalis* & *Enterococcus faecium* with high sensitivity and accuracy, directly from the color of the colony.

INTRODUCTION

The isolation of vancomycin-resistant *Enterococcus faecalis* & *Enterococcus faecium* is crucial because they are easily transmitted from person to person. In contrast, non-resistant enterococci are not easily transmitted and are not associated with outbreaks. There are two types of vancomycin-resistant enterococci. The first type of resistance (mainly VanC but also VanD, VanE, VanF, etc.) found in *E. gallinarum* and *E. casseliflavus* / *E. flavescens* shows low resistance to vancomycin. The second type of vancomycin resistance in enterococci is acquired resistance (VanA & VanB types), which is mainly observed in *E. faecium* and *E. faecalis*. Therefore, to prevent the spread of this resistance to more pathogens (e.g., *S. aureus*), it is vital to immediately detect the presence of any of these two species in the patient and accurately differentiate them from other enterococci.

PRINCIPLE OF THE METHOD

Meat peptones provide nitrogen, vitamins, minerals, and amino acids necessary for growth. Yeast extract is the source of vitamins, particularly B vitamins. Sodium chloride provides essential electrolytes for transport and regulates the osmotic balance of the medium. The bacteriological agar is the gelling agent.

FORMULA	g/litre
Chromogenic mix	27.3
Peptones and yeast extract	20.0
Sodium chloride	5,0
Agar	15,0
Vancomycin	60mg

Appearance: White opaque.

Final pH 6.8 ± 0.2 at 25 °C.

PRECAUTIONS

CHROMagar™ VRE is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The plates should be stored at **2–12 °C** within their packaging until the time of use. Prolonged storage at temperatures below **2 °C** can create excess moisture inside the material, which may lead to contamination. Freezing, even momentarily, destroys the material. Excessive heating should also be avoided. The plates can be used up until the expiration date indicated on the label. For transport, stability studies have shown that the plates can remain at **6–25 °C** for **4 days** or at **25–40 °C** for 48 hours without affecting the product's performance.

USAGE

The infectious material must be delivered to the laboratory promptly and protected from excessive heat and cold. If there is any delay in processing, the sample should be inoculated into an appropriate transport medium and kept at 2-12°C until inoculation. Allow the plates to warm to room temperature. The agar surface should be dry before inoculation. Inoculate the sample at one edge of the plate and then make successive streaks with the loop in parallel lines to create isolated colonies. Incubate the plates aerobically at 35-37°C for 24 to 48 hours.

INTERPRETATION OF RESULTS

VRE faecium & VRE faecalis will form colonies with a pink-mauve color.

VRE gallinarum & VRE casseliflavus will grow poorly and form blue colonies.

LIMITATIONS OF THE METHOD

Definitive identification may require additional testing.

Some rare varieties of *Lactobacilli* & *Pediococcus* may occasionally form mauve colonies. However, they can be differentiated with a PYR test (L-pyrrolidonyl β-naphthylamide):

PYR (+) -> VRE Enterococci

PYR (-) -> *Lactobacilli* & *Pediococcus*

After 24 hours of incubation, some rare strains of *E. gallinarum* may occasionally develop mauve colonies.



The Vancomycin resistance type (R) should be determined and confirmed by further testing.

It is recommended to use the oxidase test (Oxidase Test Strips-Code MID-61G) on colonies from pure culture for full identification.

Final identification may require additional tests such as biochemical or immunological assays (latex agglutination test).

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	Growth - Colony Color
VRE <i>Enterococcus faecalis</i> ATCC® 51299	Purple, small
<i>Enterococcus faecalis</i> ATCC® 29212	Inhibit
<i>Escherichia coli</i> ATCC® 25922	Inhibit
<i>Candida albicans</i> ATCC® 10231	Inhibit
<i>Staphylococcus aureus</i> ATCC® 25923	Inhibit

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHROMagar™ VRE - CE

ITEM	CODE	PACKAGE	KEEPING	LIFE SPAN
Plate 9cm	010360	10 pieces	2 – 12 °C	2 months
Plate 6cm	050360	10 pieces	2 – 12 °C	2 months

Produced in Greece by Bioprepare in accordance with the requirements of the European Directive 98/79/EK. FEK B2198/2-10-2009. CODE EDMA 14 01 04 02. Bioprepare company has been certified according to standards EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004.

BIBLIOGRAPHY

Evaluation of CHROMagar compared with enterococcosel broth for the isolation of vancomycin resistant enterococci. 2012. Hindley et al Department of microbiology and infectious diseases ASM 2012.

Evaluation of Three Chromogenic Media for Detection of Vancomycin-Resistant Enterococci in a tertiary-care Hospital M.L. Miller et al. 2011. Kingston General Hospital, ON, Canada. Poster P26, CACMID 2011.

Evaluation of Three Commercial Chromogenic Media and BEAA + van 6ug/mL for the Detection of Vancomycin-Resistant Enterococcus (VRE). 2011. Kornherr, Department of Microbiology Gamma Dynacare Medical Laboratories, Ottawa and Toronto, Ontario, Canada. ASM Meeting Poster 2010.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepare.gr