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CHROMagar™ *Y. enterocolitica*

Comparison of CHROMagar *Yersinia* and CIN Agar media for isolation of potentially virulent *Yersinia enterocolitica* in stools. Nicolas Renaud, a Laetitia Lecci, a René J. Courcol, a, b Michel Simonet, a, b Olivier Gaillota Lille University Medical Center and INSERM U1019, b Lille, France Journal of Clinical Microbiology p. 1184–1187 April 2013 Volume 51 Number 4 COMPLETE ARTICLE

Comparison of CHROMagar *Yersinia* and CIN Agar media for isolation of potentially virulent *Yersinia enterocolitica* in stools. Nicolas RENAUD et al. Microbiology Laboratory, Lille University Hospital, LILLE, France ASM Poster 2012.

CHROMagar™ O157

The use of selective and differential agars in the isolation of *Escherichia coli* O157 from dairy herds.

Wallace J.S. et al.

1996. Journal of Applied Bacteriology, 81 : 663-668.

Isolation of vero cytotoxin-producing *Escherichia coli* O157 from wild birds.

Wallace J.S. et al.

1997. Journal of Applied Microbiology, 82 : 399-404.

Reliability of CHROMagar O157 for the detection of enterohaemorrhagic *Escherichia coli* (EHEC) O157 but not EHEC belonging to other serogroups.

Bettelheim K.A.

1998. Journal of Applied Microbiology, 85 : 425-428.

The optimization of isolation media used in immunomagnetic separation methods for the detection of *Escherichia coli* O157 in foods.

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PRODUCT: **CHROMagar™ Y. enterocolitica**
CHROMagar™ O157

REFERENCE: **020709**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

The split plate CHROMagar™ Y. enterocolitica – CHROMagar™ O157 is used for the detection and direct differentiation of the pathogenic *Yersinia enterocolitica* (CHROMagar™ Y. enterocolitica). Additionally, it enables chromogenic differentiation and isolation of *E. coli* O157 (CHROMagar™ O157).

PRINCIPLE OF THE METHOD

In CHROMagar™ Y. enterocolitica, peptones provide the necessary nutrients for the medium. Selective agents inhibit the growth of Gram (+) cocci and certain Gram (-) bacteria. The chromogenic substrates are broken down by specific microbial enzymes, releasing insoluble colored compounds that bind to the cell membrane of the microorganisms. *Yersinia enterocolitica* uses only one of the chromogenic substrates and forms purple colonies. If none of the substrates is used, natural or white colonies will appear. Salts maintain the osmotic balance of the medium, while agar provides trace elements and electrolytes, and solidifies the medium.

COMPOSITION	g/litre
CHROMagar™ Y. enterocolitica	
Chromogenic mix	1,3
Peptones	20,0
Salts	5,0
Selective mix	0.1
Agar	15,0

Appearance: Light beige, transparent

Final pH 7.0 ± 0.2 at 25 °C.

CHROMagar™ E. coli O157 is the new solution for the isolation and identification of the enteropathogenic *E. coli* O157 from clinical samples and food samples. Almost all O157 strains (98%) form colonies with a bright pink color, while other coliforms form blue colonies. The compact pink colonies, without any diffusion of the color into the medium, the easy color differentiation of the colonies from other coliforms, and the rapid growth in the medium (18–24 hours) are the advantages of this new medium. The medium has good selectivity and inhibits the growth of Gram(+) and Gram(-) bacteria. Peptones and yeast extract provide the necessary nutrients for the medium. The chromogenic substrates are broken down by specific microbial enzymes, releasing insoluble colored compounds that bind to the cell membrane of the microorganisms. *E. coli* O157 uses two chromogenic substrates and forms colonies with a purple color. If none of the substrates is used, natural or white colonies will appear. The agar provides trace elements and electrolytes, and also solidifies the medium.

COMPOSITION	g/litre
CHROMagar™ O157	
Peptone and yeast extract	13.0
Chromogenic mix	1.2
Agar	15,0

Appearance: Light beige, transparent

Final pH 6.9 ± 0.2 at 25 °C.

PRECAUTIONS

CHROMagar™ Y. enterocolitica – CHROMagar™ O157 is an in vitro diagnostic laboratory material and must be handled only by qualified laboratory personnel. This product contains peptones and extracts of animal origin. Certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious, and standard safety precautions should be followed (they must not be ingested or inhaled). Plates should always be handled with gloves and within a Class II laminar flow hood to avoid contamination, especially from saprophytic fungi. If a plate is cracked or the packaging is damaged, do not use it. Do not use plates that show signs of microbial contamination. The agar layer should be 4–5 mm thick and the material should be free of cracks, dryness, or any other signs of deterioration. After the expiration date, the material is unsuitable for use. In case of contact with skin, wash immediately with plenty of water and soap. Positive samples must be disposed of according to hygiene regulations for the handling of infectious material.

STORAGE AND TRANSPORT CONDITIONS

The plates must be stored at **2 – 12 °C** within their packaging until use. Prolonged storage at temperatures below **2 °C** creates excess moisture in the material, which can lead to contamination. Freezing, even momentarily, damages the material. Excessive heating should also be avoided. The plates can be used up until the expiration date indicated on the label. For transport, stability studies have shown that the plates can remain at **6 – 25 °C** for **4 days** or at **25 – 40 °C** for **48 hours** without affecting the product's performance.

INSTRUCTIONS FOR USE

Allow the plates to warm to room temperature. The surface of the agar must be dry before inoculation. Inoculate the sample at one edge of the plate and then streak the sample into parallel lines with the loop to create isolated colonies. Incubate the plates aerobically at **35 – 37 °C** for 24 to 48 hours.

READING AND INTERPRETATION

CHROMagar™ *Y. enterocolitica*

After 24 hours of incubation, check the plates for possible growth of *Yersinia*. The pathogenic *Yersinia enterocolitica* forms purple colonies ranging from 2 to 3mm in size. Non-pathogenic *Yersinia enterocolitica* forms metallic blue colonies. Other bacteria that grow, but show limited growth, also present metallic blue colonies (e.g., *Citrobacter*, *Enterobacter*, *Aeromonas*, etc.).

CHROMagar™ O157

After 24 hours of incubation, check the plates for possible growth of *Escherichia coli* O157. *E. coli* O157 forms purple colonies. Other *Enterobacteriaceae* are inhibited or form colorless to blue colonies. Gram-positive bacteria do not grow.

METHOD LIMITATIONS

CHROMagar™ *Y. enterocolitica*

Some *Y. enterocolitica* may show moderate or no growth on the medium. Some rare strains of non-pathogenic *Yersinia* may form purple colonies (e.g., *Y. bercovieri*, *Y. mollareti*, *Y. kristensenii*, *Y. rohdei*, etc.). Biochemical tests are recommended on colonies from pure cultures for full identification.

CHROMagar™ O157

Final identification should be carried out using biochemical and serological tests (e.g., Microgen *E. coli*-O157 Latex Test, code: M44CE), which can be performed directly from the suspected purple colonies of *E. coli* O157.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

CHROMagar™ *Y. enterocolitica*


Microorganism	ATCC	Colony Characteristics
<i>Y. enterocolitica</i>	12982	Purple
<i>Y. enterocolitica</i>	biotype 1A	Metallic blue
<i>Escherichia coli</i>	25922	Inhibited
<i>Enterococcus faecalis</i>	29212	Inhibited
<i>Pseudomonas aeruginosa</i>	27853	Inhibited
<i>Citrobacter freundii</i>	9750	Moderate inhibition

Y. enterocolitica ATCC 12982 ($\mu\text{o}8$)

CHROMagar™ O157

Microorganism	Colony Growth / Color
<i>E. coli</i> O157 ATCC® 35150	Purple colonies
<i>Klebsiella pneumoniae</i> ATCC® 13883	Metallic blue colonies
<i>Escherichia coli</i> ATCC® 25922	Metallic blue colonies
<i>Enterococcus faecalis</i> ATCC® 29212	Inhibited

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CHROMagar™ Y. enterocolitica – CHROMagar™ O157 - 

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Split plate 9cm	020709	10 pieces	2 – 12 °C	2 months

Produced in Greece by Bioprep prepare according to Regulation (EU) 2017/746.

Basic UDI-DI: 5212037714010402WH. EDMA: (14 01 04 02) Chromogenic ID Media (Plates).

Bioprep prepare is certified according to the standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 DY86/1348/2004.

REFERENCES

CHROMagar™ Y. enterocolitica

Comparison of CHROMagar Yersinia and CIN Agar media for isolation of potentially virulent Yersinia enterocolitica in stools. Nicolas Renaud,a Laetitia Lecci,a René J. Courcol,a,b Michel Simonet,a,b Olivier Gaillota Lille University Medical Centera and INSERM U1019,b Lille, France Journal of Clinical Microbiology p. 1184–1187 April 2013 Volume 51 Number 4 COMPLETE ARTICLE

Comparison of CHROMagar Yersinia and CIN Agar media for isolation of potentially virulent Yersinia enterocolitica in stools. Nicolas RENAUD et al. Microbiology Laboratory, Lille University Hospital, LILLE, France ASM Poster 2012.

CHROMagar™ O157

The use of selective and differential agars in the isolation of Escherichia coli O157 from dairy herds.

Wallace J.S. et al.

1996. Journal of Applied Bacteriology, 81 : 663-668.

Isolation of vero cytotoxin-producing Escherichia coli O157 from wild birds.

Wallace J.S. et al.

1997. Journal of Applied Microbiology, 82 : 399-404.

Reliability of CHROMagar O157 for the detection of enterohaemorrhagic Escherichia coli (EHEC) O157 but not EHEC belonging to other serogroups.

Bettelheim K.A.

1998. Journal of Applied Microbiology, 85 : 425-428.

The optimization of isolation media used in immunomagnetic separation methods for the detection of Escherichia coli O157 in foods.

Ogden I.D. et al.

2001, Journal of Applied Microbiology, 91 : 1-7.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr www.bioprep prepare.gr