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- ATLAS, R.M., L.C. PARKS (1993) Handbook of Microbiological Media. CRC Press, Inc. London.
- BARON, E.J., L.R. PETERSON & S.M. FINEGOLD (1994) Bailey & Scott's Diagnostic Microbiology. 9th ed. Mosby-Year Book Inc. St Louis. MO. USA.
- ISENBERG, H.D. (1992) Clinical Microbiology Procedures Handbook. ASM Washington. DC. USA.
- ISO 11133:2014. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- MACKEY, J.P. & G.H. SANDYS (1966) Diagnosis of urinary tract infections. Brit. Med. J. 3:1.173.
- MURRAY, P.R., E.J. BARON, M.A. PFALLER, F.C. TENOVER & R.H. YOLKEN (1995) Manual of Clinical Microbiology 6th ed. ASM Washington. DC. USA.
- SANDYS, G.H. (1960) A new method of preventive swarming of *Proteus* sp. J. Med. Lab. Tech. 17:224.

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6th 2024

PRODUCT: CLED AGAR

(Cystine Lactose Electrolyte Deficient-Single Indicator)

REFERENCE: 010034 - 060034**DESCRIPTION**

A culture medium recommended for the detection and cultivation of urinary pathogenic microbes.

PRINCIPLE OF THE METHOD

In the presence of Bromthymol Blue, the microbes that ferment lactose produce yellow colonies, while those that do not ferment lactose give blue colonies. Bevis incorporated an additional dye of acid fucin which makes lactose positive bacteria to give bright pink colonies and lactose negatively to give colorless colonies. Due to the absence of salts in the material, the creep of the primate is inhibited.

FORMULA	g/litre
Lactose	10.00
Casein Peptone	4.00
Gelatin Peptone	4.00
Beef Extract	3.00
L-Cystine	0.128
Bromothymol Blue	0.02
Bacteriological Agar	15.00

Appearance: Green clear

Final pH 7.3 ± 0.2 at 25 °C.

PRECAUTIONS

CLED AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The petri dishes should be stored at 2 – 12 °C inside their packaging until use. Prolonged storage at temperatures below 2 °C creates excessive moisture within the medium, increasing the risk of contamination. Freezing, even momentarily, destroys the medium. Excessive heating should also be avoided. The petri dishes can be used until the expiration date printed on the label. For transportation, stability studies have shown that the petri dishes can remain at 18 - 25 °C for up to 7 days or at 30 - 40 °C for up to 48 hours without affecting the product's performance.

INSTRUCTIONS FOR USE

Urine samples should be inoculated onto the petri dish no later than one hour after collecting. Inoculate using a calibrated loop of 10µL or 1µL at the edge of the petri dish. Incubate the plates under aerobic conditions at $37 \pm 1^\circ\text{C}$ for 18-24 hours.

READING AND INTERPRETATION

After incubation, count all colonies that have appeared on the agar surface.

METHOD LIMITATIONS

The isolation of any urinary pathogen must be confirmed with further microbiological and biochemical tests.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	Growth	Collor
<i>Escherichia coli</i>	25922	Yellow.
<i>Proteus mirabilis</i>	7002	Blue.
<i>Salmonella typhimurium</i>	14028	Blue.
<i>Staphylococcus aureus</i>	25923	Yellow.

Escherichia coli ATCC 25922

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CLED AGAR - GR/CA01/GRM5/O/1 - CE

PRODUCT	CODE	PACKING	STORAGE	SELF LIFE
Plate 9cm	010034	10 pieces	2 – 12 °C	5 months
Bottle 100ml	060034	10 pieces	2 – 25 °C	12 months

Produced in Greece by Bioprep in accordance with the requirements of the European Directive 98/79/EK. FEK B2198/2-10-2009. CODE EDMA 14 01 02 01. Bioprep company has been certified according to standards EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004.

BIBLIOGRAPHY

- ATLAS, R.M., L.C. PARKS (1993) Handbook of Microbiological Media. CRC Press, Inc. London.
- BARON, E.J., L.R. PETERSON & S.M. FINEGOLD (1994) Bailey & Scott's Diagnostic Microbiology. 9th ed. Mosby-Year Book Inc. St Louis. MO. USA.
- ISENBERG, H.D. (1992) Clinical Microbiology Procedures Handbook. ASM Washington. DC. USA.
- ISO 11133:2014. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- MACKEY, J.P. & G.H. SANDYS (1966) Diagnosis of urinary tract infections. Brit. Med. J. 3:1.173.
- MURRAY, P.R., E.J. BARON, M.A. PFALLER, F.C. TENOVER & R.H. YOLKEN (1995) Manual of Clinical Microbiology 6th ed. ASM Washington. DC. USA.
- SANDYS, G H. (1960) A new method of preventive swarming of Proteus sp. J. Med. Lab. Tech. 17:224.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr