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TECHNICAL DATA SHEET

PRODUCT: **COLUMBIA C.N.A. (Horse Blood)** 
REFERENCE: **010251**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

Selective nutrient medium for the cultivation and isolation of Staphylococci and Streptococci from clinical samples.

PRINCIPLE OF THE METHOD

COLUMBIA C.N.A. (Horse Blood) contains ingredients of high nutritional value that provide vitamins, carbohydrates, and other organic elements. Sodium chloride provides the necessary minerals and maintains osmotic balance and electrolyte equilibrium. The addition of 5% horse blood enhances the nutritional properties of the medium. Additionally, hemolysis and the type of hemolysis of horse red blood cells aid in the identification of certain gram-positive cocci. The addition of colistin and nalidixic acid inhibits the growth of gram-negative bacteria.

COMPOSITION	g/litre
Peptone Mixture	20,0
Beef Extract	3,0
Sodium chloride	5.0
Corn Starch	1,0
Bacteriological Agar	15.0
Horse Blood	50ml
Colistin	10mg
Nalidixic acid	15mg

Appearance: Red – burgundy, opaque due to the addition of blood.

Final pH 7.3 ± 0.2 at 25°C.

PRECAUTIONS

COLUMBIA C.N.A. (Horse Blood) is an in-vitro diagnostic laboratory material and should only be handled by trained laboratory personnel. This material contains peptones and extracts of animal origin. Certificates regarding the source and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious and handled using standard safety precautions (avoid ingestion and inhalation). Always handle plates wearing gloves and within a Class II laminar flow cabinet to prevent contamination, especially by saprophytic fungi. Do not use plates that show signs of microbial contamination. The material should not be used after the expiration date. In case of skin contact, wash immediately with plenty of water and soap. Positive samples must be disposed of according to the appropriate hygiene regulations for infectious materials.

STORAGE AND TRANSPORT CONDITIONS

Petri dishes should be stored at **2–12 °C** in their packaging until use. Prolonged storage at temperatures below **2 °C** causes excessive moisture in the medium, increasing the risk of contamination. Freezing, even briefly, destroys the medium. Excessive heating should also be avoided. The plates can be used up to the expiration date indicated on the label. Regarding transport, our stability studies show that the plates can remain at **6–25 °C** for up to **4 days** or at **25–40 °C** for up to **48 hours** without affecting the product's performance.

INSTRUCTIONS FOR USE

Inoculate the plates using the streaking technique for isolated colonies.

Incubate at 35–37 °C for 24 hours under aerobic conditions.

INTERPRETATION OF RESULTS

Hemolytic streptococcus forms transparent or matte colonies, gray in color, surrounded by a zone of beta-hemolysis.

Pneumococcus appears as flat, transparent colonies surrounded by a small zone of green alpha-hemolysis.

Staphylococcus appears as matte, white to golden-yellow colonies with or without a zone of beta-hemolysis.

Gram-negative bacteria do not grow.

LIMITATIONS OF THE METHOD

Final identification must be carried out using biochemical and serological tests (e.g., Latex Agglutination Test), which can be performed directly from the suspected colonies. On COLUMBIA C.N.A. (Horse Blood) medium, growth of Gram-negative bacteria may be observed, although these are usually sensitive to the medium's antibiotics. This can occur when the bacterial count per ml of sample (CFU) is $\geq 10^5$. Additionally, some strains may exhibit resistance to the antibiotics in the medium.

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Streptococcus pyogenes ATCC 19615

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

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Produced in Greece by the company Bioprepere in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 5212037714010401WF. EDMA (14 01 04 01): Non-Chromogenic media (Plates). Bioprepere is certified according to the standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 / DY8d/1348/2004.

REFERENCES

- Ellner, P.D., Stoessel, C.J., Drakeford, E and Vasi, F. (1966). A new culture medium for medical bacteriology. Amer. J. Clin Pathol., 45:502-504.
- Goldberg, R.L., and Washington, J.A., (1976). Comparison of isolation of Haemophilus vaginalis (Corynebacterium vaginale) from Peptone-Starch-Dextrose Agar and Columbia Colistin-Nalidixic Acid Agar. J. Clin. Microbiol., 4:245-247.
- Thayer, D.D. and Martin, H. E. (1966). An improved medium for the cultivation of N. gonorrhoeae and N. meningitidis. Publ. Hlth. Report, 81:559-562.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepere.gr