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E-mail: bioprep1@otenet.gr www.bioprepare.gr

**TECHNICAL DATA SHEET**

PRODUCT: **COLUMBIA C.N.A. W/ 6% HORSE BLOOD – MAC CONKEY AGAR**

REFERENCE: **020261**

Date 1st Edition:

7th 2009

Date 4th Revision:

6th 2024

DESCRIPTION

Bipartite Petri Dish COLUMBIA C.N.A. W/ 6% HORSE BLOOD – MAC CONKEY AGAR. For the cultivation and isolation of aerobic Gram-positive cocci and Gram-negative bacteria from clinical samples.

PRINCIPLE OF THE METHOD

COLUMBIA C.N.A. W/ 6% HORSE BLOOD contains highly nutritious components that provide vitamins, carbohydrates, and other organic elements essential for bacterial growth. Sodium chloride supplies the necessary minerals and helps maintain osmotic and electrolyte balance. The addition of 6% horse blood enhances the medium's nutritional value. Hemolysis and the type of hemolysis observed in horse red blood cells assist in identifying certain Gram-positive cocci. The inclusion of colistin and nalidixic acid inhibits the growth of Gram-negative bacteria.

COMPOSITION	g/litre
COLUMBIA C.N.A. W/ 6% HORSE BLOOD	
Peptone Mixture	20,0
Beef Extract	3,0
Sodium chloride	5.0
Corn Starch	1,0
Bacteriological Agar	15.0
Horse Blood	60ml
Colistin	10mg
Nalidixic acid	15mg

Appearance: Red-burgundy, non-transparent due to the addition of blood.

Final pH 7.3 ± 0.2 at 25°C .

MAC CONKEY AGAR is a nutrient medium for the isolation of Enterobacteriaceae. Gram-negative enterobacteria that ferment lactose produces red or pink colonies. Bile salts No. 3 and crystal violet inhibit the growth of Gram-positive cocci. The neutral red indicator changes color with lactose fermentation.

COMPOSITION	g/litre
MAC CONKEY AGAR	
Peptone	20.0
Lactose	10.0
Bile Salts No. 3	1.5
Sodium chloride	5.0
Neutral red	0.03
Crystal violet	0.001
Agar No. 2	15.0

Appearance: pink-purple clear

Final pH 7.1 ± 0.2 at 25°C

PRECAUTIONS

The COLUMBIA C.N.A. W/ 6% HORSE BLOOD – MAC CONKEY AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The plates should be stored at $2-12^{\circ}\text{C}$ in their packaging until use. Prolonged storage below 2°C can create moisture in the medium, increasing the risk of contamination. Freezing, even momentarily, destroys the medium. Also, avoid excessive heating. The plates can be used until the expiration date indicated on the label. For transport, stability studies have shown that the plates can remain at $18-25^{\circ}\text{C}$ for 7 days or at $30-40^{\circ}\text{C}$ for 48 hours without affecting product performance.

USAGE

Inoculate the plates using the technique for single colony isolation. Incubate at 35-37°C for 24 hours in aerobic conditions. After incubation, the appearance of the colonies is as follows:

On COLUMBIA C.N.A. W/ 6% HORSE BLOOD, the hemolytic streptococcus gives clear or matte grey colonies surrounded by a beta-hemolysis zone. Pneumococcus appears as flat, clear colonies surrounded by a small green alpha-hemolysis zone. Staphylococcus appears as matte white to golden-yellow colonies with or without a beta-hemolysis zone. Gram-negative bacteria do not grow.

On MAC CONKEY AGAR, gram-positive cocci do not grow. Gram-negative bacteria that ferment lactose produce pink to red colonies due to acid production. Additionally, Escherichia coli produces pink or red colonies surrounded by a cloudy zone. Gram-negative bacteria that do not ferment lactose produce colorless colonies. Proteus does not swarm.

LIMITATIONS OF THE METHOD

In COLUMBIA C.N.A. W/ 6% HORSE BLOOD, you may observe the growth of Gram (-) bacteria that are normally sensitive to the antibiotics in the medium. This can occur when the number of bacteria per ml of sample (CFU) is $\geq 10^5$. Also, some strains may show resistance to the antibiotics in the medium.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

COLUMBIA C.N.A. W/ 6% HORSE BLOOD

Germ	ATCC	Growth	Hemolysis
<i>Escherichia coli</i>	25922	Is suspended	-
<i>Staphylococcus aureus</i>	25923	Good	Beta
<i>Streptococcus pneumoniae</i>	6303	Good	Alpha
<i>Streptococcus pyogenes</i>	19615	Excellent	Alpha

Streptococcus pyogenes ATCC 19615

GENERAL CHARACTERISTICS OF QUALITY CONTROL

MAC CONKEY AGAR

Germ	ATCC	Growth	Colonies
<i>C</i>	14028	Good	Colorless - translucent
<i>P. mirabilis</i>	12453	Good	Colorless, it does not crawl
<i>Escherichia coli</i>	25922	Good	Pink or red are surrounded by a blurred zone
<i>E. faecalis</i>	29212	Is suspended	

Escherichia coli

P. mirabilis

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

COLUMBIA W/ 6% HORSE BLOOD – MAC CONKEY AGAR - GR/CA01/GRM5/O/73 - 

ITEM	CODE	PACKAGE	STORAGE	SELF LIFE
Petri dish halved 9cm	020261	10 pieces	2 – 12 °C	2 μήνες

It is produced in Greece from the Bioprep company in accordance with the requirement of the European guidance 98/79/EK. Code by EDMA 14 01 04 90. The Bioprep company has been certified according to the standards EN ISO 9001:2015 / ΕΛΟΤ EN ISO 13485:2016 DY8d/1348/2004.

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45:502-504. Ruoff, k., I. 1995 Streptococcus.

American Public Health Association (1950). Diagnostic Procedures and Reagents. 3rd edn. A.P.H.A., New York.

American Public Health Association (1946). Standard Methods for the examination of Water and Sewage. 9th edn. A.P.H.A., New York.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr