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TECHNICAL DATA SHEET

PRODUCT: **COLUMBIA SHEEP BLOOD AGAR** 
 REFERENCE: **010250 - 050250**

Date 1st Edition:
7th 2009
 Date 4rth Revision:
6th 2024

DESCRIPTION

A medium recommended for the cultivation of a wide variety of pathogenic microorganisms. This medium is prepared with the addition of 5% sheep blood.

PRINCIPLE OF THE METHOD

COLUMBIA SHEEP BLOOD AGAR contains ingredients of high nutritional value that provide vitamins, carbohydrates, and other organic components. Sodium chloride supplies the necessary minerals and maintains osmotic and electrolyte balance. The addition of 5% sheep blood enhances the nutritive properties of the medium. Furthermore, the hemolysis and the type of hemolysis of sheep red blood cells aid in the identification of certain gram-positive cocci.

COMPOSITION	g/litre
Columbia peptone mixture	25.1
Soluble starch	1.0
Sodium chloride	5.0
Agar	12.0
Sheep Blood Defibrinated	50ml

Appearance: Red to burgundy, non-transparent, due to the addition of blood.

Final pH 7.3 ± 0.2 at 25 °C.

PRECAUTIONS

COLUMBIA SHEEP BLOOD AGAR is an in vitro diagnostic laboratory material and should only be handled by trained laboratory personnel. This material contains peptones and extracts of animal origin. Certificates regarding the source and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious and handled using standard safety precautions (avoid ingestion and inhalation). Always handle plates wearing gloves and within a Class II laminar flow cabinet to prevent contamination, especially by saprophytic fungi. Do not use plates that show signs of microbial contamination. The material should not be used after the expiration date. In case of skin contact, wash immediately with plenty of water and soap. Positive samples must be disposed of according to the appropriate hygiene regulations for infectious materials.

STORAGE AND TRANSPORT CONDITIONS

The plates must be stored at **2 – 12 °C** in their packaging until the time of use. Prolonged storage at temperatures below **2 °C** creates excessive moisture in the medium, increasing the risk of contamination. Freezing, even momentarily, destroys the medium. Excessive heating should also be avoided. The plates can be used until the expiration date indicated on the label. Stability studies have shown that for transportation, the plates can remain at **18 – 25 °C** for up to **4 days** or at **25 – 40 °C** for up to **48 hours** without affecting product performance.

INSTRUCTIONS FOR USE

Inoculate the plates using the streak plate technique to obtain isolated colonies.

Incubate at 35 – 37 °C for 24 hours under aerobic conditions.

INTERPRETATION OF RESULTS

Streptococcus pyogenes forms small, grayish-white colonies with beta-hemolysis.

Streptococcus pneumoniae forms small, flat colonies with alpha-hemolysis.

Staphylococcus aureus forms relatively large colonies (1–3 mm after 24 hours of incubation and 3–8 mm if incubation is extended up to 5 days). Their color ranges from cream-yellow to orange, depending on incubation time. It produces alpha or beta hemolysis.

LIMITATIONS OF THE METHOD

Final identification must be performed using biochemical and serological tests (e.g., Latex Agglutination Test), which may be applied directly from suspect colonies.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	ATCC	Growth	Hemolysis
<i>Escherichia coli</i>	25922	Good	-
<i>Pseudomonas aeruginosa</i>	27853	Good	-
<i>Staphylococcus aureus</i>	25923	Good	B
<i>Streptococcus pneumoniae</i>	6303	Good	A
<i>Streptococcus pyogenes</i>	19615	Excellent	B

Streptococcus pyogenes ATCC 19615

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

COLUMBIA SHEEP BLOOD AGAR - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Plate 90mm	010250	10 pieces	2 – 12 °C	2 months
Plate 60mm	050250	10 pieces	2 – 12 °C	2 months

Produced in Greece by Bioprep company in accordance with Regulation (EU) 2017/746.

Basic UDI-DI: 5212037714010490X8. EDMA: (14 01 04 90) Other Prepared Media in Plates.

Bioprep is certified according to the standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 DY86/1348/2004.

REFERENCES

Ellener, P.C., C.J. Stoessel, E. Drakeford, and F. Vassi 1966. A new culture medium for medical bacteriology. Am J. Clin Pathol. 45:502-504. Ruoff, k., I. 1995 Streptococcus.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr