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COMPARISON OF PETRIFILM AND CHROMAGAR ECC FOR ISOLATION OF E. COLI FROM CHICKEN

2005. Submitted to: International Association for Food Protection Publication Date: July 15, 2005 Citation: Bailey, J.S., Cray, P.J., Berrang, M.E., Plumlee, J. 2005. Comparison of petrifilm and chromagar ecc for isolation of e. coli from chicken [abstract]. International Association for...

Trial report : "Comparison between CHROMagar ECC method for the detection and enumeration of Escherichia coli and coliforms and FIL-IDF N°73B: 1998 method for the enumeration of coliforms

2001, National Institute of Industrial Technology.

Comparison and recovery of Escherichia coli and thermotolerant coliforms in water with a chromogenic medium incubated at 41 and 44.5°C.

1999. Alonso J. L. et al. .

1999. Applied and Environmental Microbiology, 65 : 3746-3749.

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E-mail: bioprep1@otenet.gr www.bioprepare.gr

TECHNICAL DATA SHEET

PRODUCT: **CONTACT PLATE CHROMagar™ ECC** 
(CP ECC).

REFERENCE: **20029**

Date 1st Edition:
7th 2009
Date 4rth Revision:
6th 2024

DESCRIPTION

The Contact Plate CHROMagar™ ECC is a specialized contact dish used for the detection and identification of *E. coli* and coliforms. It is recommended for microbiological monitoring of surfaces.

METHOD PRINCIPLE

In a contact plate, the agar medium can come into direct contact with the surface being tested. With the Contact Plate CHROMagar™ ECC, the user can monitor workbenches, walls, product conveyor belts in industrial settings, production machinery, tools, and more. CHROMagar™ ECC is a nutrient medium that visually identifies *E. coli* colonies by their blue color, making detection and quantification of this important hygiene indicator as straightforward as possible. Additionally, CHROMagar™ ECC reveals other *Enterobacteriaceae* as red colonies.

COMPOSITION	g/litre
Chromogenic mix	4.8
Peptone, yeast and extract	8.0
Sodium chloride	5,0
Agar	15,0

Appearance: Light beige and clear

Final pH 7.2 ± 0.2 . at 25 °C.

PRECAUTIONS

Contact Plate CHROMagar™ ECC is an in vitro diagnostic laboratory material and should only be handled by trained laboratory personnel. This material contains peptones and extracts of animal origin. Certificates regarding the source and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious and handled using standard safety precautions (avoid ingestion and inhalation). Always handle plates wearing gloves and within a Class II laminar flow cabinet to prevent contamination, especially by saprophytic fungi. Do not use plates that show signs of microbial contamination. The material should not be used after the expiration date. In case of skin contact, wash immediately with plenty of water and soap. Positive samples must be disposed of according to the appropriate hygiene regulations for infectious materials.

STORAGE AND TRANSPORT CONDITIONS

The plates must be stored at a temperature between **2–12 °C** in their packaging until the time of use. Prolonged storage below **2 °C** generates excess moisture in the medium, which increases the risk of contamination. Freezing, even momentarily, destroys the medium. Excessive heating should also be avoided. The plates can be used up to the expiration date indicated on the label. For transport, stability studies have shown that the plates can remain at **6–25 °C** for **4 days** or at **25–40 °C** for **24 hours** without affecting the product's performance.

INSTRUCTIONS FOR USE

Place the entire surface of the medium onto the surface you wish to test, which must be completely dry. Press gently for 8–12 seconds. Close the plate and place it in the incubator.

Surface	Incubation Time	Temperature
Environmental surfaces	72 ώρες	35 ± 2 °C / 24 h 25 ± 2 °C / 48 h
Human monitoring	72 ώρες	35 ± 2 °C / 24 h 25 ± 2 °C / 48 h

RESULTS

After incubation, inspect for potential growth of *E. coli* and coliforms and count the colonies, considering that the contact plate surface area is 25 cm². Colonies positive for both β-D-galactosidase and β-D-glucuronidase (dark blue to violet color) are counted as *E. coli*. Colonies positive for β-D-galactosidase only (pink to red color) are counted as coliforms. Gram-positive bacteria do not grow on this medium.

Microorganisms	Interpretation
One colony of <i>E. coli</i> or coliforms per 25 cm²	The surface must be cleaned and disinfected.


LIMITATIONS OF THE METHOD

To avoid false-positive results caused by oxidase-positive bacteria, such as *Pseudomonas spp.*, colonies that appear slightly bluish or colorless (potential false positives) should be tested for oxidase production using Oxidase Test Strips (MID-61G). The total count of coliform bacteria is the sum of oxidase-negative, β -D-galactosidase-positive colonies (pink to red) and all colonies that exhibit a dark blue to violet color (*E. coli*). Dark blue to violet colonies should be tested for indole production using Kovac's reagent (to confirm *E. coli*). It is recommended to perform biochemical testing on colonies from a pure culture for full identification.

GENERAL QUALITY CONTROL CHARACTERISTICS

Microorganism	Growth/Colony Color
<i>Escherichia coli</i> ATCC® 25922	Blue
<i>Citrobacter freundii</i> ATCC®8090	Purple-Red
<i>Enterobacter cloacae</i> ATCC® 23355	Purple-Red
<i>Klebsiella pneumoniae</i> ATCC® 13883	Purple-Red
<i>Staphylococcus aureus</i> ATCC® 25923	Inhibited
<i>Enterococcus faecalis</i> ATCC® 29212	Inhibited

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CONTACT PLATE CHROMagar™ ECC - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Petri dish 65mm	20029	10 pieces	2 – 12°C	3 months

This product is manufactured in Greece by the company Bioprep, in accordance with Regulation (EU) 2017/746. The Basic UDI-DI is 5212037714010402WH. According to the EDMA classification, it falls under code (14 01 04 02) for Chromogenic ID Media (Plates). Bioprep is certified in compliance with the standards EN ISO 9001:2015, ELOT EN ISO 13485:2016, and Ministerial Decision DY86/1348/2004.

REFERENCES

COMPARISON OF PETRIFILM AND CHROMAGAR ECC FOR ISOLATION OF E. COLI FROM CHICKEN

2005. Submitted to: International Association for Food Protection Publication Date: July 15, 2005 Citation: Bailey, J.S., Cray, P.J., Berrang, M.E., Plumblee, J. 2005. Comparison of petrifilm and chromagar ecc for isolation of e. coli from chicken [abstract].International Association for ...

Trial report : "Comparison between CHROMagar ECC method for the detection and enumeration of Escherichia coli and coliforms and FIL-IDF N°73B: 1998 method for the enumeration of coliforms

2001, National Institute of Industrial Technology.

Comparison and recovery of Escherichia coli and thermotolerant coliforms in water with a chromogenic medium incubated at 41 and 44.5°C.

1999. Alonso J. L. et al. .

1999. Applied and Environmental Microbiology, 65 : 3746-3749.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr www.bioprepare.gr