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Sandys, G. H. 1960. A new method of preventing swarming of *Proteus* spp. with a description of a new medium suitable for use in routine laboratory practice. J. Med. Lab. Technol. 17:224.

Mackey, J. P., and G. H. Sandys. 1965. Laboratory diagnosis of infections of the urinary tract in general practice by means of a dip-inoculum transport medium. Br. Med. J. 2:1286.

Mackey, J. P., and G. H. Sandys. 1966. Diagnostic of urinary tract infections. Br. Med. J. 1:1173.

MacFaddin, J. D. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. 1 Williams & Wilkins, Baltimore, MD.

Baron, E. J., L. R. Peterson, and S. M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis, MO **Isenberg, H. D. (ed.).** 1992. Clinical microbiology procedures handbook

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E-mail: bioprep1@otenet.gr www.bioprepere.gr

TECHNICAL DATA SHEET

PRODUCT: **PHENYLALANINE AGAR (PPA)** 
 REFERENCE: **080086**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

PHENYLALANINE AGAR (PPA) is used to differentiate bacteria capable of deaminating phenylalanine.

PRINCIPLE OF THE METHOD

Certain microorganisms, such as *Proteus mirabilis*, can deaminate phenylalanine, producing pyruvic acid. When ferric chloride is added after incubation, pyruvic acid reacts, forming a green complex.

Positive reaction: Blue-green ring formation

Negative reaction: No color change in the medium

FORMULA	g/litre
Yeast Extract	3.0
Dipotassium Phosphate	1.0
Sodium Chloride	5.0
DL-Phenylalanine	2.0
Agar	12.0

Appearance: Thick-colored clear.

Final pH 7.3 ± 0.2 at 25 °C.

PRECAUTIONS

PHENYLALANINE AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The tubes should be stored at **2 – 12°C** in their original packaging until use. Freezing, even momentarily, destroys the medium. Avoid excessive heating. The tubes can be used until the expiration date printed on the label. For transport, stability studies have shown that the tubes can remain at **6 - 25°C for 5 days**, or **25 - 40°C for 48 hours** without affecting product performance.

USAGE INSTRUCTIONS

Allow the PPA medium to reach room temperature before inoculation. Inoculate the tube with a single bacterial colony, streaking across the slanted surface. Incubate aerobically at 35 – 37°C for 18-24 hours. After incubation, add 4-5 drops of 10% ferric chloride solution directly onto the slanted surface. Gently shake the tube and observe color development within 1-5 minutes.

INTERPRETATION OF RESULTS

Positive reaction: Formation of a light to dark green color within 1-5 minutes after adding ferric chloride.

Negative reaction: Yellow color due to the ferric chloride reagent itself.

LIMITATIONS OF THE METHOD

The green color reaction of a positive test fades quickly. Test results should be interpreted within 5 minutes after adding ferric chloride, as delayed readings may lead to false negatives. Gentle shaking of the tube containing ferric chloride will help dislodge surface colonies and produce a faster, more intense color reaction.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

MICROORGANISMS	REACTION
<i>Proteus mirabilis</i> ATCC 7002	Green (+)
<i>Proteus vulgaris</i> ATCC 13315	Green (+)
<i>Enterobacter aerogenes</i> ATCC 13048	Yellow (-)
<i>Escherichia coli</i> ATCC 25922	Yellow (-)

Proteus mirabilis GREEN (+)

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

PHENYLALANINE AGAR (PPA) -

PRODUCT	CODE	PACKING	STORAGE	SELF LIFE
Tube 3ml	080086	20 pieces	2 – 12 °C	8 months

Produced in Greece by Bioprepure in accordance with the requirements of the European Directive 98/79/EK. FEK B2198/2-10-2009. CODE EDMA 14 01 02 01. Bioprepure company has been certified according to standards EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004

BIBLIOGRAPHY

- Sandys, G. H. 1960. A new method of preventing swarming of *Proteus* spp. with a description of a new medium suitable for use in routine laboratory practice. *J. Med. Lab. Technol.* 17:224.
- Mackey, J. P., and G. H. Sandys. 1965. Laboratory diagnosis of infections of the urinary tract in general practice by means of a dip-inoculum transport medium. *Br. Med. J.* 2:1286.
- Mackey, J. P., and G. H. Sandys. 1966. Diagnostic of urinary tract infections. *Br. Med. J.* 1:1173.
- MacFaddin, J. D. 1985. *Media for isolation-cultivation-identification-maintenance of medical bacteria*, vol. 1 Williams & Wilkins, Baltimore, MD.
- Baron, E. J., L. R. Peterson, and S. M. Finegold. 1994. *Bailey & Scott's diagnostic microbiology*, 9th ed. Mosby-Yearbook, Inc., St. Louis, MO
- Isenberg, H. D. (ed.). 1992. *Clinical microbiology procedures handbook*.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepure.gr