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New plate medium for facilitated differentiation of *Salmonella* spp. from *Proteus* spp. and other enteric bacteria. 1990. Rambach A. 1990, Applied and Environmental Microbiology, 56 : 301-303.

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TECHNICAL DATA SHEET

PRODUCT: **RAMBACH™ AGAR** 
 REFERENCE: **010092– 050092**

Date 1st Edition:
7th 2009
 Date 4th Revision:
6th 2025

DESCRIPTION

RAMBACH™ AGAR is a chromogenic medium for the selective differentiation of *Salmonella* spp., including *Salmonella Typhi*, from other Gram-negative bacteria in clinical samples and food.

INTRODUCTION

Infections caused by *Salmonella* spp., including *Salmonella Typhi*, remain a significant global health issue. In the U.S., *Salmonella* has an incidence rate of 16.47 cases per 100,000 population (CDC estimate, 2010). In Europe, it is reported as the leading cause of foodborne outbreaks (EFSA/ECDC report 2011, data from 2009). In developing countries, *Salmonella Typhi* and *Salmonella Paratyphi* are commonly encountered, with an estimated 17 million cases annually (EFSA report, 2007). Additionally, according to a recent WHO report, *Salmonella* infections are responsible for 2 million deaths annually due to diarrhea. *Salmonella* is the second most commonly reported zoonotic infection in humans (EFSA/ECDC report 2011, data from 2009).

PRINCIPLE OF THE METHOD

Peptones provide the essential nutrients for bacterial growth, serving as a primary source of organic nitrogen, especially amino acids and polypeptides. They are also rich in vitamins, minerals, and trace elements. Yeast extract serves as a vitamin source, particularly from the B-complex group. Sodium deoxycholate inhibits the growth of Gram-positive organisms. The chromogenic substrates are broken down by specific microbial enzymes, releasing insoluble colored compounds that bind to microbial cell membranes. *Salmonella* utilizes only one of the chromogenic substrates and forms red colonies. If none of the substrates are utilized, colonies appear colorless or white. Sodium chloride supplies electrolytes for transport and maintains osmotic balance. Agar provides trace elements and electrolytes and serves as the solidifying agent.

COMPOSITION	g/litre
Peptone	5.0
Yeast extract	2.0
Meat extract	1.0
Sodium chloride	5.0
Sodium deoxycholate	1.0
Chromogenic mix	1.5
Propylene glycol	10.5
Agar	15.0

Appearance: Bright pink, opaque

Final pH 7.3 ± 0.2 at 25 °C.

PRECAUTIONS

To RAMBACH™ AGAR is in vitro diagnostic laboratory material and should only be handled by trained laboratory personnel. This material contains peptones and extracts of animal origin. Certificates regarding the source and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious and handled using standard safety precautions (avoid ingestion and inhalation). Always handle plates wearing gloves and within a Class II laminar flow cabinet to prevent contamination, especially by saprophytic fungi. Do not use plates that show signs of microbial contamination. The material should not be used after the expiration date. In case of skin contact, wash immediately with plenty of water and soap. Positive samples must be disposed of according to the appropriate hygiene regulations for infectious materials.

STORAGE AND TRANSPORT CONDITIONS

Plates should be stored at **2–12 °C** in their original packaging until use. Prolonged storage below **2 °C** may lead to excess moisture within the medium, increasing the risk of contamination. Freezing, even momentarily, will destroy the medium. Excessive heating should also be avoided. Plates can be used up to the expiry date stated on the label. If the airtight packaging is accidentally opened, the plate may be stored in the refrigerator for **5–7 days**, provided it is sealed with parafilm or placed in a sealed bag. According to our stability studies, the plates can be transported at **6–25 °C** for up to **4 days**, or at **25–40 °C** for up to **24 hours**, without affecting product performance.

INSTRUCTIONS FOR USE

Clinical samples:

Place the plates in an incubator at 35–37 °C for 30–45 minutes to dry. Inoculate the sample at one end of the plate, then streak it in parallel lines using a loop to obtain isolated colonies. Incubate at 35–37 °C in an aerobic atmosphere for 18–24 hours. After 18–24 hours, examine the plates for the following color changes.

Industrial samples:

Add 25 g of the sample to 225 ml of BPW or a proportionally smaller amount. Mix thoroughly and incubate at $37 \text{ °C} \pm 1 \text{ °C}$ for 18 ± 2 hours. Then, inoculate 0.1 ml of BPW into 10 ml of Rappaport-Vassiliadis Broth (070093) and incubate at $41.5 \pm 1 \text{ °C}$ for 24 ± 3 hours. Also, inoculate 1 ml of BPW into 10 ml of MKTTn Broth (070526) and incubate at $37 \pm 1 \text{ °C}$ for 24 ± 3 hours. Then, inoculate 10 µl from each of the RVS and MKTTn broths onto two RAMBACH™ AGAR plates and incubate at $37 \pm 1 \text{ °C}$ for 24 ± 3 hours.

READING AND INTERPRETATION


After 24 hours of incubation, examine the plates for potential *Salmonella* growth. Pathogenic *Salmonella* appears as red colonies 2–3 mm in diameter. Most *Enterobacteriaceae* form blue-violet colonies. *Proteus* species typically appear as colorless colonies.

LIMITATIONS OF THE METHOD

Various studies have shown that 97% to 99% of *Salmonella* strains produce typical red colonies, except *S. Paratyphi A* and *S. Typhi*, which may form colorless colonies or none, and lactose-positive *Salmonella*, which may form blue colonies. Final identification must be confirmed using biochemical and serological tests (e.g., Microgen *Salmonella* Latex agglutination test, code F42), which can be performed directly on suspected *Salmonella* colonies. Some *Pseudomonas* strains may also form red colonies. This can be ruled out with an oxidase test (Oxidase Test Strips – Code MID-61G), which can be performed directly on the suspect colony.

GENERAL QUALITY CONTROL CHARACTERISTICS

Microorganism	Colony Characteristics
<i>Salmonella typhimurium</i> ATCC® 14028	Red
<i>Salmonella abaeetuba</i> ATCC® 35640	Red
<i>Salmonella enteritidis</i> ATCC® 13076	Red
<i>Escherichia coli</i> ATCC® 25922	Small, metallic blue
<i>Citrobacter freundii</i> ATCC® 8090	Purple-blue
<i>Staphylococcus aureus</i> ATCC® 25923	Inhibited

DISPOSAL OF MATERIALS IN WASTE

Materials that show no growth can be considered non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of in accordance with guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must manage and dispose of (or assign the management and disposal of) it according to the applicable regulations.

SPECIFICATIONS

RAMBACH™ AGAR - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Plate 9 cm	010092	10 τεμάχια	2 – 12 °C	3 months
Plate 6 cm	050092	10 τεμάχια	2 – 12 °C	3 months

Manufactured in Greece by Bioprepare in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 5212037714010402WH EDMA: (14 01 04 02) Chromogenic ID Media (Plates) The company Bioprepare is certified according to the following standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 / Ministerial Decision ΔΥ86/1348/2004.

REFERENCES

Intérêt du milieu Rambach pour l'isolement et l'identification rapide des *Salmonella*. Laudat P. et al. 1992. 1992, Revue française des laboratoires, N°233, 45-47.
Use of Rambach propylene glycol containing agar for identification of *Salmonella* spp. 1991. Gruenewald R. et al. 1991, Journal of Clinical Microbiology, 29 : 2354-2356.
New plate medium for facilitated differentiation of *Salmonella* spp. from *Proteus* spp. and other enteric bacteria. 1990. Rambach A. 1990, Applied and Environmental Microbiology, 56 : 301-303.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepare.gr