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- Isenberg, H.D., Kominos, S., and Siegel, M. (1969). Isolation of salmonellae and shigellae from an artificial mixture of fecal bacteria. *Appl. Microbiol.*, 18: 4, 656-659.
- Leifson, E. (1935). New culture media based on sodium desoxycholate for the isolation of intestinal pathogens and for the enumeration of colon bacilli in milk and water. *J. Pathol. Bacteriol.*, 40: 581-589.
- American Public Health Association (1950). *Diagnostic Procedures and Reagents*. 3rd edn. A.P.H.A., New York.
- American Public Health Association (1946). *Standard Methods for the examination of Water and Sewage*. 9th edn. A.P.H.A., New York.

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E-mail: bioprep1@otenet.gr www.bioprepare.gr

DESCRIPTION

The divided plate S.S. AGAR – MAC CONKEY AGAR is used for the isolation and identification of *Salmonella* and *Shigella* species (S.S. AGAR), as well as for the isolation and identification of Enterobacteriaceae (MAC CONKEY AGAR).

PRINCIPLE OF THE METHOD

In S.S. Agar, beef extract and balanced peptone No. 1 supply nitrogen, vitamins, minerals, and amino acids essential for bacterial growth. Lactose is the fermentable carbohydrate providing carbon and energy. The combination of Bile Salts No. 3, sodium citrate, and brilliant green inhibits Gram-positive bacteria, most coliforms, and some *Proteus* spp., while allowing the growth of *Salmonella* spp. Neutral red is the pH indicator. Sodium thiosulphate and ferric citrate enable detection of hydrogen sulfide (H₂S) production by bacteria such as *Proteus* and certain *Salmonella* strains, which form colonies with black centers and clear halos. Non-lactose fermenting bacteria (presumed pathogens) form clear, transparent, or colorless colonies, while coliforms are sufficiently inhibited and form small colonies ranging from pink to red. This medium is highly selective and not recommended for the primary isolation of *Shigella*, as some *Shigella* strains may fail to grow.

COMPOSITION S.S. AGAR	g/litre
Beef Extract	5.0
Balanced Peptone No. 1	5.0
Lactose	10.0
Bile Salts No. 3	8.5
Sodium citrate	8.5
Sodium thiosulphate	8.5
Ferric citrate	1.0
Brilliant Green	0.00033
Neutral Red	0.025
Agar No. 2	13.5

Appearance: Clear red-pink

Final pH 7.4 ± 0.2. at 25 °C

MacConkey Agar is a nutrient medium used for the isolation of *Enterobacteriaceae*. Gram-negative enteric bacteria that ferment lactose form red or pink colonies. Bile Salts No. 3 and crystal violet inhibit the growth of Gram-positive cocci. Neutral red acts as a pH indicator and changes color in response to lactose fermentation.

COMPOSITION MAC CONKEY AGAR	g/litre
Peptone	20.0
Lactose	10.0
Bile Salts No. 3	1.5
Sodium chloride	5.0
Neutral red	0.03
Crystal violet	0.001
Agar No. 2	15.0

Appearance: Clear pink to purple

Final pH 7.1 ± 0.2 at 25°C.

PRECAUTIONS

S.S. AGAR – MAC CONKEY AGAR is in vitro diagnostic laboratory material and should only be handled by trained laboratory personnel. This material contains peptones and extracts of animal origin. Certificates regarding the source and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious and handled using standard safety precautions (avoid ingestion and inhalation). Always handle plates wearing gloves and within a Class II laminar flow cabinet to prevent contamination, especially by saprophytic fungi. Do not use plates that show signs of microbial contamination. The material should not be used after the expiration date. In case of skin contact, wash immediately with plenty of water and soap. Positive samples must be disposed of according to the appropriate hygiene regulations for infectious materials.

STORAGE CONDITIONS

Plates should be stored at **2–12 °C** in their original packaging until the time of use. Prolonged storage at temperatures below **2 °C** may lead to excess moisture inside the medium, increasing the risk of contamination. Freezing, even momentarily, will destroy the medium. Excessive heating should also be avoided. Plates may be used up to the expiration date indicated on the label. If the airtight packaging is accidentally opened, the plate can be stored in the refrigerator for **5–7 days**, provided it is resealed with parafilm or placed in a sealed bag. According to our stability studies, the plates can be transported at **6–25 °C** for up to **4 days** or at **25–40 °C** for up to **48 hours** without affecting the product's performance.

INSTRUCTIONS FOR USE

Place the plates in an incubator at 35–37 °C for 30–45 minutes. Inoculate the sample as soon as possible after collection and streak using successive dilutions to obtain isolated colonies. Incubate under aerobic conditions at 35–37 °C for 24 hours. Refer to the colony color chart for interpretation.

LIMITATIONS OF THE METHOD

Proteus is a bacterium that may form black colonies on S.S. Agar, like *Salmonella*, due to its ability to produce hydrogen sulfide (H₂S) by breaking down ferric ammonium citrate.

GENERAL QUALITY CONTROL CHARACTERISTICS

S.S. AGAR

Microorganism	ATCC	Colony Morphology and Color
<i>Escherichia coli</i>	25922	Partial inhibition; red colonies
<i>Enterococcus faecalis</i>	29212	Partial inhibition
<i>Salmonella typhimurium</i>	14028	Colorless colonies with black center
<i>Shigella flexneri</i>	12022	Colorless colonies

S.S. AGAR: *Salmonella typhimurium* ATCC 14028

GENERAL QUALITY CONTROL CHARACTERISTICS

MAC CONKEY AGAR

Microorganism	ATCC	Growth	Colony Morphology
<i>S. typhimurium</i>	14028	Καλή	Colorless to semi-transparent colonies
<i>P. mirabilis</i>	12453	Καλή	Colorless, non-swarming
<i>Escherichia coli</i>	25922	Καλή	Pink or red colonies with a hazy zone
<i>E. faecalis</i>	29212	Αναστέλλεται	

Escherichia coli

P. mirabilis

DISPOSAL OF MATERIALS IN WASTE

Materials that show no growth can be considered non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of in accordance with guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must manage and dispose of (or assign the management and disposal of) it according to the applicable regulations.

SPECIFICATIONS

S.S. AGAR – MAC CONKEY AGAR - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Divided Plate 9 cm	020558	10 pieces	2 – 12 °C	3 months

Manufactured in Greece by Bioprep in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 5212037714010490X8 EDMA: (14 01 04 90) Other Prepared Media in Plates The company Bioprep is certified according to the following standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 / Ministerial Decision ΔΥ86/1348/2004

REFERENCES

Isenberg, H.D., Kominos, S., and Siegel, M. (1969). Isolation of salmonellae and shigellae from an artificial mixture of fecal bacteria. Appl. Microbiol., 18: 4, 656-659.

Leifson, E. (1935). New culture media based on sodium desoxycholate for the isolation of intestinal pathogens and for the enumeration of colon bacilli in milk and water. J. Pathol. Bacteriol., 40: 581-589.

American Public Health Association (1950). Diagnostic Procedures and Reagents. 3rd edn. A.P.H.A., New York.

American Public Health Association (1946). Standard Methods for the examination of Water and Sewage. 9th edn. A.P.H.A., New York.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr www.bioprepare.gr