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Leifson, E. (1935). New culture media based on sodium desoxycholate for the isolation of intestinal pathogens and for the enumeration of colon bacilli in milk and water. J. Pathol. Bacteriol., 40: 581-589.

Taylor, W.I., Harris, B. (1965) Isolation of shigellae. II. Comparison of plating media and enrichment broths. Am. J. Clin. Pathol. 44: 4, 476-479.

Wallace, J.S. and K. Jones. 1996. The Use of Selective and Differential Agars in the Isolation of E.coli O157 from dairy herds. J. Appl. Bacteriol. 81: 663-668

Wallace, J.S. et al. 1997. Isolation of Vero cytotoxin-producing

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## DESCRIPTION

The divided plate S.S. AGAR – RAMBACH AGAR is used for the isolation and identification of *Salmonella* and *Shigella* species (S.S. AGAR) and for the selective isolation and identification of *Salmonella* spp. (RAMBACH AGAR).

## PRINCIPLE OF THE METHOD

In S.S. Agar, beef extract and balanced peptone No. 1 supply nitrogen, vitamins, minerals, and amino acids necessary for bacterial growth. Lactose is the fermentable carbohydrate that provides carbon and energy. A combination of Bile Salts No. 3, sodium citrate, and brilliant green inhibits Gram-positive bacteria, most coliforms, and some *Proteus* spp., while allowing *Salmonella* spp. to grow. Neutral red acts as a pH indicator. Sodium thiosulphate and ferric citrate enable the detection of hydrogen sulfide (H<sub>2</sub>S) production by bacteria such as *Proteus* and some *Salmonella* strains, which appear as colonies with black centers and a clear halo. Non-lactose fermenters (presumed pathogens) form clear, transparent, or colorless colonies, while inhibited coliforms form small colonies ranging in color from pink to red. This medium is highly selective and not recommended for the primary isolation of *Shigella*, as some *Shigella* strains may not grow.

COMPOSITION S.S. AGAR	g/litre
Beef Extract	5.0
Balanced Peptone No. 1	5.0
Lactose	10.0
Bile Salts No. 3	8.5
Sodium citrate	8.5
Sodium thiosulphate	8.5
Ferric citrate	1.0
Brilliant Green	0.00033
Neutral Red	0.025
Agar No. 2	13.5

Appearance: Clear red-pink

Final pH 7.4 ± 0.2. at 25 °C

RAMBACH Agar is a chromogenic, selective nutrient medium that supports the growth of Enterobacteriaceae, with its primary use being the differentiation of *Salmonella* spp. from other enterobacteria. *Salmonella* colonies appear distinctly red—making them easily distinguishable from other species—due to a unique mechanism in this medium. Specifically, *Salmonella* metabolizes propylene glycol, leading to an immediate drop in pH. In the presence of a pH indicator, this results in red-colored colonies. To differentiate other enterobacteria from one another, the medium includes a chromogenic mix. When β-galactosidase breaks down this mix, it causes a pH shift that produces colonies in various distinct colors. Gram-positive cocci are inhibited due to the presence of sodium deoxycholate.

COMPOSITION RAMBACH AGAR	g/litre
Peptone	5
Yeast extract	2
Meat extract	1
Sodium chloride	5
Sodium deoxycholate	1
Chromogenic mix	1,5
Propylene glycol	10,5
Agar	15

Appearance: Intense pink, opaque

Final pH 7.3 ± 0.2 °C at 25 °C

## PRECAUTIONS

S.S. AGAR – RAMBACH AGAR is in vitro diagnostic laboratory material and should only be handled by trained laboratory personnel. This material contains peptones and extracts of animal origin. Certificates regarding the source and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious and handled using standard safety precautions (avoid ingestion and inhalation). Always handle plates wearing gloves and within a Class II laminar flow cabinet to prevent contamination, especially by saprophytic fungi. Do not use plates that show signs of microbial contamination. The material should not be used after the expiration date. In case of skin contact, wash immediately with plenty of water and soap. Positive samples must be disposed of according to the appropriate hygiene regulations for infectious materials.

## STORAGE CONDITIONS

Plates should be stored at 2–12 °C in their original packaging until the time of use. Prolonged storage below 2 °C may cause excessive moisture to form in the medium, increasing the risk of contamination. Freezing, even momentarily, will destroy the medium. Excessive heating should also be avoided. Plates may be used until the expiration date indicated on the label. If the airtight packaging is accidentally opened, the plate may be stored in the refrigerator for 5–7 days, provided it is sealed with parafilm or placed in a sealed bag. According to our stability studies, the plates can be transported at 6–25 °C for up to 4 days or at 25–40 °C for up to 48 hours without affecting product performance.

## INSTRUCTIONS FOR USE

Place the plates in the incubator at 35–37 °C for 30–45 minutes. Inoculate the sample as soon as possible after collection and streak using serial dilutions to obtain isolated colonies. Incubate under aerobic conditions at 35–37 °C for 24 hours. Refer to the colony color interpretation table.

## LIMITATIONS OF THE METHOD

*Proteus* is a bacterium that may form black colonies like *Salmonella*. This occurs due to its ability to produce hydrogen sulfide (H<sub>2</sub>S) by breaking down ferric ammonium citrate.

## GENERAL QUALITY CONTROL CHARACTERISTICS

### S.S. AGAR

Microorganism	ATCC	Colony Morphology and Color
<i>Escherichia coli</i>	25922	Partial inhibition; red colonies
<i>Enterococcus faecalis</i>	29212	Partial inhibition
<i>Salmonella typhimurium</i>	14028	Colorless colonies with black center
<i>Shigella flexneri</i>	12022	Colorless colonies


S.S. AGAR: *Salmonella typhimurium* ATCC 14028

## GENERAL QUALITY CONTROL CHARACTERISTICS

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## DISPOSAL OF MATERIALS IN WASTE

Materials that show no growth can be considered non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of in accordance with guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must manage and dispose of (or assign the management and disposal of) it according to the applicable regulations.

## SPECIFICATIONS

### S.S. AGAR – RAMBACH AGAR - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Bi-plate 9cm	020640	10 plates	2 – 12 °C	3 months

Produced in Greece by Bioprep in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 521203771414010405SJ EDMA: (14 01 04 05) Bi-Plates – including Non-Chromogenic & Chromogenic Media combination. Bioprep is certified according to the standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 / Ministerial Decision ΔΥ86/1348/2004.

## REFERENCES

Isenberg, H.D., Kominos, S., and Siegel, M. (1969). Isolation of salmonellae and shigellae from an artificial mixture of fecal bacteria. *Appl. Microbiol.*, 18: 4, 656-659.

Leifson, E. (1935). New culture media based on sodium desoxycholate for the isolation of intestinal pathogens and for the enumeration of colon bacilli in

milk and water. J. Pathol. Bacteriol., 40: 581-589.  
Taylor, W.I., Harris, B. (1965) Isolation of shigellae. II. Comparison of plating media and enrichment broths. Am. J. Clin. Pathol. 44: 4, 476-479.  
Wallace, J.S. and K. Jones. 1996. The Use of Selective and Differential Agars in the Isolation of E.coli O157 from dairy herds. J. Appl. Bacteriol. 81: 663-668  
Wallace, J.S. et al. 1997. Isolation of Vero cytotoxin-producing E.coli O157 from wild birds. J. Appl. Bacteriol. 82: 399-404  
Bettelheim K.A.1998. Reliability of CHROMagar O157 for the Detection of Enterohaemorrhagic E.coli (EHEC) O157 but not EHEC Belonging to Other Serogroups. J. Appl. Microbiol. 85: 425 - 428

#### IN VITRO MANUFACTURER'S DATA

**Bioprep**  
*microbiology*

#### G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: [bioprep1@otenet.gr](mailto:bioprep1@otenet.gr) [www.bioprep.gr](http://www.bioprep.gr)