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6th 2025

DESCRIPTION

The bi-plate SABOURAUD ACTIDIONE – DERMATOPHYTE TEST AGAR is used for the isolation, cultivation, and identification of pathogenic fungi. (SABOURAUD ACTIDIONE). Additionally, the Dermatophyte Test Medium base, when supplemented with chloramphenicol and cycloheximide, becomes a selective medium for the isolation of pathogenic fungi. (DERMATOPHYTE TEST AGAR).

PRINCIPLE OF THE METHOD

In SABOURAUD ACTIDIONE, peptones provide a rich source of amino acids and nitrogenous compounds required for the growth of molds and yeasts. Dextrose serves as an energy source. The pH is adjusted to approximately 5.6 to enhance fungal growth, particularly dermatophytes, and to moderately inhibit bacterial growth in clinical samples. Chloramphenicol is an antimicrobial agent that inhibits a wide range of Gram-positive and Gram-negative bacteria. Cycloheximide suppresses the growth of saprophytic fungi.

COMPOSITION SABOURAUD ACTIDIONE	g/litre
Balanced Peptone No. 1	10.0
Dextrose	40.0
Agar No. 2	12.0
Chloramphenicol	0.1
Cycloheximide	0.5

Appearance: Beige, clear

Final pH 5.6 ± 0.2 at 25 °C.

In DERMATOPHYTE TEST AGAR, soy peptone provides essential nitrogenous and carbonaceous nutrients for microbial growth. Dextrose is the energy source for fungal metabolism. The pH indicator phenol red detects acid production. Cycloheximide inhibits most saprophytic fungi. Chloramphenicol is a broad-spectrum antibiotic that inhibits a large group of Gram-positive and Gram-negative bacteria. Gentamicin primarily inhibits Gram-negative bacterial growth.

COMPOSITION DERMATOPHYTE TEST AGAR	g/litre
Soy Peptone	10.0
Dextrose	10.0
Bacteriological Agar	12.0
Phenol Red	0.2
Cycloheximide	0.5
Chloramphenicol	0.05
Gentamicin	0.1

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To SABOURAUD ACTIDIONE – DERMATOPHYTE TEST AGAR is in vitro diagnostic laboratory material and should only be handled by trained laboratory personnel. This material contains peptones and extracts of animal origin. Certificates regarding the source and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious and handled using standard safety precautions (avoid ingestion and inhalation). Always handle plates wearing gloves and within a Class II laminar flow cabinet to prevent contamination, especially by saprophytic fungi. Do not use plates that show signs of microbial contamination. The material should not be used after the expiration date. In case of skin contact, wash immediately with plenty of water and soap. Positive samples must be disposed of according to the appropriate hygiene regulations for infectious materials.

STORAGE CONDITIONS

Plates should be stored at **2 – 12 °C** in their original packaging until use. Prolonged storage below **2 °C** may cause excessive moisture inside the medium, increasing the risk of contamination. Freezing, even momentarily, damages the medium. Excessive heating should also be avoided. Plates may be used up to the expiration date indicated on the label. If the sterile package is accidentally opened, the plate can be stored in the refrigerator for **5 – 7 days** after being sealed with Parafilm or placed in a sealed bag. According to our stability studies, the plates can be transported at **6 – 25 °C** for up to **5 days** or at **25 – 40 °C** for **48 hours** without affecting product performance.

INSTRUCTIONS FOR USE

Place the plates in the incubator at **35 – 37 °C** for 30–45 minutes. Inoculate the sample as soon as possible after collection and streak using serial dilutions to obtain isolated colonies. Incubate plates at **37 °C** for 24 hours for the detection of *Candida* species. Then, incubate at **25–30 °C** under increased humidity. All cultures should be examined at least once per week for fungal growth and should be retained for 4–6 weeks before being considered negative.

READING AND INTERPRETATION**SABOURAUD ACTIDIONE**

After adequate incubation, the plates should show isolated colonies in dilution areas and mixed growth in dense inoculation areas. Examine colonies for color and morphology.

DERMATOPHYTE TEST AGAR

Check the medium at 24 hours for a color change from yellow to red and for detection of *Candida* species. Most pathogenic dermatophytes will produce a full color change within 3 to 6 days. Test interpretation after 14 days is unreliable, as false-positive results may occur. Some *Candida* species may also cause the indicator to turn red but can be distinguished by their characteristic white colonies.

METHOD LIMITATIONS

The full classification of dermatophytes requires microscopic examination of wet mounts and additional biochemical and serological procedures. Limited or normal growth of Gram-positive or Gram-negative bacteria may occur if the inoculum contains $\geq 10^6$ CFU/ml. Biochemical testing of isolated colonies is recommended for complete identification.

GENERAL QUALITY CONTROL CHARACTERISTICS

SABOURAUD ACTIDIONE

Microorganism	ATCC	Growth
<i>Candida albicans</i>	10231	Good
<i>Trichophyton mentagrophytes</i>	28185	Good
<i>Aspergillus brasiliensis</i>	16404	Inhibited

Trichophyton mentagrophytes ATCC 28185

GENERAL QUALITY CONTROL CHARACTERISTICS

DERMATOPHYTE TEST AGAR

Microorganism	ATCC	Growth	Colonies
<i>Microsporum audouinii</i>	9079	Good	Pink, red coloration in the medium
<i>Trichophyton mentagrophytes</i>	9533	Good	Pink, red coloration in the medium
<i>Aspergillus niger</i>	16404	Inhibited	
<i>Pseudomonas aeruginosa</i>	10145	Inhibited	

Trichophyton mentagrophytes ATCC 9533

DISPOSAL OF MATERIALS IN WASTE

Materials that show no growth can be considered non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of in accordance with guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must manage and dispose of (or assign the management and disposal of) it according to the applicable regulations.

SPECIFICATIONS

SABOURAUD ACTIDIONE – DERMATOPHYTE TEST AGAR - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Bi-plate 9cm	020256	10 pieces	2 – 12 °C	5 months

Manufactured in Greece by Bioprep in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 5212037714030201WT EDMA Code: (14 03 02 01) PPM for Yeasts and Fungi Bioprep is certified according to the standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 / DY86/1348/2004.

REFERENCES

SABOURAUD ACTIDIONE

Sabouraud, R. (1910). Les Teignes Paris. Pagano, J., Levin, J.D. and Trejo, W. (1957-8). Diagnostic medium for the differentiation of species of *Candida*. Antibiotics Annual, 137-143.

DERMATOPHYTE TEST AGAR

Allen, A.M., Drewry, R.A., A. Weaver, R.E.: Evaluation Of A New Color Indicator Media For Diagnosis Of Dermatophytosis. -Arch. Derm., 102; 68-70 (1970).

Mertz, W.G., Berger, C.L., A. Silva-Hutner, M.: Media With Ph-Indicator For The Isolation Of Dermatophytes. - Arch. Derm., 99; 203-209 (1969).

Taplin, D., Allen, A.M., A. Mertz, P.M.: Experience With A New Indicator Medium (Dtm) For The Isolation Of Dermatophyte Fungi, In "Proceedings Of The International Symposium Of Mycoses", Scientific Publication 205. Washington, D.C. Pan American Health Organization, 55-58 (1970).

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr www.bioprepare.gr