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E-mail: bioprep1@otenet.gr www.bioprepare.gr

TECHNICAL DATA SHEET

PRODUCT: **SABOURAUD DEXTROSE AGAR** 
W/ CHLORAMPHENICOL

REFERENCE: **010098 – 050098 – 060098 – 070098**

Date 1st Edition:
7th 2009
Date 4th Revision:
6th 2024

DESCRIPTION

SABOURAUD DEXTROSE AGAR WITH CHLORAMPHENICOL is a medium used for isolation, cultivation, and identification of pathogenic and non-pathogenic fungi.

PRINCIPLE OF THE METHOD

Peptones provide a nutrient source of amino acids and nitrogen compounds necessary for the growth of fungi and yeasts. Dextrose serves as an energy source for fungal development. The pH is adjusted to approximately 5.6 to enhance fungal growth, particularly dermatophytes, while slightly inhibiting bacterial growth in clinical samples. Chloramphenicol is an antimicrobial agent that inhibits the growth of a broad spectrum of Gram-positive and Gram-negative bacteria.

COMPOSITION	g/l
Balanced Peptone No. 1	10.0
Dextrose	40.0
Agar No. 2	12.0
Chloramphenicol	0.1

Appearance: Beige, clear

Final pH 5.6 ± 0.2 at 25 °C.

PRECAUTIONS

SABOURAUD DEXTROSE AGAR W/ CHLORAMPHENICOL is an in vitro diagnostic laboratory material and should only be handled by trained laboratory personnel. This material contains peptones and extracts of animal origin. Certificates regarding the source and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious and handled using standard safety precautions (avoid ingestion and inhalation). Always handle plates wearing gloves and within a Class II laminar flow cabinet to prevent contamination, especially by saprophytic fungi. Do not use plates that show signs of microbial contamination. The material should not be used after the expiration date. In case of skin contact, wash immediately with plenty of water and soap. Positive samples must be disposed of according to the appropriate hygiene regulations for infectious materials.

STORAGE AND TRANSPORT CONDITIONS

Plates should be stored at **2 – 12 °C** in their original packaging until the time of use. Prolonged storage at temperatures below **2 °C** may create excess moisture in the medium, increasing the risk of contamination. Freezing, even briefly, will damage the product. Excessive heat should also be avoided. Plates can be used up to the expiration date indicated on the label. Stability studies have shown that during transport, the plates can remain at **18 – 25 °C** for up to **10 days** or at **25 – 40 °C** for up to **48 hours** without affecting product performance.

INSTRUCTIONS FOR USE

Incubate the plates at 37 °C for 24 hours to detect *Candida* fungi. Then continue incubation at 25–30 °C with increased humidity. All cultures should be examined at least once a week for fungal growth and kept for 4–6 weeks before being considered negative.

RESULT INTERPRETATION

After sufficient incubation, plates should show isolated colonies in streaked areas and mixed growth in areas with heavy inoculation. Examine colonies for color and morphology.

METHOD LIMITATIONS

It is recommended to perform biochemical tests on individual colonies for complete identification.

There may be weak or normal growth of Gram-positive and Gram-negative bacteria when the inoculated sample contains $\geq 10^6$ CFU/ml.

GENERAL QUALITY CONTROL CHARACTERISTICS

Microorganism	ATCC	Growth
<i>Candida albicans</i>	10231	Good
<i>Trichophyton mentagrophytes</i>	28185	Good
<i>Sacharomyces cerevisiae</i>	CECT 1045	Good
<i>Escherichia coli</i>	25922	Moderately Inhibited

Candida albicans ATCC 10231

DISPOSAL OF MATERIALS IN WASTE

Materials that show no growth can be considered non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of in accordance with guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must manage and dispose of (or assign the management and disposal of) it according to the applicable regulations.

SPECIFICATIONS

SABOURAUD DEXTROSE AGAR W/ CHLORAMPHENICOL - C€

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
PETRI DISH 9CM	010098	10 pieces	2 – 12 °C	5 months
PETRI DISH 6CM	050098	10 pieces	2 – 12 °C	5 months
TUBE 10ML	070098	20 pieces	2 – 12 °C	8 months
VIAL 100ML	060098	10 pieces	2 – 25 °C	12 months

Produced in Greece by Bioprepare in accordance with Regulation (EU) 2017/746. BASIC UDI-DI: 5212037714030201WT EDMA: (14 03 02 01) PPM for Yeasts and Fungi Bioprepare is certified according to the standards: EN ISO 9001:2015 / ELOT EN ISO 13485:2016 / DY86/1348/2004

REFERENCES

Sabouraud, R. (1910). Les Teignes Paris. Pagano, J., Levin, J.D. and Trejo, W. (1957-8). Diagnostic medium for the differentiation of species of *Candida*. Antibiotics Annual, 137-143.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepare.gr