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Date 1st Edition:
7th 2009
 Date 5th Revision:
7th 2025

DESCRIPTION

SABOURAUD DEXTROSE AGAR is a medium used for the isolation, cultivation, and identification of pathogenic and non-pathogenic fungi.

PRINCIPLE OF THE METHOD

Peptones serve as a source of nitrogenous growth factors. Dextrose provides an energy source for fungal growth. Agar solidifies the medium.

FORMULA	g/litre
Balanced Peptone No. 1	10.0
Dextrose	40.0
Agar	15.0

Appearance: Beige clear.

Final pH $5,6 \pm 0,2$ at 25 °C.

PRECAUTIONS

SABOURAUD DEXTROSE AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The plates should be stored at **2–12 °C** in their packaging until ready for use. Prolonged storage at temperatures below **2 °C** creates significant moisture within the medium, increasing the risk of contamination. Freezing, even momentarily, destroys the medium. Excessive heating should also be avoided. The plates can be used until the expiration date indicated on the label. If the airtight packaging of a plate is accidentally opened, it can be stored in the refrigerator for **5–7 days** if resealed with parafilm or a bag. Stability studies have shown that the plates can remain at **18–25 °C** for up to **10 days** or at **25–40 °C** for up to **48 hours** without compromising product performance. Vials should be stored at **8–25 °C**. Freezing, even briefly, destroys the material.

USAGE

Spread the material on the plates as soon as possible after receiving it in the laboratory, using a sterile swab to obtain isolated colonies. For isolating fungi from contaminated samples, inoculate both a selective and a non-selective medium. Incubate the plates at 37 °C for 24 hours for the detection of *Candida* fungi. Then, incubate at 25–30 °C with increased humidity. All cultures should be examined at least once a week for fungal growth and should be kept for 4–6 weeks before being considered negative.

READING AND INTERPRETATION

After sufficient incubation, the plates should display isolated colonies in the dilution areas and mixed growth in the areas of dense inoculation. Examine the colonies for their colour and morphology.

LIMITATIONS OF THE METHOD

It is recommended to perform biochemical tests on the isolated colonies for complete identification. Weak or normal growth of Gram-positive and Gram-negative bacteria may occur when their concentration in the inoculum is $\geq 10^6$ CFU/ml.

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	ATCC	Growth
<i>Candida albicans</i>	10231	Good
<i>Trichophyton mentagrophytes</i>	28185	Good
<i>Sacharomyces cerevisiae</i>	CECT 1045	Good
<i>Aspergillus brasiliensis</i>	16404	Good
<i>Escherichia coli</i>	25922	Moderate inhibition

Candida albicans ATCC 10231

Aspergillus brasiliensis ATCC 16404

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

SABOURAUD DEXTROSE AGAR (USP/EP) - 

PRODUCT	CODE	PACKING	STORE	SELF LIFE
Plate 9cm	010292	10 pieces	2 – 12 °C	5 months
Plate 6cm	050292	10 pieces	2 – 12 °C	5 months
Bottle 100ml	060292	10 pieces	2 – 25 °C	12 months
Tube 100 ml	070292	40 pieces	2 – 25 °C	12 months

Produced in Greece by Bioprep in accordance with the requirements of the European Directive 98/79/EK. FEK B2198/2-10-2009. CODE EDMA 14 03 02 01. Bioprep company has been certified according to standards EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004.

BIBLIOGRAPHY

Sabouraud, R. (1910). Les Teignes Paris. Pagano, J., Levin, J.D. and Trejo, W. (1957-8). Diagnostic medium for the differentiation of species of *Candida*. Antibiotics Annual, 137-143.
European Pharmacopoeia 9th Edition (2017)

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 19001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprep.gr