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Sabouraud, R. 1892. Ann. Dermatol. Syphilol. 3:1061

Jarett, L., and A.C. Sonnenwirth (ed) 1980. Gradwohl's clinical laboratory methods and diagnosis, 8th ed. CV Mosby

Davidson, A.M., E.S Dowding, and A.H.R Buller. 1932 Hyphal fusions in dermatophytes. Can J. Res. 6:1.

Association of Official Analytical Chemists. 1995. Bacteriological analytical manual, 8th ed. AOAC International, Gaithersburg, MD
European Pharmacopoeia. 7.0

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E-mail: bioprep1@otenet.gr www.bioprepare.gr

TECHNICAL DATA SHEET

PRODUCT: **SABOURAUD DEXTROSE BROTH** 
REFERENCE: **070605 – 060605**

Date 1st Edition:
7th 2009
Date Revision:
9th 2025

DESCRIPTION

SABOURAUD DEXTROSE BROTH is used for the cultivation of yeasts and fungi.

PRINCIPLE OF THE METHOD

The high concentration of dextrose and the acidic pH make this medium selective for fungi from clinical specimens and other sources. It is employed for the cultivation of yeasts and pathogenic fungi, particularly those associated with dermal infections. It is also used in sterility testing. The formulation complies with the European Pharmacopoeia. Dextrose provides carbon and energy, while the peptone mixture supplies nitrogen, vitamins, minerals, and amino acids essential for growth.

COMPOSITION	g/litre
Dextrose	20.00
Mixture of Peptic of Animal Tissue and Pancreatic Digest of Casein	10.00

Appearance: Beige to light orange

Final pH 5,6 ± 0.2 at 25 °C.

PRECAUTIONS

To SABOURAUD DEXTROSE BROTH is in vitro diagnostic laboratory material and should only be handled by trained laboratory personnel. This material contains peptones and extracts of animal origin. Certificates regarding the source and health status of the animals do not fully guarantee the absence of transmissible pathogens. Therefore, these materials should be treated as potentially infectious and handled using standard safety precautions (avoid ingestion and inhalation). Always handle plates wearing gloves and within a Class II laminar flow cabinet to prevent contamination, especially by saprophytic fungi. Do not use plates that show signs of microbial contamination. The material should not be used after the expiration date. In case of skin contact, wash immediately with plenty of water and soap. Positive samples must be disposed of according to the appropriate hygiene regulations for infectious materials.

STORAGE CONDITIONS

The tubes should be stored at **2 – 12 °C** in their original packaging until the time of use. Even momentary freezing will damage the medium. Excessive heating should also be avoided. The tubes may be used until the expiry date indicated on the label. Based on our stability studies, the tubes can be transported at **6 – 25 °C** for up to **7 days** or at **25 – 40 °C** for up to **48 hours** without compromising product performance.

INSTRUCTIONS FOR USE

Inoculate the broth by collecting an isolated colony from a pure culture, e.g., *Candida albicans*, using a loop. Incubate the tubes with loosely fitted caps at 30 – 35 °C for 3 to 5 days. After incubation, subculture on a SABOURAUD DEXTROSE AGAR plate (Code: 010292) at 30 – 35 °C for 24 to 48 hours.

GENERAL QUALITY CONTROL CHARACTERISTICS

MICROORGANISM	GROWTH
<i>Aspergillus brasiliensis</i> ATCC 16404	Good
<i>Candida albicans</i> ATCC 26790	Good
* <i>Candida albicans</i> ATCC 10231	Good (*)
<i>Escherichia coli</i> ATCC 25922	Partial inhibition
<i>Lactobacillus casei</i> ATCC 9595	Good
<i>Saccharomyces cerevisiae</i> ATCC 9763	Good

(*) According to the European harmonized method:

Incubate at 20–25 °C for 2–3 days and for specific microorganisms at 30–35 °C for 3–5 days.

DISPOSAL OF MATERIALS IN WASTE

Materials that show no growth can be considered non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of in accordance with guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must manage and dispose of (or assign the management and disposal of) it according to the applicable regulations.

SPECIFICATIONS

SABOURAUD DEXTROSE BROTH - 

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Tube 7 ml	070605	20 pcs	2 – 12 °C	12 months
Bottle 100 ml	060605	10 pcs	2 – 25 °C	12 months

Produced in Greece by Bioprepare in accordance with Regulation (EU) 2017/746. Basic UDI-DI: 5212037714030301WY EDMA Code: (14 03 03 01) Tubes for Yeasts and Fungi Bioprepare is certified according to the following standards: EN ISO 9001:2015 / EAOT EN ISO 13485:2016 / ΔΥ86/1348/2004.

REFERENCES

Sabouraud, R. 1892. Ann. Dermatol. Syphilol. 3:1061
Jarett, L., and A.C. Sonnenwirth (ed) 1980. Gradwohl's clinical laboratory methods and diagnosis, 8th ed. CV Mosby
Davidson, A.M., E.S Dowding, and A.H.R Buller. 1932 Hyphal fusions in dermatophytes. Can J. Res. 6:1.
Association of Official Analytical Chemists. 1995. Bacteriological analytical manual, 8th ed. AOAC International, Gaithersburg, MD

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr www.bioprepare.gr