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**TECHNICAL DATA SHEET**

PRODUCT: **T.C.B.S. AGAR**
 REFERENCE: **010109**

Date 1st Edition:
 7th 2009
 Date 4th Revision:
 6th 2024

DESCRIPTION

TCBS (Thiosulfate Citrate Bile Salt Agar) is a highly selective medium used for the selective isolation of *Vibrio cholerae* and *Vibrio parahaemolyticus* from a variety of clinical and non-clinical samples.

PRINCIPLE OF THE METHOD

Enterobacteriaceae and *Proteus* species that may be present in stool samples are inhibited. TCBS can also be used for the isolation of *Vibrio fluvialis* and *Vibrio vulnificus*. Gram-positive bacteria are inhibited by the presence of bile salts in the medium. Sucrose is included as a carbohydrate source for the metabolism of *Vibrio* species. The alkaline pH enhances the growth of *Vibrio cholerae*. Thymol blue and bromothymol blue function as pH indicators.

COMPOSITION	g/litre
Yeast extract	5
Peptone mix	10
Sodium Citrate	10
Sodium Thiosulfate	10
Bile salts	9
Sucrose	17
Sodium chloride	10
Ferric Ammonium citrate	1
Bromthymol blue	0.04
Thymol blue	0.04
Agar No. 1	14

Appearance: Green clear.

Final pH 8.6 ± 0.2 in 25 °C

PRECAUTIONS

T.C.B.S. AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE AND TRANSPORTS CONDITIONS

The plates must be stored at **2 – 12 °C** in their packaging until the time of use. Prolonged storage at temperatures below **2 °C** creates excessive moisture inside the medium, increasing the risk of contamination. Freezing, even momentarily, destroys the medium. Excessive heating should also be avoided. The plates can be used until the expiration date indicated on the label. Stability studies have shown that during transport, the plates can remain at **6 - 25 °C** for **4 days** or at **25 - 40 °C** for **24 hours** without affecting product performance.

INSTRUCTIONS FOR USE

The infectious material must reach the laboratory quickly without delay and must be protected from excessive heat or cold. If processing is delayed, the sample should be inoculated into a suitable transport medium and kept at 2 - 12 °C until inoculation. Allow the plates to reach room temperature before use. The agar surface must be dry before inoculation. Inoculate the sample at one edge of the plate and then streak with a loop in parallel lines to obtain isolated colonies. Incubate the plates aerobically at 35 - 37 °C for 24 to 48 hours.

INTERPRETATION OF RESULTS

Vibrio parahaemolyticus forms blue colonies with a green center.

LIMITATIONS OF THE METHOD

An oxidase test (Oxidase Test Strips - Code MID-61G) is recommended on colonies from a pure culture for complete identification. Final identification may require additional tests, such as biochemical assays or immunological tests (latex agglutination).

GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	Growth - Colony Color
<i>Vibrio cholerae</i> ATCC® 9459	Good - Yellow
<i>Vibrio parahaemolyticus</i> ATCC® 17802	Good – Blue with blue center
<i>Escherichia coli</i> ATCC® 25922	Inhibited
<i>Enterococcus faecalis</i> ATCC® 13124	Inhibited

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

T.C.B.S. AGAR - CE

PRODUCT	CODE	PACKING	STORE	SELF LIFE
Plate 9cm	010109	10 pieces	2 – 12 °C	2 months

Produced in Greece by Bioprepere in accordance with the requirements of the European directive 98/79/EK. FEK B2198/2-10-2009. Code EDMA 14 01 04 01. Bioprepere company has been certified according to standards EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004

BIBLIOGRAPHY

Kobayashi, T., Enomoto, S., Sakazaki, R. and Kuwahara, S. (1963). Jap. Bacteriol 18: 10-11, 387-391.
Furniss, A.L., Lee J.V. and Donovan, T.J. (1978). The Vibrios, PHLS Monograph No. 11.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepere.gr