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##### BLOOD AGAR

Ellner, P.D., Stoessel, C.J., Drakeford, E. and Vasi, F. (1966). A new culture medium for medical bacteriology. Amer J. Clin Pathol. 45. 502-504.  
American Public Health Association (1950). Diagnostic Procedures and Reagents. 3rd edn. A.P.H.A., New York.

##### MAC CONKEY AGAR

American Public Health Association (1950). Diagnostic Procedures and Reagents. 3rd edn. A.P.H.A., New York.

American Public Health Association (1946). Standard Methods for the examination of Water and Sewage. 9th edn. A.P.H.A., New York.

MANNITOL SALT AGAR (CHAPMAN)

American Public Health Association (1966). Recommended Methods for Microbiological Examination of Foods, 2nd Edn. (ed. J.M. Sharf) A.P.H.A. Washington.

Davis, J.G., (1959). Milk Testing 2nd edn, Dairy Industries, London.

European Pharmacopeia 6th Edition Supplement 6.3. Council of Europe European (COE) - European Directorate for the Quality of Medicines (EDQM). June 2008

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PRODUCT: **BLOOD AGAR – MAC CONKEY AGAR –  
MANNITOL SALT AGAR (CHAPMAN)**

REFERENCE: **030155**

Date 1st Edition:  
7th 2009  
Date 4th Revision:  
6th 2024

**DESCRIPTION**

A tri-sectioned plate containing BLOOD AGAR – MAC CONKEY AGAR – MANNITOL SALT AGAR for the cultivation of all aerobic bacteria and the identification of hemolytic Gram-positive cocci (BLOOD AGAR). It is also used for the cultivation, isolation, and identification of Enterobacteriaceae (MAC CONKEY AGAR), where Enterococcus does not grow. Additionally, it supports the cultivation, isolation, and identification of pathogenic Staphylococcus (MANNITOL SALT AGAR).

**PRINCIPLE OF THE METHOD**

The Blood Agar Base contains highly nutritious components that provide essential vitamins, carbohydrates, and other organic elements. Sodium chloride supplies necessary minerals and helps maintain osmotic balance and electrolyte equilibrium. The addition of 6% horse blood enhances the nutritive properties of the medium. Furthermore, hemolysis and the type of hemolysis observed in horse red blood cells assist in the identification of specific Gram-positive cocci.

FORMULA BLOOD AGAR	g/litre
Columbia Peptone Mixture	25.1
Corn Starch	1.0
Sodium chloride	5.0
Agar No. 2	12.0
Horse Blood	60ml

Appearance: Red to burgundy, non-transparent due to the addition of blood.

Final pH  $7.3 \pm 0.2$  at 25 °C.

MAC CONKEY AGAR is a nutrient medium used for the isolation of Enterobacteriaceae. Gram-negative enterobacteria that ferment lactose produces red or pink colonies. Bile salts No. 3 and crystal violet inhibit the growth of Gram-positive cocci. The neutral red indicator changes color in response to lactose fermentation.

FORMULA MAC CONKEY AGAR	g/litre
Peptone	20.0
Lactose	10.0
Bile Salts No. 3	1.5
Sodium chloride	5.0
Neutral red	0.03
Crystal violet	0.001
Agar No. 2	15.0

Appearance: Pink to purple, transparent.

Final pH  $7.1 \pm 0.2$  at 25°C.

In Mannitol Salt Agar (CHAPMAN), peptones and beef extract provide essential growth factors such as nitrates, carbon, and trace elements. The high sodium chloride concentration partially or completely inhibits the growth of all bacteria except Staphylococcus. Mannitol fermentation, indicated by a color change in the phenol red indicator, aids in the differentiation of Staphylococcus species.

FORMULA MANNITOL SALT AGAR (CHAPMAN)	g/litre
Beef Extract	1.0
Balanced Peptone No1	10.0
Sodium chloride	75.0
D-Mannitol	10.0
Agar No 2	12.0
Phenol Red	0.025

Appearance: Red, transparent.

Final pH  $7.4 \pm 0.2$  at 25 °C.

**PRECAUTIONS**

To BLOOD AGAR – MAC CONKEY AGAR – MANNITOL SALT AGAR (CHAPMAN) is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the

skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

### STORAGE CONDITIONS

The plates should be stored at 2 – 12 °C within their packaging until use. Prolonged storage at temperatures below 2 °C may cause excessive moisture buildup inside the material, risking contamination. Freezing, even momentarily, damages the material. Additionally, excessive heating should be avoided. The plates can be used until the expiration date indicated on the label. For transport, our stability studies have shown that the plates can remain at 6 – 25 °C for 4 days or at 25 – 40 °C for 48 hours without affecting product performance.

### METHOD OF USE

Inoculate the plates using the technique for obtaining isolated colonies. Incubate at 35 – 37 °C for 24 hours under aerobic conditions.

### INTERPRETATION OF RESULTS

#### BLOOD AGAR:

*Streptococcus pyogenes* forms small, grayish-white colonies with β-hemolysis.

*Streptococcus pneumoniae* forms small, flat colonies with α-hemolysis.

*Staphylococcus aureus* forms relatively large colonies (1-3mm after 24 hours of incubation, 3-8mm if incubation is extended to 5 days). The color ranges from cream-yellow to orange depending on the incubation time. It causes either α or β hemolysis.

#### MAC CONKEY AGAR:

Gram-positive cocci do not grow.

Gram-negative bacteria that ferment lactose produce pink to red colonies due to acid production. Additionally, *Escherichia coli* forms pink or red colonies surrounded by a cloudy zone due to the precipitation of bile salts because of pH decrease.

Gram-negative bacteria that do not ferment lactose form colorless colonies.

*Proteus* does not swarm.

#### MANNITOL SALT AGAR (CHAPMAN):

*Staphylococcus aureus* colonies appear white yellow, with a yellow zone in the medium due to mannitol fermentation.

Other *Staphylococcus* species form red colonies without any change in the color of the medium since they do not ferment mannitol.


### LIMITATIONS OF THE METHOD

Final identification should be confirmed with biochemical and serological tests (e.g., Latex Agglutination Test), which can be performed directly from the suspect colonies.

### QUALITY CONTROL

#### BLOOD AGAR

Microorganism	ATCC	Growth	Haemolysis
<i>Escherichia coli</i>	25922	Good	-
<i>Pseudomonas aeruginosa</i>	27853	Good	-
<i>Staphylococcus aureus</i>	25923	Good	B
<i>Streptococcus pneumoniae</i>	6303	Good	A
<i>Streptococcus pyogenes</i>	19615	Excellent	B


*Streptococcus pyogenes* ATCC 19615

#### MAC CONKEY AGAR

Microorganism	ATCC	Growth	Colonies
<i>S. typhimurium</i>	14028	Good	Colorless - translucent
<i>P. mirabilis</i>	12453	Good	Colorless, does not swarm
<i>Escherichia coli</i>	25922	Good	Pink or red colonies surrounded by a cloudy zone
<i>E. faecalis</i>	29212	Inhibited	


*Escherichia coli*

*P. mirabilis*

**MANNITOL SALT AGAR (CHAPMAN)**

Microorganism	Growth	Colonies
<i>Staphylococcus aureus</i> ATCC 25923	Good	Yellow colonies, with a yellow zone in the medium.
<i>Staphylococcus epidermidis</i> ATCC 12228	Good	Red colonies without a zone or color change.
<i>Enterobacter aerogenes</i> ATCC 13048	Partial to complete inhibition	
<i>Escherichia coli</i> ATCC 25922	Partial to complete inhibition	
<i>Proteus mirabilis</i> ATCC 7002	Partial inhibition	

*Staphylococcus aureus* ATCC 25923

**WASTE DISPOSAL OF WASTE**

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

**SPECIFICATIONS**

BLOOD AGAR – MAC CONKEY AGAR – MANNITOL SALT AGAR (CHAPMAN) - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Tri-sectioned plate 9cm	030155	10 pieces	2 – 12 °C	2 months

Produced in Greece by the company Bioprepare in accordance with the requirements of the European Directive 2017/746.

BASIC UDI-DI: 5212037714010401WF. EDMA (14 01 04 01) non-chromogenic media (Plates).

The Bioprepare company has been certified according to the standards: EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004

**LITERATURE REFERENCES**

**BLOOD AGAR**

Ellner, P.D., Stoessel, C.J., Drakeford, E. and Vasi, F. (1966). A new culture medium for medical bacteriology. Amer J. Clin Pathol. 45. 502-504.  
American Public Health Association (1950). Diagnostic Procedures and Reagents. 3rd edn. A.P.H.A., New York.

**MAC CONKEY AGAR**

American Public Health Association (1950). Diagnostic Procedures and Reagents. 3rd edn. A.P.H.A., New York.  
American Public Health Association (1946). Standard Methods for the examination of Water and Sewage. 9th edn. A.P.H.A., New York.

**MANNITOL SALT AGAR (CHAPMAN)**

American Public Health Association (1966). Recommended Methods for Microbiological Examination of Foods, 2nd Edn. (ed. J.M. Sharf) A.P.H.A. Washington.  
Davis, J.G., (1959). Milk Testing 2nd edn, Dairy Industries, London.  
European Pharmacopoeia 6th Edition Supplement 6.3. Council of Europe European (COE) - European Directorate for the Quality of Medicines (EDQM). June 2008

IN VITRO MANUFACTURER'S DATA

**G. PAPANIKOLAOU & CO**

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: [bioprep1@otenet.gr](mailto:bioprep1@otenet.gr) [www.bioprepare.gr](http://www.bioprepare.gr)