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6th 2024

**DESCRIPTION**

Bordet Gengou Agar Base with 15% sheep blood is used to isolate *Bordetella pertussis* and other *Bordetella* strains.

**PRINCIPLE OF THE METHOD**

Potato infusion and Proteose peptone provide nitrogen, vitamins, minerals and amino acids essential for growth. Glycerol is a carbon source. Sodium chloride maintains osmotic balance in the material. The addition of 15% sheep's blood makes the material particularly nutritious for *Bordetella* strains. Many factors can inhibit the growth of *Bordetella pertussis*, such as fatty acids, or cotton from the stylet. Starch from the potato extract absorbs fatty acids. The colonies of *Bordetella pertussis* appear small, white, with a characteristic haemolysis zone. This zone is not obvious if the material contains more than 20% blood. Methicillin inhibits the growth of *Staphylococcus* which produces a toxin that inhibits the growth of *Bordetella pertussis*.

COMPOSITION	g/litre
Potato Infusion	4.5 g
Proteose Peptone	10,0 g
Sodium Chloride	5.5 g
Bacteriological Agar	16.0 g
Glycerol	10 ml
Sheep Blood	150 ml

Appearance: Crystal clear due to the addition of blood.

Final pH 6,7 ± 0.2 in 25 °C.

**PRECAUTIONS**

BORDET GENGOU AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

**STORAGE AND TRANSPORTS CONDITIONS**

Petri dishes must be stored in **2 - 12°C** in their packaging until the time of their use. Prolonged storage at a temperature below **2 °C** creates enough moisture in the material at risk of contamination. Freezing even instantly destroys the material. Avoid overheating also. Petri dishes must be used until the expiration date which is shown on the label. When you open the airtight package, if you have any leftover plates, store them in the bag until the expiration date. For transport, our stability studies showed that the plates may remain at **17-25°C** for **4 days** or at **27-40°C** for **48 hours**, without affecting the attribution of the product.

**USAGE**

Roll the swab on the surface of the plate and inoculate for isolation

Incubate the plates at 35 ° C for 5-7 days in a humid environment. The presence of CO2 is not recommended. Isolates suspected isolates for *Bordetella pertussis* should be confirmed using antiserum or staining procedure.

**INTERPRETATION OF RESULTS**

The development of *Bordetella pertussis* will occur in 3-5 days. Other strains can grow in 1-3 days. Pharyngeal specimens may contain *Staphylococcus* which produces a substance that inhibits the growth of *Bordetella pertussis*. For such samples, use a material with methicillin 2.2µg / ml or cephalexin 40µg / ml.

**LIMITATIONS OF THE METHOD**

Isolates suspected isolates for *Bordetella pertussis* should be confirmed using antiserum or staining procedure, or by a staining process.

The final identification must be done by biochemical and serological tests. (e.g., a Latex Test welding test and can be run directly from suspicious colonies.

**GENERAL CHARACTERISTICS OF QUALITY CONTROL**

Microorganism	ATCC	Growth	Haemolysis
<i>Bordetella bronchiseptica</i>	19395	Good	No
<i>Bordetella pertussis</i>	8467	Medium to good	Yes

*Bordetella bronchiseptica*

#### WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

#### SPECIFICATIONS

##### BORDET GENGOU AGAR - CE

PRODUCT	CODE	PACKAGE	STORAGE	SELF LIFE
Petri dish 90mm-20ml	010014	10 items	2 – 12 °C	2 months

Produce in Greece by BIOPREPARE Company in accordance with regulations of the European committee 98/79/EK. Code EDMA 14 01 04 90. Bioprepare company has been certified according to standards EN ISO 9001:2015 / EAOT EN ISO 13485:2016.

#### BIBLIOGRAPHY

1. Hoppe. 1999. In Murray, Baron, Pfaller, Tenover and Tenover (ed.), Manual of clinical microbiology 7th ed. American Society for Microbiology, Washington, D.C.
2. Bordet and Gengou. 1906. Ann. Inst. Pasteur 20:731.
3. Kendrick and Eldering. 1934. Am. J. Public Health 24:309

#### IN VITRO MANUFACTURER'S DATA

#### G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: [bioprep1@otenet.gr](mailto:bioprep1@otenet.gr)

[www.bioprep.gr](http://www.bioprep.gr)