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Herbert G.A., Hollis D.G., Weaver R.E., Karmali M.A., Simor A.E., Roscoe M., Fleming P.C., Smith, S.S. Lane J. Evaluation of a Blood-Free, Charcoal-Based, Selective Medium for the Isolation of *Campylobacter* organisms from Faeces. J. Clin. Micro. (1986) 23 p456-459.

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6th 2024

DESCRIPTION

A selective enrichment broth for the isolation of *Campylobacter spp.* from food, environmental samples, and feces according to ISO 10272 (food) and ISO 17995:2019 (water).

METHOD PRINCIPLE

Peptone and beef extract provide nitrogen, vitamins, and essential minerals necessary for the growth of *Campylobacter*. Sodium chloride supplies essential ions and regulates osmotic balance. Sodium pyruvate, sodium metabisulphite, and ferrous sulphate enhance microaerophilic growth. Polymyxin B inhibits Gram-negative bacteria except for *Proteus*, as some strains exhibit resistance. Trimethoprim mainly inhibits *Proteus* species as well as other Gram-negative bacteria. Rifampicin also acts against Gram-negative bacteria. Amphotericin B functions as an antifungal agent. Lysed horse blood enriches the broth with hemin and other components essential for *Campylobacter* growth.

FORMULA	g/litre
Beef Extract	10.0
Peptone	10.0
Sodium chloride	5.0
Sodium pyruvate	0.25 g
Sodium metabisulphite	0.25 g
Ferrous sulphate	0.25 g
Rifampicin	10.0mg
Polymyxin B	5000 IU
Trimethoprim	10.0mg
Amphotericin B	10.0mg
HORCE BLOOD LYSED	50ml

Appearance: A red (wine-colored) broth with a fine black suspension.

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CAMPYLOBACTER BROTH PRESTON is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

STORAGE CONDITIONS

The tubes should be stored at 2–12°C in their packaging until use. Freezing, even momentarily, damages the material. Excessive heating should also be avoided. The tubes can be used until the expiration date indicated on the label. For transportation, stability studies have shown that the tubes can remain at 6–25°C for up to 4 days or at 25–40°C for up to 48 hours without affecting product performance.

INSTRUCTIONS FOR USE

Add 0.5 g or 0.5 ml of fresh feces or another sample into the Preston broth and mix well. Incubate the broth under microaerophilic conditions at 35–37°C or 42°C (for better isolation) for 18–24 hours. Using a 10 µl loop, inoculate the edge of the CAMPYLOBACTER AGAR (SELECTIVE MEDIUM) (010020) or CAMPYLOBACTER BLOOD-FREE AGAR (mCCD) (010021) plates. Perform serial dilutions using a sterile loop. Incubate at 41.5°C ± 1°C for 44 ± 4 hours under microaerophilic conditions optimized for *Campylobacter* growth (5–6% oxygen, 10% CO₂, and 84–85% nitrogen).

READING AND INTERPRETATION

Campylobacter jejuni on the medium will appear as small, mucous-like colonies, usually grayish, flat with irregular edges, and non-hemolytic. Another colony morphology includes round colonies with a diameter of 1–2 mm, which are convex and shiny. A small percentage of colonies may appear pink. The colonies tend to spread, especially when isolated from a fresh clinical sample. If plates need to be examined after 24 hours of incubation, they should be checked quickly and immediately returned to a reduced-oxygen atmosphere after examination.

METHOD LIMITATIONS

It is recommended to perform biochemical tests on individual colonies for complete identification.

QUALITY CONTROL

Microorganism	ATCC	Growth
<i>Campylobacter jejuni</i>	33291	Good
<i>Escherichia coli</i>	25922	Inhibited

WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

SPECIFICATIONS

CAMPYLOBACTER BROTH PRESTON - CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Tube 3ml	080022	10 pieces	2 – 12 °C	3 months
Vial 90ml	040022	10 pieces	2 – 12 °C	3 months

Produced in Greece by the company Bioprepate in accordance with the requirements of the European Directive 2017/746.

BASIC UDI-DI: 5212037714010401WF. EDMA (14 01 04 01) non-chromogenic media (Plates).

The Bioprepate company has been certified according to the standards: EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004

LITERATURE REFERENCES

Bolton F.J. Hutchinson D.N., Parker G. Reassessment of Selective Agars and Filtration Techniques for Isolation of *Campylobacter* Species from Feces. *Eur.J. Clin. Microbiol. Infects. Dis.* (1988) 7 p 155-160.

Bolton F. J. (1988) Personal Communication.

Bolton F.J. Hutchinson D.N., Parker G. Isolation of *Campylobacter*: What are we missing? *J.Clin.Path.* (1987) 40 p 702-703.

Goosens H., De Boeck M., Coignau H., Vlaes L., Van Den Borre C., Butzler J.P. Modified Selective Medium for Isolation of *Campylobacter* spp from Feces: Comparison with Preston Medium, a Blood Free Medium, and a Filtration System. *J.Clin. Micro.* (1986) 24 p 840-843.

Gun-Munro J., Rennie R.P., Thornley J.H. Richardson H.L., Hodge D., Lynch J. Laboratory and Clinical Evaluation of Isolation Media for *Campylobacter jejuni* *J. Clin Micro.* (1987). 25 p2274-2277.

Herbert G.A., Hollis D.G., Weaver R.E., Karmali M.A., Simor A.E., Roscoe M., Fleming P.C., Smith, S.S. Lane J. Evaluation of a Blood-Free, Charcoal-Based, Selective Medium for the Isolation of *Campylobacter* organisms from Faeces. *J. Clin. Micro.* (1986) 23 p456-459.

IN VITRO MANUFACTURER'S DATA

G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: bioprep1@otenet.gr

www.bioprepate.gr