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6th 2024

**DESCRIPTION**

Split plate CHOCOLATE AGAR – SABOURAUD DEXT. AGAR for the cultivation of microorganisms with specific nutritional requirements, such as Haemophilus, Neisseria. (CHOCOLATE AGAR) – and for the cultivation of fungi. It contains chloramphenicol to inhibit Gram-positive and Gram-negative bacteria. (SABOURAUD DEXT. AGAR).

**METHOD PRINCIPLE**

Chocolate Agar is an enriched medium for qualitative procedures, used for the isolation and cultivation of pathogenic microorganisms, particularly Neisseria species and Haemophilus species, from a variety of clinical samples.

FORMULA	g/litre
CHOCOLATE AGAR	
Columbia Peptone Mixture	23.0
Corn Starch	1.0
Sodium chloride	5.0
Agar No. 2	12.0
Horse Blood	60ml

Appearance: Brown-chocolate, non-transparent, due to the addition of blood.

Final pH  $7.3 \pm 0.2$  at 25 °C

In SABOURAUD DEXT. AGAR, peptones provide a nutritional source of amino acids and nitrogenous compounds for the growth of fungi and yeasts. Dextrose serves as an energy source for fungal growth. The pH is adjusted to around 5.6 to enhance the growth of fungi, particularly dermatophytes, while slightly inhibiting bacterial growth in clinical samples. Chloramphenicol is an antimicrobial agent that inhibits the growth of a broad spectrum of Gram-positive and Gram-negative bacteria.

FORMULA SABOURAUD DEXT. AGAR	g/litre
Balanced Peptone No. 1	10.0
Dextrose	40.0
Agar No. 2	12.0
Chloramphenicol	0.1

Appearance: Beige, transparent,

Final pH  $5.6 \pm 0.2$  at 25 °C.

**PRECAUTIONS**

CHOCOLATE AGAR – SABOURAUD DEXT. AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

**STORAGE AND TRANSPORT CONDITIONS**

The plates should be stored at **2 – 12°C** inside their packaging until use. Prolonged storage at temperatures below **2°C** may cause excessive moisture buildup inside the material, risking contamination. Freezing, even momentarily, will destroy the material. Excessive heating should also be avoided. The plates can be used until the expiration date printed on the label. For transport, stability studies have shown that the plates can be kept at **6 – 25°C** for **5 days** or at **25 – 40°C** for **48 hours** without affecting the product's performance.

**METHOD OF USE**

Inoculate and spread the sample as soon as possible after receiving it in the laboratory. Alternatively, if the sample is to be cultured after collection with a swab, the following procedure can be followed:

1. Pass the swab over the agar, drawing a large "Z" on its surface. This ensures the swab is exposed to the medium for enough time, allowing better transfer of microorganisms.
2. Spread the material using a sterile loop, crossing in a perpendicular direction over the initial "Z." This process is preferably done during the sample collection. If not, it should be performed in the laboratory.
3. Place cultures in an aerobic environment enriched with 10% CO<sub>2</sub> as soon as possible.
4. Incubate at 35-37°C and examine the cultures after overnight incubation as an initial phase and again after approximately 48 hours in the final phase.

## INTERPRETATION OF RESULTS

### CHOCOLATE AGAR:

*Haemophilus influenzae* forms small, moist colonies with a faint odor.

*Streptococcus pneumoniae* forms small, flat colonies that exhibit a green color.

*Neisseria meningitidis* forms large colonies, gray to blue in color.

*Neisseria gonorrhoeae* forms small, gray, white, colorless, mucoid colonies.

### SABOURAUD DEXT. AGAR:

After sufficient incubation, the plates should show isolated colonies of *Candida albicans* in the dilution areas. The colonies are white, circular, creamy, and their size varies depending on the incubation time, ranging from 2 to 10 mm.

## LIMITATIONS OF THE METHOD

Final identification should be performed with biochemical and serological tests (e.g., latex agglutination test), which can be conducted directly from the suspected colonies.

## QUALITY CONTROL

### CHOCOLATE AGAR

Microorganism	ATCC	Colonies
<i>Haemophilus influenzae</i>	10211	Small, moist, with a faint odor.
<i>Streptococcus pneumoniae</i>	6305	Small, flat, exhibiting a green color.
<i>Neisseria meningitidis</i>	13090	Large colonies, gray to blue in color.
<i>Neisseria gonorrhoeae</i>	43069	Small, gray, white, colorless, mucoid.

*Haemophilus influenzae* ATCC 10211

### SABOURAUD DEXT. AGAR

Microorganism	ATCC	Growth
<i>Candida albicans</i>	10231	Good
<i>Trichophyton mentagrophytes</i>	28185	Good
<i>Sacharomyces cerevisiae</i>	CECT 1045	Good
<i>Escherichia coli</i>	25922	Moderately inhibited.

*Candida albicans* ATCC 10231

## WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

## SPECIFICATIONS

### CHOCOLATE AGAR – SABOURAUD DEXT. AGAR – CE

ITEM	CODE	PACKAGE	STORAGE	SHELF LIFE
Bi - plate 9cm	020139	10 pieces	2 – 12 °C	3 months

Produced in Greece by the company Bioprepate in accordance with the requirements of the European Directive 2017/746.

BASIC UDI-DI: 5212037714010401WF. EDMA (14 01 04 01) non-chromogenic media (Plates).

The Bioprepate company has been certified according to the standards: EN ISO 9001:2015 / EAOT EN ISO 13485:2016 DY8d/1348/2004

## LITERATURE REFERENCES

### CHOCOLATE AGAR

Sng, E.H., Rajan, V.S., and Lim, A.L. (1977) Simplified media for isolating *Neisseria gonorrhoeae* J. Clin. Microbiol., 5 (4), 387.

Sottnek, F. O., Biddle, J. W., Kraus, S. J., Weaver, R. E., and Stewart, J. A. (1980) Isolation and Identification of *Haemophilus ducreyi* in a clinical study. J. Clin. Microbiol., 12 (2), 170.

Lewis, J. S., and Wiesner, P.J. (1980) Gonorrhea: Current laboratory methods. Lab. Management, Sept., p.33.

### SABOURAUD DEXT. AGAR

Sabouraud, R. (1910). Les Teignes Paris. Pagano, J., Levin, J.D. and Trejo, W. (1957-8). Diagnostic medium for the differentiation of species of *Candida*. Antibiotics Annual, 137-143.

## IN VITRO MANUFACTURER'S DATA

### G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: [bioprep1@otenet.gr](mailto:bioprep1@otenet.gr)

[www.bioprepare.gr](http://www.bioprepare.gr)