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Bergey's Manual of Systematic Bacteriology, Vol. 1, (1984). Williams and Wilkins, Baltimore/London.

MacFadden, J.F. (1983). Biochemical Tests for the Identification of Medical Bacteria, 2nd edn. Williams and Wilkins, Baltimore/London.

Maslen L.G.C. (1952). Routine use of liquid urea medium for identifying Salmonella and Shigella organisms. J. Brit. Med. 2: 545-546.

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## DESCRIPTION

Identifying agar used to differentiate bacteria that can break down urea.

## PRINCIPLE OF THE METHOD

The addition of agar, peptone and dextrose to the material allows for better growth of gram-negative grains with greater urease production. Urea decomposition causes the release of ammonia which increases the pH of the material to alkaline.

FORMULA	g/litre
Gelatin Peptone	1.0
Dextrose	1.0
Sodium chloride	5.0
Monopotassium phosphate	2.0
Phenol red	0.012
Urea	20.0
Bacteriological Agar	15.0

Appearance: yellow-pale pink clear.

Final pH  $6.8 \pm 0.2$  at 25 °C.

## PRECAUTIONS

UREA AGAR is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

## STORAGE CONDITIONS

The tubes should be stored at **2 – 12 °C** inside their packaging until use. Freezing, even momentarily, destroys the material. Excessive heating should also be avoided. The tubes can be used until the expiration date indicated on the label. For transportation, stability studies have shown that the tubes can remain at **6 – 25 °C** for **7 days** or at **25 – 40 °C** for **48 hours** without affecting the performance of the product.

## INSTRUCTIONS FOR USE

Inoculate the slanted tubes with a pure culture. Incubate the tubes, with loose caps, for 24 to 48 hours or up to 4 days at 35 – 37 °C in an aerobic atmosphere. A positive reaction appears with growth and a purple color in the material. A negative reaction shows no growth or color change in the material.

## GENERAL CHARACTERISTICS OF QUALITY CONTROL

Microorganism	ATCC	Urease reaction
<i>Enterobacter aerogenes</i>	13048	-
<i>Escherichia coli</i>	25922	-
<i>Klebsiella pneumonia</i>	13883	+ Purple material
<i>Proteus vulgaris</i>	13315	+ Purple material
<i>Salmonella typhimurium</i>	14028	-

Proteus vulgaris (+)

#### WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

#### SPECIFICATIONS

UREA AGAR (CHRISTENSEN MEDIUM) - 

PRODUCT	CODE	PACKING	STORE	SELF LIFE
Tube 7ml	070124	40 pieces	2 – 12 °C	8 months
Tube 3ml	080124	20 pieces	2 – 12 °C	8 months

Produced in Greece by Bioprep in accordance with the requirements of the European Directive 98/79/EK. FEK B2198/2-10-2009. CODE EDMA 14 01 02 90. Bioprep company has been certified according to standards EN ISO 9001:2015 / ΕΛΟΤ EN ISO 13485:2016 DY8d/1348/2004

#### BIBLIOGRAPHY

Christensen, W.B. (1946). Urea decomposition as a means of differentiating Proteus and paracolon cultures from each other and from Salmonella and Shigella types. J. Bacteriol. 52: 461-466.

#### IN VITRO MANUFACTURER'S DATA

**Bioprep**  
microbiology

#### G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: [bioprep1@otenet.gr](mailto:bioprep1@otenet.gr)

[www.bioprep.gr](http://www.bioprep.gr)