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## TECHNICAL DATA SHEET

PRODUCT: **UREA BROTH**   
 REFERENCE: **080125**

Date 1st Edition:  
7th 2009  
Date 4th Revision:  
6th 2024

### DESCRIPTION

Urea Broth is used for differentiating organisms, especially Enterobacteriaceae, based on urease production.

### PRINCIPLE OF THE METHOD

When organisms break down urea, ammonia is produced, which causes an alkaline reaction in the medium, resulting in a pink-red color. Therefore, urease is detected by the change in color of the phenol red indicator.

FORMULA	g/litre
Yeast Extract	0.1
Disodium phosphate	9.5
Monopotassium Phosphate	9.1
Phenol red	0.01
Urea	20.0

Appearance: Orange clear.

Final pH  $6.8 \pm 0.2$  at 25 °C.

### PRECAUTIONS

UREA BROTH is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

### STORAGE CONDITIONS

The tubes must be stored at **2 – 12 °C** inside their packaging until use. Freezing, even momentarily, will destroy the material. Excessive heating should also be avoided. The tubes can be used until the expiration date indicated on the label. For transport, stability studies have shown that the tubes can remain at **6 - 25 °C** for **7 days** or at **25 - 40 °C** for **48 hours** without affecting the product's performance.

### USAGE METHOD

UREA BROTH is not intended for primary isolation. Inoculate the material by taking a single colony (a small amount of the colony). Inoculate into the broth and gently mix until the colony dissolves. Leave the cap loosely on the tube. Incubate aerobically for 18-24 hours at 35 - 37 °C.

### INTERPRETATION OF RESULTS

Observe for color changes at 4, 8, 12, and 24 hours. In a positive reaction, the color of the medium turns purple. In a negative reaction, the color of the medium does not change.

### LIMITATIONS OF THE METHOD

If not stored at 6 - 12 °C, urea has been shown to undergo auto-hydrolysis, which could affect the urease reaction.

### GENERAL CHARACTERISTICS OF QUALITY CONTROL

MICROORGANISMS	ATCC	Urease
<i>Enterobacter aerogenes</i>	13048	-
<i>Escherichia coli</i>	25922	-
<i>Klebsiella pneumonia</i>	13883	+ Mωβ
<i>Proteus vulgaris</i>	13315	+ Mωβ
<i>Salmonella typhimurium</i>	14028	-

*Proteus vulgaris*

Pink (+)

#### WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

#### SPECIFICATIONS

##### UREA BROTH - CE

PRODUCT	CODE	PACKING	STORE	SELF LIFE
Tube 2ml	080125	20 pieces	2 – 12 °C	8 months

Produced in Greece by Bioprep in accordance with the requirements of the European Directive 98/79/EK. FEK B2198/2-10-2009. CODE EDMA 14 01 02 01. Bioprep company has been certified according to standards EN ISO 9001:2015 / ΕΛΟΤ EN ISO 13485:2016 DY8d/1348/2004

#### BIBLIOGRAPHY

Rustigian and Stuart. Proc. Soc. Exp. Biol. and Med. 47:109, 1941. Stuart, Van Stratum and Rustigian. J. Bact. 48:437. 1945. McKay, Edwards and Leonar A. J. Clin. Path. 17:479, 1947. Gordon and Mihn. J. Gen. Microbiol., 21:736. 1959. Goldsmith and Latlief. Applied Microbiol., 3:195. 1955.

#### IN VITRO MANUFACTURER'S DATA

#### G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: [bioprep1@otenet.gr](mailto:bioprep1@otenet.gr)

[www.bioprep.gr](http://www.bioprep.gr)