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Bergey's Manual of Systematic Bacteriology, Vol. 1, (1984). Williams and Wilkins, Baltimore/London.

MacFadden, J.F. (1983). Biochemical Tests for the Identification of Medical Bacteria, 2nd edn. Williams and Wilkins, Baltimore/London.

Maslen L.G.C. (1952). Routine use of liquid urea medium for identifying Salmonella and Shigella organisms. J. Brit. Med. 2: 545-546.

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**TECHNICAL DATA SHEET**

PRODUCT: **UREA INDOLE BROTH**  
 REFERENCE: **080183, 070183**

Date 1st Edition:  
7th 2009  
 Date 4th Revision:  
6th 2024

**DESCRIPTION**

UREA INDOLE BROTH is used for testing the breakdown of urea and the production of indole.

**PRINCIPLE OF THE METHOD**

The medium contains 5% tryptophan, which some bacteria (e.g., *E. coli*) can break down to produce indoles. It also contains urea, which some bacteria (e.g., *P. mirabilis*) can break down. The production of indole is visible when 3-4 drops of Kovac's reagent are added to the medium, forming a ring on the surface of the broth. Initially, the ring is green and yellow, but with the presence of indole, it changes color to red violet within 5-10 minutes. When the urea is broken down, the pH of the medium becomes alkaline, and due to the phenol red indicator, the color changes from pale pink to bright pink. Peptone enhances the performance of the medium, while sodium chloride regulates the osmotic balance of the medium.

COMPOSITION	g/litre
Tryptone	5.0
Peptone	5.0
Sodium chloride	5.0
Urea	20.0
Phenol red	0,004

Appearance: Beige.

Final pH  $6.9 \pm 0.2$  at 25 °C.

**PRECAUTIONS**

UREA INDOLE BROTH is an in-vitro laboratory diagnostic material and should only be handled by qualified people in the laboratory. This material contains peptones and extracts of animal origin. The certificates regarding the origin and health status of the animals do not fully guarantee the absence of transmissible pathogens. For this reason, it is recommended that these materials be treated as potentially infectious, with the usual safety precautions (avoiding ingestion or inhalation). Plates should always be handled with gloves and in Laminar flow Class II, to avoid contamination mainly by saprophytic fungi. If the plate is cracked or the bag has a hole, do not use it. Do not use petri dishes if there are signs of microbial contamination. The thickness of the agar must be 4 - 5 mm and the material without cracks, dryness or other signs of deterioration. After the expiry date the material is unfit for use. In case of contact with the skin, wash immediately with plenty of water and soap. Positive samples must be destroyed according to the hygienic rules prescribed for the management of contaminated samples.

**STORAGE CONDITIONS**

The tubes should be stored at **2 – 12 °C** in their packaging until they are used. Freezing, even momentarily, destroys the material. Excessive heating should also be avoided. The tubes can be used until the expiration date printed on the label. For transport, stability studies have shown that the tubes can remain at **6 – 25 °C** for **7 days** or at **25 – 40 °C** for **48 hours** without affecting the product's performance.

**USAGE INSTRUCTIONS**


Inoculate the broth by taking enough of a single colony from a pure culture that has been incubated for at least 24 hours. Shake the tubes well to dissolve the bacteria. Incubate the tubes with loose caps at 35 – 37 °C for 18 to 24 hours.

**INTERPRETATION AND RESULTS**

If the color of the medium changes from beige to bright pink, the urea test is positive (+). For indole testing, add 3 – 4 drops of Kovac's reagent. If the color of the ring changes from green-yellow to red-violet (within 5 – 10 minutes), the indole test is positive (+). If the color of the ring remains green-yellow, the indole test is negative (-).

**GENERAL CHARACTERISTICS OF QUALITY CONTROL**

GERM	INDOLE	UREA
<i>Escherichia coli</i> ATCC 25922	(+) Red ring	(-) Beige material
<i>Klebsiella pneumoniae</i> ATCC 13883	(-) Yellow ring	(+) Purple material
<i>Salmonella typhimurium</i> ATCC 14028	(-) Yellow ring	(-) Beige material
<i>Proteus mirabilis</i> ATCC12473	(-) Yellow ring	(+)Purple material

*E. coli* (Indole +)

*P. mirabilis* (Urea +)

#### WASTE DISPOSAL OF WASTE

Materials that show no growth can be considered as non-hazardous waste and disposed of accordingly. Materials that show colony growth must be disposed of according to the guidelines for infectious or potentially infectious waste. The laboratory is responsible for the proper management of infectious waste according to its nature and level of risk and must handle and dispose of it (or assign its management and disposal) in compliance with the applicable regulations.

#### SPECIFICATIONS

UREA INDOLE BROTH 

ITEM	CODE	PACKAGE	KEEPING	SELF LIFE
Tube 2ml	080183	10 pieces	2 – 12 °C	8 months

It is produced in Greece from the Bioprep Company in accordance with the requirements of the European guidance 98/79/EK. ΦΕΚ Β2198/2-10-2009. Code by EDMA 14 01 02 01 Bioprep Company has been certified according to the standards EN ISO 9001:2015 / ΕΛΟΤ EN ISO 13485:2016 DY8d/1348/2004

#### BIBLIOGRAPHY

Bergey's Manual of Systematic Bacteriology, Vol. 1, (1984). Williams and Wilkins, Baltimore/London.  
 MacFadden, J.F. (1983). Biochemical Tests for the Identification of Medical Bacteria, 2nd edn. Williams and Wilkins, Baltimore/London.  
 Maslen L.G.C. (1952). Routine use of liquid urea medium for identifying Salmonella and Shigella organisms. J. Brit. Med. 2: 545-546.

#### IN VITRO MANUFACTURER'S DATA

#### G. PAPANIKOLAOU & CO

PRODUCTION LABORATORIES OF CULTURE MEDIA

Potamou 5, Industrial Area Keratea, Attica

P.O. Box: 4893, Postal Code: 9001 - Tel: +30 2299066113. Fax: +30 2299066112

E-mail: [bioprep1@otenet.gr](mailto:bioprep1@otenet.gr)

[www.bioprep.gr](http://www.bioprep.gr)